



AUG 31, 2023

## 🌐 Qiagen QIAmp PowerFecal Pro extraction kit

Michael Dan                      Dilan  
Siemon<sup>1</sup>,                      Joseph<sup>1</sup>,                      Jessica Pardy<sup>1</sup>,  
Richard Gibson<sup>1</sup>, christopher.degroot<sup>1</sup>

<sup>1</sup>Western University



Michael Dan Siemon  
Western University

### ABSTRACT

This protocol is used by Western University to process wastewater samples for wastewater-based epidemiology. The protocol is modified from the protocol described in the RNeasy PowerFecal Pro Kit Handbook provided by Qiagen with the PowerFecal Pro kit.

### MATERIALS

This protocol requires the Qiagen QIAmp PowerFecal Pro extraction kit

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.j8nlko8owv5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlko8owv5r/v1)

**Protocol Citation:** Michael Dan Siemon, Dilan Joseph, Jessica Pardy, Richard Gibson, christopher.degroot 2023. Qiagen QIAmp PowerFecal Pro extraction kit.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlko8owv5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

This protocol is used for wastewater RNA extraction for wastewater-based epidemiology at Western University

**Created:** Aug 30, 2023

**Last Modified:** Aug 31, 2023

**PROTOCOL integer ID:**  
87125

**Keywords:** wastewater-based epidemiology, RNA extraction, Qiagen, ImPaKT

## Sample Preparation

- 1 Thoroughly mix wastewater sample then aliquot 40 mL into 50 mL Falcon tube. Centrifuge at 4500 x g for 120 min. Decant supernatant, assume 280 µl pellet.
- 2 Spin the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom. Add 100 µl of phenol-chloroform-isoamyl alcohol to the PowerBead Pro Tube.
- 3 Add 650 µl of Solution CD1 to the wastewater pellet and transfer solution to the PowerBead Pro Tube.
- 4 As a substitute for vortexing described in the kit protocol, bead-beating is used. Bead-beating is conducted 4x for 30s at 4 m/s.
- 5 Centrifuge the PowerBead Pro Tube at 15,000 x g for 1 min.
- 6 Transfer the supernatant to a clean 2 ml Microcentrifuge Tube. Add 200 µl of Solution CD2 and vortex for 5 s.

- 7** Centrifuge at 15,000 x g for 1 min. Avoiding the pellet, transfer up to 650 µl of supernatant to a clean 2 ml Microcentrifuge Tube.
- 8** Add 650 µl of Solution EA. Vortex briefly.
- 9** Load 650 µl supernatant-EA mix into an MB RNA Spin Column and centrifuge at 15,000 x g for 1 min. Discard the flow-through and repeat until all of the solution has been passed through the Spin Column.
- 10** Add 650 µl Solution EA to the Spin Column and centrifuge at 15,000 x g for 1 min.
- 11** Mix 45 µl DNase Digestion Solution and 5 µl DNase I stock enzyme. Place the MB RNA Spin Column into a clean 2 ml Collection Tube and add the DNase Digestion Solution to the center of the filter. Incubate at room temperature for 15 minutes.
- 12** Add 650 µl Solution EA to the Spin Column and centrifuge at 15,000 x g for 1 min.
- 13** Discard the flow-through. Add 500 µl Solution C5. Centrifuge at 15,000 x g for 1 min.
- 14** Replace the collection tube and dry the Spin Column by centrifuging at 20,000 x g for 1 min.
- 15** Elute the RNA by placing the spin column in a 1.5 ml collection tube. Add 100 µl of RNase-free water and centrifuge at 15,000 x g for 1 min.

