

Oct 22, 2025

Molecular cloning of SHIP164 plasmids for expression in mammalian cells

DOI

dx.doi.org/10.17504/protocols.io.8epv5z5kjb/v1

michael.hanna¹, Michael G Hanna^{2,3}, Pietri De Camilli^{2,3}

¹Yale University;

²Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Michael G Hanna

Yale University

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5z5kjb/v1>

Protocol Citation: michael.hanna , Michael G Hanna, Pietri De Camilli 2025. Molecular cloning of SHIP164 plasmids for expression in mammalian cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5z5kjb/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: July 18, 2021

Last Modified: October 22, 2025

Protocol Integer ID: 51653

Keywords: Molecular cloning of SHIP164 plasmids, expression in mammalian cells, ASAPCRN, molecular cloning of ship164 plasmid, molecular cloning of ship164, ship164 plasmid, molecular cloning, expression in mammalian cell, ship164, mammalian cell, sequences of other protein, other protein, protein, cell

Abstract

This protocol is to help with the molecular cloning of SHIP164 and the sequences of other proteins.

Attachments




Molecular cloning of...

46KB

Troubleshooting

Safety warnings

 For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Before start

We experienced difficulty in amplifying SHIP164 from a cDNA library, therefore we opted to generate a codon optimized SHIP164 open reading frame synthesized. The codon optimized cDNA was cloned into pUC57 and was subsequently subcloned into multiple other mammalian expression vectors.

The general cloning strategy is described below.



Cloning strategy

- 1 Linearize destination vector using restriction enzyme digest following manufacturers protocol (New England Biolabs).
- 2 Amplify desired insert by PCR making sure to include 15-20 bp overhang homology with destination vector.
- 3 Run both cut vector and PCR product on an agarose gel (e.g., 1% supplemented with Midori Green or Ethidium Bromide for visualization) and purify the DNA from gel using a NucleoSpin Gel and PCR Clean-up kit (Takara).
- 4 Ligate PCR insert into vector using NEB HiFi assembly mixture following manufacturers protocol (<https://www.neb.com/protocols/2014/11/26/nebuilder-hifi-dna-assembly-reaction-protocol>).