



1 ▼

Oct 06, 2021

# DNA Ethanol Precipitation (SOP009.v1.1) V.1

Rory Kruithoff<sup>1</sup>, Douglas Shepherd<sup>1</sup><sup>1</sup>Arizona State University

1



protocol .

Human Cell Atlas Method Development Community

Rory Kruithoff  
Arizona State University

**Document Summary:** This document, SOP002 - DNA Ethanol Precipitation, describes a method to concentrate or dry DNA using ethanol to reduce the solubility of dissolved DNA causing it to precipitate out of solution. Alternatively, a vacuum concentrator can be used if one is available.

[DNA Ethanol Precipitation,  
SOP009.v1.1.pdf](#)

Rory Kruithoff, Douglas Shepherd 2021. DNA Ethanol Precipitation  
(SOP009.v1.1). **protocols.io**  
<https://protocols.io/view/dna-ethanol-precipitation-sop009-v1-1-byc7pszn>



protocol

From Qiagen: <https://www.qiagen.com/it/resources/faq?id=5d591b8b-968a-4a17-849f-9d0f719b40af&lang=en&Print=1>

DNA, ethanol, ethyl alcohol, precipitation

protocol ,

Sep 19, 2021

Oct 06, 2021

Sep 19, 2021



zahara93.zs



Sep 28, 2021



Rory Kruithoff Arizona State University

53375

### Required Reagents

-  [Sodium](#)
- [Acetate Sigma Catalog #S2889-250g](#)
-  [Ethanol 99.5% ACS](#)
- [Reagent Thermofisher Catalog #AC615090010](#)
- Nuclease-free water

### Required Equipment

- Microcentrifuge
- Sterile Eppendorf tubes
- -20°C Freezer

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

### Quick Overview:

Step 1 – Precipitate DNA

Step 2 – Pellet DNA and dry

Step 3 – Resuspend DNA

### v1.1 revision notes





1. Updated SOP to standard formatting
2. Updated document description

#### Precipitate DNA 2h

- 1 Add 1/10 volume of 3 M Na-Acetate [pH 5.2](#), and 2 to 2.5 volumes of ice-cold 100% ethanol to the DNA sample.




2h



Mix, and store at  **-20 °C** for at least  **01:00:00** to precipitate the DNA. Typically, DNA will be left at  **-20 °C**  **Overnight** before proceeding to part 2.



#### Pellet DNA and dry 35m



20m

3 


Recover the precipitated DNA by centrifugation at full speed in a microcentrifuge for  **00:15:00** -  **00:20:00** to form a pellet.

4 Pour off the ethanol and wash the pellet twice with  **Room temperature**  
 **70 % (v/v) ethanol** .

5 Allow the DNA pellet to air-dry.

#### Resuspend DNA

6 Re-suspend the DNA in a suitable volume of sterile TE buffer or nuclease-free water.

7 Store at  **-20 °C** and avoid any unnecessary freeze-thaw cycles.