



Aug 12, 2022

DNA Purification (NEB)

Brian Teague¹¹University of Wisconsin - Stout1 *Works for me* Share

This protocol is published without a DOI.

Yeast ORFans CURE**Brian Teague**
University of Wisconsin - Stout

ABSTRACT

This protocol purifies DNA from an enzymatic reaction (like a PCR or a digestion with a restriction enzyme.) It is NOT an miniprep -- if you are trying to purify DNA from a culture of E. coli, you are in the wrong place!

PROTOCOL CITATION

Brian Teague 2022. DNA Purification (NEB). **protocols.io**
<https://protocols.io/view/dna-purification-neb-ce6sthee>



KEYWORDS

pcr, purification, monarch, spin column

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 11, 2022

LAST MODIFIED

Aug 12, 2022

PROTOCOL INTEGER ID

68530

MATERIALS TEXT

- Binding buffer
- Wash buffer
- Elution buffer
- One spin column and collection tube per cleanup
- Chemical waste container (50 ml conical)

The buffers and spin column are from the

[☒ Monarch PCR and DNA Cleanup Kit - 50 preps New England](#)

Biolabs Catalog #T1030S

SAFETY WARNINGS

The binding buffer may cause irritation to skin and eyes. Additionally, we are shedding nucleases -- enzymes that degrade DNA -- all the time. Wear lab coats, gloves and safety glasses.

The flow-through (containing the binding buffer or wash buffer) cannot go down the drain. Dispose of it as chemical waste, per the directions of your instructor.


1 

Make sure you are using the "Monarch PCR & DNA Cleanup Kit", not the Monarch Plasmid Miniprep Kit. They come in identical boxes -- read the label!






2 Mix the ENTIRE DNA sample with 5 times its volume of Binding Buffer. Pipette up and down several times to mix the sample thoroughly.

le, if you had  **200 µL** of sample, mix it with  **1000 µL** of binding buffer.

3 Insert the column into the collection tube and load the sample onto the column.

4 Centrifuge  **16000 x g, 00:01:00** . Discard the flow-through in the **chemical waste container**.

1m

- 5 Re-insert the column into the collection tube. Add  **200 µL** of Wash Buffer.
- 6 Centrifuge  **16000 x g, 00:01:00** . Discard the flow-through in the **chemical waste container**. 1m
- 7 Repeat steps 5 and 6 once.
- 8 Transfer the column to a clean 1.7 ml microcentrifuge tube.
- 9 Add  **10 µL** of Elution Buffer to the center of the silica matrix. Wait  **00:01:00** , then 2m
centrifuge  **16000 x g, 00:01:00**
- 10 Quantify your purified DNA on the Nanodrop.