



MAR 19, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.5qpvox3bl4o/v1

Protocol Citation: Nisha Mohd Rafiq, Pietro De Camilli 2024. Cell culture, transfection, immunocytochemistry, and imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvox3bl4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 12, 2024

🌐 Cell culture, transfection, immunocytochemistry, and imaging

Nisha Mohd Rafiq¹, Pietro De Camilli^{2,3,4,5}

¹University of Tuebingen;

²Department of Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

³Department of Cell biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

⁴Program in Cellular Neuroscience, Neurodegeneration and Repair.;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA.

ASAP Collaborative Research Network



Berrak Ugur
Yale University

ABSTRACT



This protocol describes the maintenance, transfection, immunocytochemistry, and imaging of RPE1 and also transfection, immunocytochemistry, and imaging of iPSCs, i³ Neurons and DA neurons.

ATTACHMENTS

[Rafiq_transfection_imaging.docx](#)

MATERIALS

Reagents:

-  Lipofectamine™ 2000 Transfection Reagent **Thermo Fisher Scientific Catalog #11668019**
-  Lipofectamine™ Stem Transfection Reagent **Thermo Fisher Scientific Catalog #STEM00008**

Last Modified: Mar 19, 2024

PROTOCOL integer ID: 96916


Keywords: ASAPCRN

Funders Acknowledgement:

ASAP

Grant ID: ASAP-000580

General cell culture for RPE1








- 1 Grow hTERT-RPE1 cells in DMEM/F12 (Thermo Fisher Scientific) supplemented with 10% FBS (Thermo Fisher Scientific), 1% glutaMAX, 1% penicillin-streptomycin. Keep cells at  37 °C with 5% CO2 in an enclosed incubator.

Note

For general maintenance, when cells reach 80-90% confluency, detach them from the dish with Trypsin and dilute 1:10-20 in a new dish.

Cell transfection for RPE1





5d

- 2 For live-cell imaging experiments, seed the cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $1-2 \times 10^5$ cells.
- 3 For RPE1: allow the cells to adhere for  08:00:00 -  24:00:00 before being transiently transfected using  4 μ L Lipofectamine™ 2000 Transfection Reagent (Invitrogen) in Opti-MEM media, mix them with the respective plasmids ( 1 μ g -  2 μ g) and visualize after  48:00:00 . 3d
- 4 For cilia generation, serum-starve the cells in DMEM/F12 media (without FBS) for  48:00:00 . 2d

Cell transfection for iPSCs, i³Neurons and DA Neurons

2d

5 For live-cell imaging experiments, seed cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $3-5 \times 10^5$ cells.

6 Use  2.4 μ L of Lipofectamine™ Stem Transfection Reagent (Invitrogen) in Opti-MEM media with respective plasmids ( 1 μ g -  3 μ g) and visualize at least  48:00:00 later.


2d

Immunocytochemistry of RPE1, iPSCs, i³Neurons and DA Neurons

2h 7m


7 For fixed-cell imaging experiments, seed cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $1-3 \times 10^5$ cells.

20m

- For immunofluorescence visualization, fix the cells with 4% (v/v) paraformaldehyde (Electron Microscopy Sciences) in 1× phosphate-buffered saline (PBS) for  00:20:00 .



8 Wash cells thrice in PBS.



9 Perform cell extraction with 0.25-0.5% (v/v) Triton X-100 in PBS for  00:10:00 .

10m

10 Wash cells thrice in PBS.

11 For removal of free aldehyde groups, quench the cells with fresh  1 undetermined sodium borohydride (Sigma-Aldrich) in PBS for  00:07:00 , and then wash thrice in PBS.

7m

12 Further block the cells for 00:30:00 in 5% bovine serum albumin (BSA, Sigma-Aldrich) in PBS and then incubate Overnight at 4 °C with the respective antibodies listed in Table S1. 30m

13 Wash the cells with PBS thrice the following day and incubate with Alexa Fluor-conjugated secondary antibodies (Thermo Fisher Scientific) for 01:00:00 at Room temperature . 1h

14 Wash the cells thrice in 1×PBS.

15 Use DAPI (Thermo Fisher Scientific) for nuclear staining, when necessary.

Imaging

16 For live imaging, maintain cells in a caged incubator with humidified atmosphere (5% CO₂) at 37 °C .

Note

The Yokogawa spinning disk field scanning confocal system with microlensing (CSU-W1 SoRa, Nikon) controlled by NIS elements (Nikon) software was used for imaging. Excitation wavelengths between 405-640 nm, CFI SR Plan ApoIR 60XC WI objective lens and SoRa lens-switched light path at 1x, 2.8x or 4x were used. SoRa images were deconvolved using the Batch Deconvolution (Nikon) software.