



APR 12, 2024

3,3-Diaminobenzidine Tetrahydrochloride (DAB) Immunohistochemistry on Mouse Brain Tissue Sections

madalynn.erb Erb¹

¹Van Andel Research Institute

ASAP Collaborative Research Network

Team Lee



Jane Balster
ASAP - Team Lee

OPEN ACCESS



ABSTRACT

This protocol details the staining for free floating mouse brain sections.

DOI:

dx.doi.org/10.17504/protocols.io.3byl497zjgo5/v1

Protocol Citation: madalynn.erb Erb 2024. 3,3-Diaminobenzidine Tetrahydrochloride (DAB) Immunohistochemistry on Mouse Brain Tissue Sections.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.3byl497zjgo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Apr 08, 2024

Last Modified: Apr 12, 2024

PROTOCOL integer ID: 98141

Keywords: ASAPCRN

MATERIALS

Pyrex® Staining Dish **Ted Pella Inc. Catalog #36754-60**

8-Section Staining Nets **Ted Pella Inc. Catalog #36154-64**

VECTASTAIN® Elite® ABC-HRP Kit, Peroxidase (Standard) **Vector Laboratories Catalog #PK-6100**

Pierce™ DAB Substrate Kit **Thermo Scientific Catalog #34002**

Fisherbrand™ Superfrost™ Plus Microscope Slides **Fisher Scientific Catalog #12-550-15**

Entellan **Electron Microscopy Sciences Catalog #14802**

Recipes

10X TBS

A	B
1M Tris	500 ml
5M NaCl	300 ml
H2O	200 mL

1X TBS + 0.1% triton → 1 liter

A	B
10X TBS	100 ml
Triton X	1 ml
Water	899 ml

TBS → 1 liter

A	B
1M Tris	50 ml
5M NaCl	30ml
Water	919 ml

TBST (0.1% Tween-20) → 1 liter

A	B
10X TBS	100 ml
Tween-20	1ml
Water	899 ml

0.2M Phosphate Buffer

A	B
Na ₂ HPO ₄	22.72 g
NaH ₂ PO ₄	5.52 g
H ₂ O	1000 mL

0.1M Phosphate Buffer

A	B
0.2M PB	500 mL
H ₂ O	500 mL

Cresyl Violet (500 mL)

A	B
H ₂ O	500 mL
Cresyl Violet Acetate	0.5g
100% Glacial Acetic Acid	1.25mL

Warm to 60°C and stir until completely dissolved filter with whatman paper (wrap in foil to protect from light)

Quenching solution→100ml

A	B
30% hydrogen peroxide (Sigma H1009)	1 ml
Methanol	99 ml

ABC mix→60ml

A	B
TBS	60 ml
Reagent A	24 drops
Reagent B	24 drops

DAB Solution→30ml

A	B
PBS	30 ml
(DAB reagent 1)	12 drops buffer pH = 7.5
(DAB reagent 2)	24 drops DAB
(DAB reagent 3)	12 drops hydrogen peroxide

vortex between reagents

Blocking solution→100ml

A	B
NGS	10 ml
TBS + triton	90 ml + 0.1%

Primary Antibodies:

A	B	C	D
Target	Species	Conc	Manufacturer
Iba1	Rabbit	1:1000	WAKO 019-19741
GFAP	Mouse	1:1000	Sigma G3893
pSer129-Synuclein	Rabbit	1:1000	Abcam ab51253
TH	Rabbit	1:2000	Novus Biological N300109
pSer202, pThr205 Tau (AT8)	Mouse	1:1000	Thermofisher MN1020
GFP	Rabbit	1:1000	Thermofisher A-11122

Secondary Antibodies:

A	B	C	D	E
Target	Species	Conc	Manufacturer	
Mouse	Goat	1:500	BA-9200-1.5	Biotinylated
Rabbit	Goat	1:500	BA-1000-1.5	Biotinylated

Day 1

- 1 Staining protocol for 35µm free floating mouse brain sections.
- 2 Staining is performed in glass staining dishes (Pyrex 36754-60) using 8-section staining nets (Ted Pella 36154-64) .

2.1 The sections can be transferred between wells using a paint brush.

2.2 Volume of solution:



1. 20 mL - 30 mL for antibodies or ABC mix .
2. 50 mL for washing.

3 Gently rock the plates during the washing and incubation steps.

4 Wash sections: 3 times (5 min per wash) in Tris Buffer Saline (TBS) at Room temperature .



4.1



Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (1/3).

5m

4.2



Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (2/3).

5m

4.3



Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (3/3).

5m

5



Quenching: Incubate sections in Quenching Solution (0.3% H_2O_2 in Methanol) for 00:10:00 at 4 °C .

10m

6



Wash sections: 3 short rinses in TBS then 2 times 5min in TBS at Room temperature .

6.1



Then wash for 00:05:00 in TBS at Room temperature (1/2).

5m



6.2

Then wash for 00:05:00 in TBS at Room temperature (2/2).

5m



7 Block sections in blocking solution (10% Normal Goat Serum (NGS) + 0.1% Triton-X in TBS) for

1h

01:00:00 at Room temperature .



8 Primary Antibodies: Dilute antibodies in 5% NGS + 0.1% Triton-X in TBS.

8.1

Incubate in primary antibody solution for 48:00:00 at 4 °C .

2d



8.2

Cover dishes with parafilm for this step.



Day 4

1d 3h 47m

9 Wash 3 times (10 min each wash) in TBST (0.1% Tween 20 in TBS) at Room temperature .





9.1

Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (1/3).

10m



9.2

Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (2/3).

10m

9.3

Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (3/3).

10m



10

Secondary Antibodies:







10.1

Dilute secondary antibodies in TBST.

10.2

Biotinylated goat anti-rabbit IgG (1:500 dilution) or Biotinylated goat anti-mouse IgG (1:500 dilution).


10.3

Incubate for  02:00:00 at  Room temperature or  24:00:00 at  4 °C .

1d 2h









11 Wash sections 3 times (10min each wash) in TBST at  Room temperature .

11.1 Wash for  00:10:00 in TBST at  Room temperature (1/3).

10m



11.2 Wash for  00:10:00 in TBST at  Room temperature (2/3).

10m

11.3 Wash for  00:10:00 in TBST at  Room temperature (3/3).

10m

12 Prepare ABC solution using VECTASTAIN® Elite® ABC-HRP Kit, Peroxidase (Standard) (PK-6100).


12.1 ABC mix:  5 mL TBS + 2 drops of reagent A (vortex) + 2 drops reagent B (vortex) let stand  30m



 Room temperature for  00:30:00 before using.

12.2 Vortex between reagents (after adding reagent A and after adding reagent B).



12.3 Make  30 mL ABC solution per dish.

13

Incubate sections in ABC solution for 01:00:00 at Room temperature .

1h



14

Wash sections 3 times (10min each wash) in TBST at Room temperature .



14.1

Wash for 00:10:00 in TBST at Room temperature (1/3).

10m



14.2

Wash for 00:10:00 in TBST at Room temperature (2/3).

10m



14.3



Wash for 00:10:00 in TBST at Room temperature (3/3).

10m





15

Wash sections 2 times (5 min each wash) in TBS at Room temperature .

15.1 Wash sections for  00:05:00 in TBS at  Room temperature (1/2).

5m





15.2 Wash sections for  00:05:00 in TBS at  Room temperature (2/2).

5m



16 Prepare DAB solution using DAB Substrate Kit, Peroxidase (HRP), with Nickel, (3,3'-diaminobenzidine) (SK-4100).



1.  5 mL TBS + 2 drops of Reagent 1 (Buffer pH=7.5) + 4 drops Reagent 2 (DAB) + 2 drops Reagent 3 (H2O2).
2. Vortex between reagents (after reagent 1, after reagent 2 and after reagent 3).
3. Make  30 mL DAB solution per dish.

17 Revelation with DAB chromogen (until cells are visible):

2m

Note

typically about  00:02:00 for TH.

- All dishes must be developed for the same amount of time – do them together, ask a friend for help if you have more than 2 dishes of tissue.

18 Perform 3 short rinses in TBS at  Room temperature .



18.1 Perform rinse in TBS at Room temperature (1/3).



18.2 Perform rinse in TBS at Room temperature (2/3).



18.3 Perform rinse in TBS at Room temperature (3/3).



19 Wash sections 5 time (5 min each wash) in TBS at Room temperature .



19.1 Wash sections for 00:05:00 in TBS at Room temperature (1/5).





5m

19.2 Wash sections for 00:05:00 in TBS at Room temperature (2/5).



5m



19.3

Wash sections for  00:05:00 in TBS at  Room temperature (3/5).

5m



19.4

Wash sections for  00:05:00 in TBS at  Room temperature (4/5).

5m



19.5

Wash sections for  00:05:00 in TBS at  Room temperature (5/5).

5m



20





Mount sections:

20.1

Use a petri dish filled with PBS to mount the sections on superfrost plus slides (Fisher Scientific 12-550-15).



- This is easiest to do using a medium sized paint brush.

20.2

After mounting let slides dry at  Room temperature ( 00:02:00 -  00:10:00) then  12m in miliQ H₂O.



20.3

Dry slides  Overnight at  Room temperature before moving to the next step.

5m





Day 5 or 6

9h 18m

21 Rehydrate slides in H₂O for 00:02:00 - 00:03:00 at Room temperature .

5m



22 Dehydrate slides in sequential EtOH washes (70%, 95%, 100% EtOH) for 00:05:00 each wash at Room temperature .

5m



23 Clear tissue for 00:05:00 in fresh Xylene I then incubate in fresh Xylene for ≥ 00:05:00 .

10m





- Keep slides in Xylene until coverslipping.
- Perform Xylene steps and coverslipping in a chemical fume hood.



24 Use Entellan mounting medium (Electron Microscopy Sciences 14802) to coverslip slides.



24.1 Apply one line of Entellan across the middle of the slide.

24.2 Cover with glass coverslip and gently press down with forceps.

24.3 Let slides dry in the fume hood for  01:00:00 – then dry on the bench  Overnight at **1h 5m**




 Room temperature .

24.4 Gently remove excess Entellen with a razor blade.

Optional Cresyl Violet Protocol (NISSL stain)

10m

25 After rehydration and before ethanol steps (dehydration).

26 Incubate slides in fresh cresyl violet for  00:10:00 at  Room temperature .

10m

27 Proceed with dehydration steps normally.

28 Must stain with cresyl violet for TH stereology in substantia nigra.

Secondary Antibodies

29

	A	B	C	D	E
	Target	Species	Conc	Manufacturer	
	Mouse	Goat	1:500	BA-9200-1.5	Biotinylated
	Rabbit	Goat	1:500	BA-1000-1.5	Biotinylated