



AUG 01, 2023

## KAPP-Sen TMC: Dissociation of Pancreatic Islets (recovered) Protocol

Jessica

Juliana Alcoforado Diniz<sup>1</sup>, Dylan Baker<sup>1</sup>, Garofalo<sup>1</sup>, Paul Robson<sup>2</sup>

<sup>1</sup>The Jackson Laboratory; <sup>2</sup>The Jackson Laboratory for Genomic Medicine

Cellular Senescence Network (SenNet) Method Development Community

OPEN  ACCESS



Ashley M Raynock

UConn Health, UConn Center on Aging

### ABSTRACT

The dispersed samples were shipped cold from PRODOLABS. Prior to scRNA-seq dispersed samples from brain dead donor's pancreatic islets were recovered and dissociated as follows.

**Protocol Citation:** Juliana Alcoforado Diniz, Dylan Baker, Jessica Garofalo, Paul Robson 2023. KAPP-Sen TMC: Dissociation of Pancreatic Islets (recovered) Protocol.

**protocols.io**

<https://protocols.io/view/kapp-sen-tmc-dissociation-of-pancreatic-islets-rec-cxz6xp9e>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Aug 01, 2023

**Last Modified:** Aug 01, 2023

## Abstract

- 1 The dispersed samples were shipped from **PRODOLABS**. Prior to scRNA-seq dispersed samples from brain dead donor's pancreatic islets were recovered and dissociated as follows.

## Cell Dissociation with Accutase

- 2 

NOTE: Before beginning cell dissociation coat all the materials (pipettes, tubes, etc.) with PIM-S001GMP media to prevent sticking.

  1. Transfer cell suspension of pure islets to a new 50ml tube. Use additional media to rinse original container.
  2. Centrifuge at room temperature 130g for 5 mins.
  3. Aspirate the supernatant and add media to the appropriate concentration of 1,000 islets/1 ml.  
--> Ex: If receiving 7,000 islets, use 7 ml media.
  4. To recover cells, add half of the islets to a coated flask and incubate at 37C overnight.
  5. Coat pipette and add islets from flask to a 50ml tube. Wash flask with media to make sure all islets have been collected.
  6. Centrifuge at 130 g, room temperature for 2 mins.
  7. Aspirate media and resuspend in 4 ml accutase. Incubate at 37C for 8 mins, mixing with pipette every 2 mins. Check at 6 mins.  
-->1 mL accutase/1,000 islets
  8. Add CMRL 1066 (Cat. 11530037) media to 25 ml then centrifuge at 230 g.
  9. Aspirate supernatant and resuspend in 1.5 ml of CMRL
  10. Filter through a 40 µm Flowmi.
  11. Count cells using AO/PI (acridine orange/propidium iodide) Cell Viability Kit for Luna-FL automated cell counter.
  12. Proceed to cell fixation.

## Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling

- 3 Cells were fixated prior to scRNAseq according to [https://dx.doi.org/10.17504/protocols.io.\[...\]/v1](https://dx.doi.org/10.17504/protocols.io.[...]/v1)