



APR 06, 2023

OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.36wgqjerovk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqjerovk5/v1)

**Protocol Citation:** Hannah M. Alt, David G. Brandner, Gregg T. Beckham, Kelsey J. Ramirez 2023. Lignin Reductive Catalytic Fractionation (RCF) Monomers Analysis by Gas Chromatography Flame Ionization Detection (GC-FID).

**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.36wgqjerovk5/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Feb 07, 2023

**Last Modified:** Apr 06, 2023

**PROTOCOL integer ID:**  
 76618

## Lignin Reductive Catalytic Fractionation (RCF) Monomers Analysis by Gas Chromatography Flame Ionization Detection (GC-FID)

Hannah M. Alt<sup>1</sup>, David G. Brandner<sup>1</sup>, Gregg T. Beckham<sup>1</sup>, Kelsey J. Ramirez<sup>1</sup>

<sup>1</sup>Renewable Resources and Enabling Sciences Center, National Renewable Energy Laboratory

NREL

Tech. support email: [stefan.haugen@nrel.gov](mailto:stefan.haugen@nrel.gov)



Hannah M. Alt

### ABSTRACT

A gas chromatography with flame ionization detection (GC-FID) method was developed to quantify reductive catalytic fractionation (RCF) monomers.

### GUIDELINES

#### NOTICE

This work was authored by the National Renewable Energy Laboratory, operated by Alliance for Sustainable Energy, LLC, for the U.S. Department of Energy (DOE) under Contract No. DE-AC36-08GO28308. Funding provided by U.S. Department of Energy Office of Energy Efficiency and Renewable Energy Bioenergy Technologies Office. The views expressed herein do not necessarily represent the views of the DOE or the U.S. Government.

### MATERIALS

#### Woody RCF Standards

 Phenol **Contributed by users**

 Guaiacol **Contributed by users**

 4-Ethylphenol **Contributed by users**

 4-Methylguaiacol **Contributed by users**

 4-Ethylguaiacol **Contributed by users**

 Syringol **Contributed by users**

 4-Propylguaiacol **Contributed by users**

**Keywords:** Reductive  
Catalytic Fractionation, Gas  
Chromatography

- ⊗ Isoeugenol **Contributed by users**
- ⊗ Methylparaben **Contributed by users**
- ⊗ 4-Ethylsyringol **Contributed by users**
- ⊗ 4-Propylsyringol **Contributed by users**
- ⊗ 4-Propanolguaiacol **Contributed by users**
- ⊗ Propenylsyringol **Contributed by users**
- ⊗ 4-Propanolsyringol **Contributed by users**

**Non-Woody RCF Standards**

- ⊗ Methyl-3-(4-hydroxyphenyl) propionate **Contributed by users**
- ⊗ Methyl-3-(4-hydroxy-3-methoxyphenyl) propionate **Contributed by users**
- ⊗ Methylcoumarate **Contributed by users**
- ⊗ Methylferulate **Contributed by users**

**Internal Standard**

- ⊗ Tri-tert-butylbenzene **Contributed by users**

**Reagents**

- ⊗ Methanol **P212121 Catalog #PA-33900HPLCCS4L**


**Equipment**

Equipment	
Gas Chromatograph	NAME
8890 GC System	TYPE
Agilent	BRAND
8890 GC System	SKU
<a href="https://www.agilent.com/en/product/gas-chromatography/gc-systems/8890-gc-system">https://www.agilent.com/en/product/gas-chromatography/gc-systems/8890-gc-system</a>	LINK

Equipment	
<b>HP5-MS</b>	NAME
GC Column	TYPE
Agilent	BRAND
19091S-433UI	SKU
<a href="https://www.agilent.com/store/productDetail.jsp?catalogId=19091S-433UI">https://www.agilent.com/store/productDetail.jsp?catalogId=19091S-433UI</a>	LINK

Equipment	
<b>Inlet Liner</b>	NAME
Agilent	BRAND
5183-4711	SKU
<a href="https://www.agilent.com/store/productDetail.jsp?catalogId=5183-4711">https://www.agilent.com/store/productDetail.jsp?catalogId=5183-4711</a>	LINK
Inlet liner, split, single taper, glass wool, deactivated	SPECIFICATIONS

## SAFETY WARNINGS

-  All chemicals used for this procedure are hazardous. Read the Safety Data Sheet (SDS) for all chemicals and follow all applicable chemical handling and waste disposal procedures. Manufacturer specific SDS information can be found by following the CAS numbers of compounds in 'Materials' list.

## BEFORE START INSTRUCTIONS

All solvents, analytes, and chemicals used in this protocol are listed in the 'Materials' section. They are excluded from in-line referencing to keep steps clear and concise.

## Internal Standard Preparation

- 1 This analysis uses 1,3,5-tri-tert-butylbenzene (TTB) as an internal standard. Prepare adequate volume to allow for addition of 1  $\mu\text{L}$  of TTB per 100  $\mu\text{L}$  of sample or standard volume.
- 2 Create a 10,000  $\mu\text{g/mL}$  internal standard working solution by weight of TTB using methanol as the diluent.

## Preparation of Standards

- 3 By weight, create individual 20,000  $\mu\text{g/mL}$  stock solutions of all monomers (listed in Materials section) and use methanol as the diluent.

### Note

If compounds have not previously been analyzed on the column specified in this method, it is necessary to analyze each monomer separately (at a concentration of approximately 100  $\mu\text{g/mL}$ ) to determine retention time for each before analyzing a mixture containing all compounds.

- 4 Combine the stock solutions to create a 1000  $\mu\text{g/mL}$  mixed standard working solution in methanol.
  - 4.1 For example, to prepare a 10 mL mixed standard working solution of woody RCF analytes (distinguished in Materials section), add 500  $\mu\text{L}$  of each of the 20,000  $\mu\text{g/mL}$  stock solutions (14 analytes) and add 3000  $\mu\text{L}$  methanol.
- 5 Using the mixed standard working solution at 1000  $\mu\text{g/mL}$ , create a calibration curve from 10  $\mu\text{g/mL}$  to 1000  $\mu\text{g/mL}$  with a minimum of five calibration points using methanol as the diluent.

- 6** Add 1  $\mu\text{L}$  of 10,000  $\mu\text{g/mL}$  TTB internal standard working solution per every 100  $\mu\text{L}$  of calibration standard.  
Example: 10  $\mu\text{L}$  into 1000  $\mu\text{L}$  or 5  $\mu\text{L}$  into 500  $\mu\text{L}$

## 6.1

Calibration Level	Concentration ( $\mu\text{g/mL}$ )	Volume of Mixed Standard Working Solution ( $\mu\text{L}$ )	Volume of Methanol ( $\mu\text{L}$ )	Volume of Internal Standard ( $\mu\text{L}$ )
9	1000	1000	0	10
8	750	750	250	10
7	500	500	500	10
6	250	250	750	10
5	100	100	900	10
4	75	75	925	10
3	50	50	950	10
2	25	25	975	10
1	10	10	990	10

Example Calibration

## Sample Preparation

- 7 Oil Samples:**
1. Weigh approximately 5-10 mg of oil into a GC vial and record weight.
  2. Add a known volume of methanol to each vial containing sample (for example 1 mL).
  3. Add 1  $\mu\text{L}$  of 10,000  $\mu\text{g/mL}$  TTB for every 100  $\mu\text{L}$  of methanol added to the sample vial.
- Example: 10  $\mu\text{L}$  into 1000  $\mu\text{L}$  or 5  $\mu\text{L}$  into 500  $\mu\text{L}$

### Liquid Samples:

1. Aliquot a known volume of sample into a GC vial (for example 1 mL)
  2. Add 1  $\mu\text{L}$  of 10,000  $\mu\text{g/mL}$  TTB for every 100  $\mu\text{L}$  of methanol added to the sample vial.
- Example: 10  $\mu\text{L}$  into 1000  $\mu\text{L}$  or 5  $\mu\text{L}$  into 500  $\mu\text{L}$

## GC-FID Analysis

- 8** Analyze samples using an 8890 Agilent Gas Chromatograph (GC) or equivalent equipped with a flame ionization detector (FID) per the method parameters below:

## 8.1 Method Parameters

### *Split/Splitless Inlet Parameters:*

Inlet Temperature: 280 °C

Injection Volume: 1 µL

Split Ratio: 2:1

Inlet Liner: split, single taper, glass wool, deactivated (see Materials)

Syringe: P/N 5181-8809

Wash Solvent: Methanol

### *Column:*

HP5-MS (see Materials)

Carrier Gas: Helium

Flow Rate: 1 mL/min (constant flow)

### *Oven Parameters:*

Maximum Oven Temperature: 280 °C

	Rate °C/min	Value °C	Hold Time min	Run Time min
Initial		70	2	2
Ramp	10	280	2	25

### *Detector Parameters:*

Detector: FID

Detector Temperature: 300 °C

Air Flow: 400 mL/min

H2 Fuel Flow: 40 mL/min

Makeup Flow: 10 mL/min

## Analytical Quality Control

- 9 Several strategies are utilized when performing this analysis to ensure instrument stability and reproducibility.

## 9.1 Calibration Curves

All compounds must have a correlation coefficient ( $r^2$ ) of 0.995 or greater using a linear calibration fit.

## 9.2 Calibration Verification Standards (CVS)

A calibration verification standard (CVS) is a standard from the calibration curve that is re-analyzed every 20 or fewer samples to ensure instrument drift remains within the determined acceptance criteria. Acceptable CVS recoveries for this analysis are within 15% of the expected amount. Acceptance criteria may differ between instruments and should be determined experimentally.