



VERSION 3

NOV 14, 2022



IN DEVELOPMENT

 Goga Lab RT-qPCR protocol: QuantStudio6 Machine V.3**This protocol is published without a DOI.**Jeremy.williams<sup>1</sup><sup>1</sup>UCSF

Jeremy.williams

COMMENTS 0

## ABSTRACT

Guidelines for preparing RT-qPCR samples for QuantStudio 6 located in HSW7 lab space.

## ATTACHMENTS

[GogaLab-RTqPCR-Preparation.xlsx](#)

## PROTOCOL CITATION

Jeremy.williams 2022. Goga Lab RT-qPCR protocol: QuantStudio6 Machine. **protocols.io**  
<https://protocols.io/view/goga-lab-rt-qpcr-protocol-quantstudio6-machine-ci95uh86>  
Version created by Jeremy.williams



## LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 14, 2022

## LAST MODIFIED

Nov 14, 2022

## GUIDELINES

- \*start with full tip boxes and use tips in coordination with your plate map, so you never get lost
- \*watch 1uL volumes in the pipet tip like a hawk - major source of variability
- \*cap and gently vortex your \*mixture per primer set\* every 12 replicates or so, sometimes things separate

## MATERIALS TEXT

PowerUp SYBR Green Master Mix  
Your generated cDNA samples  
100uM single primer stocks

## BEFORE STARTING

\*Thaw on ice (leave time!) and keep all reagents on ice through preparation. Prepare the plate on ice.

## ATTACHMENTS

[GogaLab-RTqPCR-Preparation.xlsx](#)

### Prior to plate preparation:

- 1 Dilute stock IDT primers to 100uM (see note in Guidelines).
- 2 Use 'GogaLab-RTqPCR-Preparation' Excel spreadsheet to input your sample number and calculate reagent volumes.

### Prepare your plate:

- 3 Mix \*total reagent volumes required\* (see spreadsheet, green) for DI and PowerUP
- 4 Mix your forward and reverse primer pairs together, to a final dilution of 10uM forward and 10uM reverse. For example, add 10uL each of forward and reverse primers to 80uL PCR-quality DI for 100uL final volume.

- 5 Mix \*mixtures per primer set\* (blue) volumes together.
- 6 Pipet 1uL diluted cDNA into respective qPCR well, aiming for the sidewall of each respective well.
- 7 Add 19uL \*mixture per primer set\* into each respective well.
- 8 Seal plate, and spin down 1000rpm for 1 minute. Use bacterial, not tissue culture, centrifuge.
- 9 Load plate into centrifuge and proceed using QuantStudio software suite.