



JAN 23, 2024

🌐 DNA damage assessment in the adult Drosophila brain via pH2Av (S137) immunostaining

Mel Feany^{1,2}

¹Department of Pathology, Brigham and Women's Hospital, Harvard Medical School, Boston, Massachusetts 02115;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815

ASAP Collaborative Research Network

OPEN  ACCESS



Beatrice Weykopf

ABSTRACT

This protocol describes how to stain adult Drosophila brains with pH2Av to assess DNA damage.

DOI:

dx.doi.org/10.17504/protocols.io.dm6gp36q5vzp/v1



Protocol Citation: Mel Feany 2024. DNA damage assessment in the adult Drosophila brain via pH2Av (S137) immunostaining. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp36q5vzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 17, 2024

Fly fixation and paraffin embedding




- 1 Fix flies of the desired genotype  Overnight in 4% formalin, embedded heads in paraffin, cut 4-micron sections, and dry at 42°C  Overnight

Deparaffinization

- 2 Deparaffinize and bring through ethanol to water (xylene x 2, 100% ethanol x 2, 95% ethanol x 2, H2O x 2)


Antigen-retrieval

35m

- 3 Microwave slides in 10 mM sodium citrate for  00:15:00 . Cool for  00:20:00 . 35m
Stock: 100 mM sodium citrate, pH 6.0. Use at least  1 L of citrate solution in a large glass box to avoid drying

Blocking

1h 30m

- 4 Block in blocking solution (PBS with 0.3% Triton and 2% dry milk) for 30min to  01:00:00 1h

Antibody incubation

- 5 Incubate with primary antibody (Rockland pH2AV, cat # 600-401-914, dilution 1:100) overnight at room temperature

- 6 Wash 3 x in PBST (PBS with 0.3% Triton)
- 7 Incubate with fluorophore-coupled secondary antibody (Thermo Fisher) diluted 1:200 in blocking solution at room temperature
- 8 Wash 2 x in PBST and 1 x in PBS

Mounting

- 9 Mount in medium containing DAPI (Southern Biotech)