

May 23, 2022

🌐 Microglia FACS staining after isolation (from UCSD)

[rabdelha](#)¹

¹California Institute of Technology

1



dx.doi.org/10.17504/protocols.io.dm6gpbxp8lzp/v1

 [rabdelha](#)

This protocol details about microglia FACS staining after isolation (from UCSD).

[437-925.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.dm6gpbxp8lzp/v1

[rabdelha](#) 2022. Microglia FACS staining after isolation (from UCSD). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpbxp8lzp/v1>



Microglia FACS staining, antibodies (CD11b, CD45)

_____ protocol ,

May 17, 2022

May 23, 2022

May 17,
2022

 [renuka.s](#)

May 23,
2022

 [rabdelha](#)

62713

Gating Notes

- live vs dead
- FSC vs SSC
- singles vs. doublets.
- CD11B (high) vs. CD45 (low).
- CX3CR1 (optional) ◇ gate on positive cells.



Misc. Notes

- **don't need to compensate with only 2-3 colors.
- CD45, CD11B, and CX3CR1 are all surface markers.
- Usually get 20-50% live cells after staining.
- Better to do staining right after microglia prep rather than waiting O/N for staining.

Processing Notes

Can do RNA seq or Atac Seq

RNA Seq

- Spin down sample in Eppendorf (not in facs tubes)
- Remove supernatant
- Add  **150 µL** of trizol → resuspend
- Store in  **-80 °C** (can keep stored for months).




ATAC Seq

- Do transposase reaction and then freeze for processing.

Materials


use 5ml Polypropylene eppendorf tubes

Antibodies

-  [CD11b Monoclonal Antibody \(M1/70\)](#)
[APC eBioscience/Invitrogen Catalog #17-0112-82](#)
-  [Alexa Fluor® 488 anti-mouse CD45](#)
[Antibody BioLegend Catalog #103122](#)
-  [CD16/CD32 Monoclonal Antibody](#)
[\(93\) eBioscience/Invitrogen Catalog #14-0161-82](#)
- DAPI 1:10,000, for 45 sec, wash twice with PBS; dilute in H2O

Protocol

25m 45s

- 1 Resuspend cells staining buffer ( **300 µL** of HBSS+EDTA+BSA).
 - Resuspend in a 15-ml falcon tube.

2



15m

Add Fc Block (~1:100 dilution) and incubate for ⌚ 00:15:00 at 🌡 4 °C or in the fridge.

- Add this to cells resuspended in the HBSS.

3

Take ~5% of cells for unstained control, put 🌡 On ice .

Already did this during the isolation during the percoll separation step.

4



Add antibodies (CD11b, CD45, 1:100 dilution) for 20-30 min at 🌡 4 °C .

4.1



Add directly to cells in the FC block.

4.2



Mix by tapping tube; do not vortex.

Don't add to unstained cells!!

5



45s

Add DAPI, ⌚ 00:00:45 , then dilute with HBSS+BSA+EDTA.

5.1 Add 1:1000 dilution.

5.2 Add to unstained cells.

- 6 Put 70-µm filter on top of a FACS test tube and filter resuspended pellet into the tube.

Do this for unstained samples as well.


- 6.1 Push filter hard onto tube.



10m

Centrifuge for  **400 x g, 00:10:00** .

- 7.1 Prepare collection tubes during spin time.

- Coat FACS tubes with  **1 mL** of staining buffer to prevent cells from sticking to walls of tubes.
- Invert tube several times.

- 8 Vacuum out supernatant.
- Remove most of supernatant; not all.

- 9 Resuspend pellet in  **500 µL** to  **1000 µL** staining buffer.

Depending on machine.

Keep cells on ice the whole time.