



SOLUTION- 04 - Wash solution (RPMI/FBS) for PBMC V.2

Marco Cosentino¹, Elisa Storelli¹, Alessandra Luini¹, Massimiliano LM Legnaro¹, Emanuela Rasini¹, Marco Ferrari¹, Franca Marino¹

¹Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

1 Works for me dx.doi.org/10.17504/protocols.io.biaikace



Elisa Storelli
Center for Research in Medical Pharmacology, University of I...

DOI

dx.doi.org/10.17504/protocols.io.biaikace

DOCUMENT CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 04 - Wash solution (RPMI/FBS) for PBMC. **protocols.io**
dx.doi.org/10.17504/protocols.io.biaikace

LICENSE

————— This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 07, 2020

LAST MODIFIED

Jul 23, 2020

DOCUMENT INTEGER ID

38954

Sterile RPMI/FBS 10% (v/v)

Prepare with sterile RPMI and FBS (Fetal Bovine Serum). Both can be found in Fridge 2 - (Room PS03).

Prepare 50 ml of RPMI/FBS by adding 5 ml of FBS to 45 ml of RPMI, work under a laminar flow hood (Room PS03).

Storage: **4 °C** Fridge 2 - (Room PS03).

RPMI: catalog number ECM0495L, Euroclone, Italy.

FBS: catalog number ECS0180L, Euroclone, Italy.