



Version 3

Jul 29, 2021

# Simple electroporation for efficient CRISPR/Cas9 genome editing in murine zygotes V.3

Simon E. Tröder<sup>1</sup>, Branko Zevnik<sup>1</sup><sup>1</sup>CECAD, University of Cologne, Cologne, Germany

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bwzhpf36](https://doi.org/10.17504/protocols.io.bwzhpf36)

Simon Tröder

## ABSTRACT

Electroporation of zygotes represents a rapid alternative to the elaborate pronuclear injection procedure for CRISPR/Cas9-mediated genome editing in mice. However, current protocols for electroporation either require the investment in specialized electroporators or corrosive pre-treatment of zygotes which compromises embryo viability. Here, we describe an easily adaptable approach for the introduction of specific mutations in C57BL/6N mice by electroporation of intact zygotes using a common electroporator with synthetic CRISPR/Cas9 components and minimal technical requirement.

## DOI

[dx.doi.org/10.17504/protocols.io.bwzhpf36](https://doi.org/10.17504/protocols.io.bwzhpf36)

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0196891>

## PROTOCOL CITATION

Simon E. Tröder, Branko Zevnik 2021. Simple electroporation for efficient CRISPR/Cas9 genome editing in murine zygotes. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bwzhpf36>

Version created by Simon Tröder

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Tröder SE, Ebert LK, Butt L, Assenmacher S, Schermer B, Zevnik B (2018) An optimized electroporation approach for efficient CRISPR/Cas9 genome editing in murine zygotes. PLoS ONE 13(5): e0196891. doi: 10.1371/journal.pone.0196891

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 29, 2021

## LAST MODIFIED

Jul 29, 2021

## PROTOCOL INTEGER ID

51977

## GUIDELINES

Use only embryo-grade reagents.

M2 and M16 media are prepared as published (Behringer, R., et al., (2014) Manipulating the mouse embryo : a laboratory manual, Fourth edition. ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York). Commercial M2 (Sigma, #M7167) and KSOM/AA (Merck Millipore, #MR-106-D) can be used instead. Mouse pre-implantation embryos are incubated in at least 4h pre-equilibrated M16 or KSOM/AA in a CO<sub>2</sub> incubator (5% CO<sub>2</sub>, 37°C, 95% humidity).

Embryo-grade T<sub>10</sub>E<sub>0.1</sub> (10 mM Tris-HCl, 0.1 mM EDTA, pH 7.4) buffer is prepared as described (Chu, V.T., et al., (2016) Efficient generation of Rosa26 knock-in mice using CRISPR/Cas9 in C57BL/6 zygotes. BMC Biotechnol 16, 4.).

Mouse zygotes are collected from the oviducts of superovulated females as described in published protocols (Behringer, R., et al., (2014) Manipulating the mouse embryo : a laboratory manual, Fourth edition. ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York).

If you use this protocol please cite our publication:

Tröder SE, Ebert LK, Butt L, Assenmacher S, Schermer B, Zevnik B (2018) An optimized electroporation approach for efficient CRISPR/Cas9 genome editing in murine zygotes. PLoS ONE 13(5): e0196891. doi: 10.1371/journal.pone.0196891

## MATERIALS TEXT

### MATERIALS

 [Cas9 Nuclease Integrated DNA](#)

**Technologies Catalog #1074181**

 [tracrRNA Integrated DNA](#)

**Technologies Catalog #1072532**

 [crRNA Integrated DNA Technologies](#)

 [ssODN \(Ultramer DNA Oligonucleotides\) Integrated DNA Technologies](#)

 [Opti-MEM \(Reduced Serum Medium\) Thermo Fisher](#)

**Scientific Catalog #31985062**

 [M2 Medium Contributed by users](#)

 [M16 Medium Contributed by users](#)

 [T10E0.1 Buffer Contributed by users](#)

 [Electroporation Cuvette 1mm BioRad](#)

**Sciences Catalog #1652089**

 [60 mm Center Well Organ Culture](#)

**Dish Corning Catalog #353037**

### guide RNA annealing

- 1 Resuspend lyophilized crRNA, tracrRNA and ssODN in T<sub>10</sub>E<sub>0.1</sub> buffer to 100 µM (e.g. 5 nmol in 50 µl) (Store at -80 °C until use)
- 2 Combine 5 µl crRNA (100 µM) and 5 µl tracrRNA (100 µM) in a nuclease-free PCR tube to yield an equimolar crRNA:tracrRNA duplex solution of 50 µM

- Heat to 95 °C for 5 min and cool down at 5 °C/ min in a thermocycler (crRNA:tracrRNA duplex can be stored for weeks at -80 °C)

#### Preparation of the electroporation mix

- Add 1.6 µl crRNA:tracrRNA duplex and 1.3 µl Cas9 nuclease to 15.1 µl Opti-MEM in a nuclease-free tube and vortex
- Incubate mix at room temperature for 10 min  
🕒 00:10:00
- Place tube on ice, add 2 µl ssODN and vortex
- Quick-spin at 4 °C and keep tube on ice until use
- Summary of the 20 µl electroporation mix:

Reagent	Stock concentration	Final concentration	Volume
crRNA:tracrRNA duplex	50 µM	4 µM	1.6 µl
Cas9 nuclease	61 µM (10 µg/µl)	4 µM	1.3 µl
ssODN	100 µM	10 µM	2.0 µl
Opti-MEM	-	-	15.1 µl

#### Electroporation of zygotes

- Collect zygotes from the oviducts of superovulated females as described in published protocols
- Wash the zygotes in five drops of M2
- Wash up to 50 zygotes in one drop of Opti-MEM
- Transfer zygotes with as little media as possible to the 20 µl electroporation mix
- Using a 20 µl pipette transfer the entire drop including the zygotes into a pre-warmed (37 °C) 1 mm electroporation cuvette  
(Ensure retrieving all zygotes by quickly aspirating the entire drop. Slow aspiration will leave zygotes behind)
- Insert cuvette into a standard electroporator (e.g. BioRad Gene Pulser Xcell electroporator)

- 15 Apply two square wave pulses at 30 V and 3 ms duration with a 100 ms interval
- 16 Retrieve the zygotes by flushing the cuvette with 100 µl pre-incubated M16 using a 100 µl pipette into a culture dish (e.g. 60 mm Center Well Organ Culture Dish)
- 17 Wash the cuvette with 100 µl pre-incubated M16
- 18 Transfer all zygotes to a new culture dish with 500 µl pre-incubated M16
- 19 Incubate zygotes in M16 until the two cell stage and transfer the developed embryos into pseudopregnant foster mice (Embryos may also be transferred at the one cell stage)