



Nov 25, 2020

OG1RF lipid extraction protocol

Elizabeth Fozo¹¹In-house protocol

1

Works for me

This protocol is published without a DOI.



Eadewunm

ABSTRACT

Lipid Extraction for *E. faecalis* OG1RF

PROTOCOL CITATION

Elizabeth Fozo 2020. OG1RF lipid extraction protocol. [protocols.io](https://protocols.io/view/og1rf-lipid-extraction-protocol-bp4smqwe)
<https://protocols.io/view/og1rf-lipid-extraction-protocol-bp4smqwe>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 25, 2020

LAST MODIFIED

Nov 25, 2020

PROTOCOL INTEGER ID

44914

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Lipid Extraction for *E. faecalis* OG1RF

Steps

- 1 Inoculate 50mL of BHI at O.D. of 0.01 from overnight culture. Spike supplement at OD600 0.25. Incubate for 30 minutes and harvest cells. Spin at 2,739 G for 10 minutes.
- 2 Wash twice in 25mL 1X PBS. Spin for 5 minutes at 3,000 rpm.

- 3 Re-suspend pellets in 1mL of PBS
- 4 Add 100uL of stock lysozyme (100mg/mL – stock made in water) pre-warmed to 37oC. Incubate cells for 20 minutes at 37oC in 50 mL conical.

Meanwhile, add 0.5g glass beads to the screw cap microcentrifuge tube.
- 5 Add the lysozyme treated cells to screw cap microcentrifuge tube containing 0.5g of glass beads (<100um). Bead beat for two intervals of 1 minute each. Leave on ice at least 30 seconds between cycles.
- 6 Transfer mixture to 15 mL conical tube (vortex the bead beating tube then pour contents into 15mL conical). Add 1.5mL methanol and 3mL chloroform. (You can add this to the tubes ahead of time.)
- 7 Then, gently vortex mixture. Centrifuge for 5 minutes.
There will be two layers, separated by a white interface, and glass beads at the bottom of the tube.
- 8 Using a glass Pasteur pipet – transfer the liquid above the white interface and the liquid below the white interface (leaving behind the beads at the bottom and as much of the interface as possible) to a new 15mL conical.
- 9 Add 1.5mL of 0.9% NaCl to the new 15mL conical.
- 10 Gently vortex mixture and centrifuge for 5 minutes.
- 11 Remove the top layer. Pour interface and bottom layer to glass screw lid tube.
- 12 Store at -20 until delivery to Mass Spec facility.