



Version 2

May 18, 2020

Acropora DNA extraction with Qiagen DNeasy tissue kit V.2

Iliana Baums¹, Sheila Kitchen¹¹Pennsylvania State University

1

Works for me

This protocol may be deleted by the owner

Astrangia

Coral Genotyping Protocols



Iliana Baums

ABSTRACT

DNA extraction protocol for Acropora or other coral tissue based on Qiagen DNeasy kit. This extraction protocol works well for the Acropora SNPchip and other coral genotyping applications (such as microsatellite genotyping). It preferentially extracts coral host DNA but some *Symbiodiniaceae* DNA will be present.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Baums IB, Hughes CR, Hellberg MH (2005) Mendelian microsatellite loci for the Caribbean coral *Acropora palmata*. Mar Ecol Prog Ser 288:115-127. doi:10.3354/meps288115.

ATTACHMENTS

[Acroporid tissue DNA extraction.docx](#)

GUIDELINES

This protocol works well on coral tissue preserved in 95% non-denatured Ethanol. Do not overload the columns. A few polyps is enough. Make sure your centrifuge can spin high enough to complete this protocol before you start. Both the 96-well plate version and the single 1.5 ml tube version of the QIAGEN kit work. For the 1.5ml tube version, we like the electronic multi-channel MATRIX pipettor with flexible spacing to go from 1.5 ml tubes to PCR tubes/plates. Avoid defrost cycles of samples and extracted DNA.

MATERIALS

| NAME | CATALOG # | VENDOR |
|--|-----------|----------------------|
| QIAGEN DNeasy Blood and Tissue Kit, 50 rxn | 69504 | Qiagen |
| Proteinase K, 2mL | 19131 | Qiagen |
| Incubator | | |
| Razor blade | | |
| Centrifuge | 5415 D | Eppendorf Centrifuge |

Sample preparation

- 1 Use sterile technique to cut 3-4 polyps (skeleton with tissue) off the coral fragment with a razor blade.
- 2 Transfer the pieces to a 1.5 ml tube/96-well plate.
- 3 Add 180 ul of Buffer ATL and 20 ul of proteinase K (from Qiagen kit) to each tube/well.
- 4 If using the tubes, invert the tubes several times to gently mix the reagents. *Do not vortex*

- 5 Incubate the mixture overnight at 56 °C.
- 6 a. The mixture will have a slight yellow hue after lysis from the algal symbionts.
- 7 b. The tubes or plate can be placed on a rocking platform for gentle agitation.

Qiagen DNeasy Blood and Tissue Kit protocol

- 8 Follow Qiagen's protocol with the following changes:
Following step 6 of the spin-column protocol, we perform the optional centrifugation for 1 min at 20,000 xg to dry the column membrane.

Elution 2x

- 9 Elution 1- Add 100 ul of low TE buffer directly to the column membrane, incubate for 1 min and centrifuge as recommended.
a. This elution typically contains smaller fragmented DNA.
- 10 Elution 2 – Move the column to a new 1.5 ml tube. Repeat elution as above, but increase incubation of 100 ul of low TE to 5 min.
a. This elution contains high molecular weight DNA and is used for downstream analyses.
b. Store at -20 °C.

Check extraction quality

- 11 Run DNA on a 1.5% agarose gel to check extraction quality of elution 1 and 2. Look for one high molecular weight band at the top of the gel with little smearing.