



Oct 06, 2020

DNA Extraction from FANS sorted nuclei

Stefania S Policicchio¹, Jonathan P Davies¹, Barry Chioza¹, Aaron Jeffries¹, Joe Burrage¹, Jonathan Mill¹, Emma L Dempster¹

¹University of Exeter Medical School, Exeter, UK

1 Works for me dx.doi.org/10.17504/protocols.io.bmpmk5k6

Complex Disease Epigenetics Group



Stefania Policicchio

ABSTRACT

Here we describe an optimised protocol for the extraction of genomic DNA from frozen nuclei samples collected from different neural cells using Fluorescence-assisted Nuclei Sorting (FANS) of post mortem human brain. The same protocol is also suitable for DNA extraction from cells allowing optimal recovery in terms of DNA purity and yield.

DOI

dx.doi.org/10.17504/protocols.io.bmpmk5k6

PROTOCOL CITATION

Stefania S Policicchio, Jonathan P Davies, Barry Chioza, Aaron Jeffries, Joe Burrage, Jonathan Mill, Emma L Dempster 2020. DNA Extraction from FANS sorted nuclei. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bmpmk5k6>

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 24, 2020

LAST MODIFIED

Oct 06, 2020

PROTOCOL INTEGER ID

42445

MATERIALS TEXT

Reagents required

- **Slagboom buffer (50mL):**
 - 10x STE buffer (5mL)
 - 5% SDS (5mL)
 - Water (40 mL)
- **RNAase A** (stock concentration: 10mg/ml)
- **Proteinase K** (stock concentration: 20 mg/mL)
- **Majiik mix** (1:1 ratio yeast Reagent 3 (Autogen Bioclear, Caine, Wiltshire, UK) + 100% ethanol)
- **100% Isopropanol**
- **80% Ethanol**
- **TE or water**



NOTE- Slagboom buffer can be made in advance while 80% Ethanol should be prepared fresh the same day.

ABSTRACT

Here we describe an optimised protocol for the extraction of genomic DNA from frozen nuclei samples collected from different neural cells using Fluorescence-assisted Nuclei Sorting (FANS) of post mortem human brain. The same protocol is also suitable for DNA extraction from cells allowing optimal recovery in terms of DNA purity and yield.

DNA extraction from FANS purified nuclei

- 1 Defrost frozen FANS sorted nuclei contained in a 1.5ml eppendorf tube **On ice**













NOTE 1-2 nuclei aliquots (~200,000 nuclei/tube) per each nuclei population should yield between 300-500ng DNA).

- 2 To each sample add either 250µL / 500µL / 1 mL of Slagboom buffer (SB) depending on size of cell pellet. Normally 250µL is sufficient for 200,000 nuclei aliquot or 1×10^6 cells
- 3 Add 500µL of SB to each nuclei sample (if nuclei pellet)*



* If nuclei are stored in running buffer: calculate the volume of 10x STE and 5% SDS to add to each sample to be consistent with the composition of SB (1x STE; 0.5% SDS).
e.g. If you have collected a 400µL sample you would add 50µL 10x STE and 50µL 5% SDS to make the sample up to 500µL








- 4 Add 1µL of RNase A per 500uL buffer
- 5 Incubate at **37 °C** for **00:45:00** (using a heat block)
- 6 Add proteinase K to a final concentration 2mg/mL (e.g 5µL for every 500µL of SB)
- 7 Mix by inverting 10 times (do not vortex or pipette mix)

- 8 Incubate at  **60 °C** for  **01:00:00** (water bath)
- 9 Move samples from the water bath and leave at  **Room temperature** for  **00:05:00**
- 10 Briefly spin down tubes (pulse spin)
- 11 For every 500µL of SB used, add 100 µL of Majiik Mix (e.g. 200 µL for 1 mL of SB)
- 12 Mix by vigorous inversions (**DO NOT** vortex)
- 13 Centrifuge at  **17000 x g** for  **00:10:00** at  **Room temperature**
- 14 Carefully recover supernatant and transfer it to a new labelled 1.5ml tube (leaving ~50µL at the bottom of the tube)
- 15 Repeat step above by adding another 100µL of Majiik Mix to the new tube
- 16 Mix by vigorous inversions (**DO NOT** vortex)
- 17 Centrifuge at  **17000 x g** for  **00:10:00**  **Room temperature**
- 18 Recover supernatant and transfer to a new labelled 1.5ml tube leaving ~50ul at the bottom of the old tube.



NOTE - If supernatant exceeds 1mL, split into 2 tubes

- 19 Add an equal volume of 100% Isopropanol to each tube (e.g. 1mL supernatant + 1 mL 100% Isopropanol)

- 20 After adding the Isopropanol, add 0.7-1 µL Glycogen Blue per tube (this step is optional – but advised as it helps to visualize DNA pellet at the end)
- 21 Mix by inversion (10 times)
- 22 Centrifuge at  **17000 x g** for  **00:15:00** at  **Room temperature** (hinges of tubes facing upwards)
- 23 Carefully remove supernatant and discard
-  **NOTE** - Care must be taken when pouring supernatant since the DNA pellet only weakly adheres to the side of the tube
- 24 Proceed by adding 500µL of 80% ethanol to each tube
- 25 Mix gently by pipetting
- 26 Centrifuge at  **17000 x g** for  **00:05:00** at  **Room temperature**
- 27 Carefully discard supernatant
- 28 Pulse centrifuge to collect remaining liquid at the bottom of the tube
- 29 Recover and discard residual liquid from the bottom of the tube making sure **not to** disturb the DNA pellet
- 30 Leave DNA pellets to air dry for 5-15 minutes (leave lids open)

 **OPTIONAL** - Leave to dry for additional  **00:10:00** at  **37 °C** but be careful not to over-dry pellets

- 31 Resuspend pellet in 20uL of ddH₂O or 1xTE (15μL per tube if expected very low yield). Avoid pipette mixing, only gently flicking
- 32 Pulse spin tubes to help pellet resuspend/ dislodge from the tube wall
- 33 Leave tubes at **4 °C** overnight to fully resuspend before quantifying DNA samples (Nanodrop or Qubit measurements).