



OCT 23, 2023

🌐 JMN-MSMP Notexin Muscle scRNAseq

ccherry¹

¹JMN-MSMP

Cellular Senescence Network (SenNet) Method Development Community

SenNet JMN-MSM



ccherry

ABSTRACT

Murine notexin muscle injury followed by single cell RNA sequencing

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.e6nvwd8r9lmk/v1

Protocol Citation: ccherry
2023. JMN-MSMP Notexin
Muscle scRNAseq.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.e6nvwd8r9lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Oct 23, 2023

Last Modified: Oct 23, 2023

PROTOCOL integer ID:
89762

Keywords: SenNet JMN-MSMP

Skeletal muscle injury and collection

- 1 In a given mouse, both tibialis anterior (TA) muscles were injected with 10 µl of notexin (10 µg/ml; Latoxan, France).
- 2 At various days post injury, the mouse was sacrificed and the TAs were collected. Each TA was processed independently to generate single cell suspensions.

Generate single-cell suspensions

- 3 Muscles were enzymatically digested with 8 mg/ml Collagenase D (Roche, Basel, Switzerland) and 10 U/ml Dispase II (Roche, Basel, Switzerland) and then manually dissociated to generate cell suspensions.
- 4 Myofiber debris was removed by filtering the cell suspensions through a 100 µm and then a 40 µm filter (Corning Cellgro # 431752 and # 431750).
- 5 After filtration, erythrocytes were removed by incubating the cell suspension in erythrocyte lysis buffer (IBI Scientific # 89135-030).
- 6 After digestion, the single-cell suspensions were washed and resuspended in 0.04% BSA in PBS at a concentration of 10^6 cells/ml. A hemocytometer was used to manually count the cells to determine the concentration of the suspension.

Single-cell RNA library preparation and sequencing

- 7 Single-cell RNA-sequencing libraries were prepared using the Chromium Single Cell 3' reagent kit v3 (10x Genomics, Pleasanton, CA) following the manufacturer's protocol.
- 8 Cells were diluted into the Chromium Single Cell A Chip to yield a recovery of 6,000 single-cell transcriptomes with <5% doublet rate.
- 9 Libraries were sequenced on the NextSeq 500 (Illumina, San Diego, CA).

