



APR 05, 2024

Brain Homogenization and MSD Protocol for Mouse Brain and Serum

Scott Vermilyea¹

¹University of Minnesota

ASAP Collaborative Research Network

Team Lee



Jane Balster
ASAP - Team Lee

OPEN ACCESS



ABSTRACT

This protocol details about the homogenization and MSD protocol of the mouse brain using various panel kits.

ATTACHMENTS

[320-681.pdf](#)

DOI:

dx.doi.org/10.17504/protocols.io.kxygxprqkl8j/v1

Protocol Citation: Scott Vermilyea 2024. Brain Homogenization and MSD Protocol for Mouse Brain and Serum. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxprqkl8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 12, 2021

Last Modified: Apr 05, 2024

PROTOCOL integer ID: 55080


Keywords: Homogenization, MSD Protocol, Mouse Brain and Serum, ASAPCRN

MATERIALS


Reagent, Kits and Equipment List


1.


Equipment	
gentleMACS™ C Tubes	NAME
C Tubes	TYPE
gentleMACS™	BRAND
130-093-237	SKU
https://www.miltenyibiotec.com/US-en/products/gentlemacs-c-tubes.html#grep	LINK


2.  RIPA Lysis and Extraction Buffer **Thermo Fisher Catalog #89901**
3. Halt PIC (ThermoFisher 1861279)
4.

Equipment	
Direct Detect® Assay-free Cards	NAME
Assay-free Cards	TYPE
Direct Detect®	BRAND
DDAC00010-GR	SKU
https://www.emdmillipore.com/US/en/product/Direct-Detect-Assay-free-Cards,MM_NF-DDAC00010-GR?bd=1	LINK
Membrane cards used for protein quantitation in the Direct Detect® Infrared Spectrometer	SPECIFICATIONS

- 5. HCl (Millipore Sigma H1758-500ML)
- 6. NaOH (Millipore Sigma 415413-500ML)
- 7.  HEPES **Sigma Aldrich Catalog #83264-500ML-F**
- 8.

 Proinflammatory Panel 1 mouse **MESO SCALE DIAGNOSTICS, LLC. Catalog #MSD K15048D**

- 9.
-  Cytokine Panel 1 Mouse Kit **MESO SCALE DIAGNOSTICS, LLC. Catalog #K15245D**

- 10.
-  U-PLEX TGF- β Combo mouse **MESO SCALE DIAGNOSTICS, LLC. Catalog #MSD K15242K-2**

11.

Equipment	
Octo Dissociator with Heaters	NAME
Octo Dissociator	TYPE
gentleMACS™	BRAND
130-096-427	SKU
https://www.miltenyibiotec.com/US-en/products/gentlemacs-octo-dissociator-with-heaters.html	LIN K

12.

Equipment	
Direct Detect Spectrometer	NAME
Spectrometer	TYPE
Direct Detect®	BRAND
DDHW00010-WW	SKU
https://www.emdmillipore.com/US/en/product/Direct-Detect-Spectrometer,MM_NF-DDHW00010-WW?bd=1	LINK

13. Meso Sector S 600 (MSD 1201)

Brain Homogenization Using Miltenyi gentleMACS Octodissociator

43m



- 1
- Weight out brain using pre-chilled C-tube.
- 2
- Use

🧪 2 mL

 RIPA (1x Halt PIC) per whole brain,

🧪 1 mL

 RIPA (1x Halt PIC) per half brain.
- 3
- Homogenize brain using “protein_01_01” routine on Miltenyi gentle MACS Octodissociator.
- 3.1
- Repeat homogenization routine if not fully homogenized.

4 Vortex briefly and incubate  00:30:00  On ice .


30m




5 Centrifuge at  3500 rpm, 4°C, 00:03:00 .

3m




6 Transfer homogenate to two  1.5 mL eppendorf tubes.



7 Vortex briefly and centrifuge at  14000 x g, 4°C, 00:10:00 .

10m




8 Transfer supernatant to new  1.5 mL eppendorf tubes







9 Perform 1:10 serial dilution.










9.1 Dilute using  10 µL sample to  90 µL RIPA (1x Halt PIC) to achieve 1:100 dilution. (1/2)

9.2 Dilute using  10 µL sample to  90 µL RIPA (1x Halt PIC) to achieve 1:100 dilution. (2/2)

- 10 Blot  2 μ L diluted sample onto Direct Detect assay card (run duplicates).
- 11 Use  2 μ L RIPA (1x Halt PIC) as blank.
- 12 Read assay card on Direct Detect Spectrometer.
- 13 Detectable range is  0.2 mg/mL to  5.0 mg/mL , dilute or concentrate as needed.
- 14 Average all readings for each sample to determine concentration.

Brain MSD Protocol Proinflammatory Panel 1 mouse and Cytokine Panel 1...

- 15 Using concentration determined, calculate volume needed for  200 μ g protein.
- 16 Use RIPA(1x PIC) to dilute  200 μ g protein to  100 μ L ( 2 undetermined).
- 
- 17 Combine  100 μ L sample and  100 μ L MSD diluent in U-bottom plate.

18

Load  50 μ L diluted sample to MSD plate per duplicate.


19

Proceed with MSD assay protocol from manufacturer.



Brain MSD Protocol U-PLEX TGF- β Combo mouse

10m

20

Using concentration determined, calculate volume needed for  200 μ g protein.





21

Use RIPA(1x PIC) to dilute  200 μ g protein to  100 μ L ( 2 undetermined).



22



Load  100 μ L sample into U-bottom plate.

23


Add  20 μ L  1 Molarity (M) HCl into U-bottom plate and shake for  00:10:00 at  25 $^{\circ}$ C . 10m

24

Neutralize sample with  14 μ L  1.2 Molarity (M) NaOH in  0.5 Molarity (M) HEPES.

25 Combine  100 μ L treated sample and  100 μ L MSD diluent in new U-bottom plate.




26 Load  50 μ L diluted sample to MSD plate per duplicate.




27 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol Proinflammatory Panel 1 mouse

28 Add  25 μ L MSD Diluent to MSD Plate per duplicate.




29 Add  25 μ L serum sample to MSD Plate per duplicate.



30 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol Cytokine Panel 1 mouse

31 Add  37.5 μ L MSD Diluent to MSD Plate per duplicate.



32 Add 12.5 μ L serum sample to MSD Plate per duplicate.



33 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol U-PLEX TGF- β Combo mouse

10m

34 Add 50 μ L serum sample to U-bottom plate.



35 Add 10 μ L 1 Molarity (M) HCl into U-bottom plate and shake for 00:10:00 at 25 $^{\circ}$ C. 10m



36 Neutralize sample with 7 μ L 1.2 Molarity (M) NaOH in 0.5 Molarity (M) HEPES.



37 Add 25 μ L MSD Diluent to MSD Plate per duplicate.



38 Add 25 μ L treated serum sample to MSD Plate per duplicate.



39 Proceed with MSD assay protocol from manufacturer.

