



Apr 09, 2020

# Adapter ligation with AMII

**Josh Quick**<sup>1</sup><sup>1</sup>University of Birmingham**1** Works for me [dx.doi.org/10.17504/protocols.io.bdp9i5r6](https://dx.doi.org/10.17504/protocols.io.bdp9i5r6)**ARTIC** **Coronavirus Method Development Community** [view 1 more group](#) **Josh Quick**

## ABSTRACT

This is a subprotocol for performing adapter ligation with AMII

## EXTERNAL LINK

<http://lab.loman.net/protocols/>

## ATTACHMENTS

[One-pot native barcoding protocol \(1\).pdf](#)

## SAFETY WARNINGS

See SDS (Safety Data Sheet) for safety warnings and hazards.


- 1 Set up the following AMII adapter ligation reaction:

Component	Volume
End-repaired amplicon pools	30 µl
NEBNext Quick Ligation Reaction Buffer (5X)	10 µl
Adapter Mix (AMII)	5 µl
Quick T4 DNA Ligase	5 µl
<b>Total</b>	50 µl






There will be some variation in clean-up efficiencies but expect to carry around 80% through a clean-up.

- 2 Incubate at room temperature for 00:20:00

- 3 Add  **50 µl** (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting.



Vortex SPRI beads thoroughly before use to ensure they are well resuspended, the solution should be a homogenous brown colour.

- 4 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 5 Incubate for  **00:05:00** at room temperature.
- 6 Place on magnetic rack and incubate for  **00:02:00** or until the beads have pelleted and the supernatant is completely clear.
- 7 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 8 Add  **250 µl** SFB and resuspend beads completely by pipette mixing.





SFB will remove excess adapter without damaging the adapter-protein complexes. Do not use 70% ethanol as in early clean-ups.

- 9 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 10 Remove supernatant and discard.
- 11 Repeat steps 14-16 to perform a second SFB wash.
- 12 Pulse centrifuge and remove any residual SFB.



You do not need to allow to air dry with SFB washes.

- 13 Add  **15 µl** EB and resuspend beads by pipette mixing.
- 14 Incubate at room temperature for  **00:02:00** .
- 15 Place on magnetic rack.
- 16 Transfer final library to a new 1.5mL Eppendorf tube.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits

