

Sep 26, 2021

## Cyanobacteria Yolk-Shell Preparation

Jorge Fernández Méndez<sup>1</sup>, Celia Martin Morales<sup>1</sup>

<sup>1</sup>iGEM MADRID\_UCM 2021

1 Works for me Share

dx.doi.org/10.17504/protocols.io.byirpud6

Jorge Fernández Méndez

## ABSTRACT

Cyanobacteria as well as other photosynthetic microorganisms can be encapsulated within biohybrid nanomaterials. Cyanobacterial Yolk-shell are biohybrid nanostructures with a size of 1 to 6 µm, where the cell is encapsulated within a non-contacting silica nanocapsule. Yolk-shell encapsulation provides long-term cell viability, higher resistance against harsh environments and superior photosynthetic activiy. In addition, engineering the colloidal packing allows tunable shell-pore diameter for size-dependent permeability and introduction of new functionalities for specific molecular recognition.

Yolk-Shell structures can be easily synthetized via self-assebly of colloidal suspensions of silica nanoparticles as described in this protocol

DOI

dx.doi.org/10.17504/protocols.io.byirpud6

PROTOCOL CITATION

Jorge Fernández Méndez, Celia Martin Morales 2021. Cyanobacteria Yolk-Shell Preparation. **protocols.io** https://dx.doi.org/10.17504/protocols.io.byirpud6

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

L. Wang et al., "Single-cell yolk-shell nanoencapsulation for long-term viability with size-dependent permeability and molecular recognition", National Science Review, vol. 8, no. 4, 2020. Available: 10.1093/nsr/nwaa097 [Accessed 26 September 2021].

**KEYWORDS** 

Cyanobacteria, Encapsulation, Immobilization, Yolk-Shell, Nanomaterials, Biohybrids, Photosynthesis, 4C\_Fuels, iGEM

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 26, 2021

LAST MODIFIED

Sep 26, 2021

PROTOCOL INTEGER ID

53553

Citation: Jorge FernÃÅ¡ndez MÃÅ®ndez, Celia Martin Morales (09/26/2021). Cyanobacteria Yolk-Shell Preparation. <a href="https://dx.doi.org/10.17504/protocols.io.byirpud6">https://dx.doi.org/10.17504/protocols.io.byirpud6</a>

• Fresh Cyanobacterial Culture

**⊗**LUDOX® TMA colloidal silica **Sigma** 

- Aldrich Catalog #420859
- Protamine Sulfate (Biochemistry Grade)
- PBS x1 Buffer pH = 7.4
- BG-11 Media (HEPES 10 mM pH = 8)

## Cells Pretreatment

10m

1 @2800 rpm, 00:10:00

In sterility, centrifugue a sufficient ammount of exponentially growing cyanobacteria culture (20 to 100 mL). OD<sub>720</sub> should be between 0.5 to 1. It is recommended to split the culture in two different sterile centrifugue tubes.

- In sterility, discard the supernatant and add 1/10<sup>th</sup> of the initial volume with sterile PBS buffer.

  Resuspend the cell pellet gently, gently flicking the tube or slowly pipetting up and down with a wide open pipette tip.
- 3 Take 1000 μl sample from the resuspended cells and measure the Optical Density at 720 nm using a Spectrophotometer. OD<sub>720</sub> should be 1.3. If OD<sub>720</sub> is greater than 1.3, ogo to step #2 and add more sterile PBS till reaching the required OD<sub>720</sub> value.
- 4 In Sterility, add 1/50 part (v/v) of the available resuspended cells volume of 5% protamine sulfate in sterile PBS. It is recommended to flick the tube up and down a couple of times and incubate the cells with protamine sulfate for © 00:05:00. (No agitation required)

10m

(3) 2800 rpm, 00:10:00 Centrifugue the cells and discard the supernatant.

**IMPORTANT!** Write the final volume of cell suspension in the tube before discarding the supernatant.

After centrifugation, resuspend the cell pellet in the same volume of cell suspension formerly used. Final OD<sub>720</sub> should still be around 1.3

Yolk-Shell Self Assembly

5

22m

Add the **necessary ammount of Ludox(R) TMA** colloidal suspension to reach a final concentration of 1 mg/mL in the cell suspension.

$$v_{Ludox-TMA} = VLudox(\mu L) = 418/VcellSuspension(mL) * 1000$$

**NOTE.** Ludox TMA has a  $\rho$  = 1.23 g/mL and it's a colloidal suspension at 34%  $_{Wt}$ , then each mL of the stock will carry 418 mg of silica nanoparticles. To know the required ammount of Ludox, just divide your cell suspension volume by 418.

 8 Incubate the suspension for  $\bigcirc$  **00:20:00** with mild agitation using a rocker shaker at minimum speed.

Storing Yolk-Shell 10m

- 9 **3600 rpm, 00:04:00** After incubation, centrifugue the former suspension and **discard the supernatant**. Wash the yolk-shell pellet with fresh **BG-11 media** using the same ammount of volume used during incubation (Other media could also be used for marine cyanobacteria strains).
- Repeat the wash step  $\circ$  go to step #9, and resuspend the cells with fresh BG-11 media to a final OD<sub>720</sub> = 1.

In this conditions yolk shell can be stored at room temperature and mild illumination for long periods of time.