

Sep 17, 2020

# COVID Airway Processing for scRNAseq

Peter Szabo<sup>1</sup>, Steven Wells<sup>1</sup>, Peter A. Sims<sup>1</sup>, Donna Farber<sup>1</sup><sup>1</sup>Columbia University

1

Works for me

[dx.doi.org/10.17504/protocols.io.bjj8kkw](https://dx.doi.org/10.17504/protocols.io.bjj8kkw)

Peter Sims

## ABSTRACT

This protocol describes the isolation of lymphocytes and pan-mononuclear cells from airway samples.

## ATTACHMENTS

[Columbia\\_Wells\\_COVID\\_Airway\\_Ficoll\\_Processing\\_for\\_scRNAseq.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.bjj8kkw](https://dx.doi.org/10.17504/protocols.io.bjj8kkw)

## PROTOCOL CITATION

Peter Szabo, Steven Wells, Peter A. Sims, Donna Farber 2020. COVID Airway Processing for scRNAseq.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjj8kkw>



## KEYWORDS

SARS-CoV-2, COVID-19, lymphocytes, isolation, pan-mononuclear cells

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 11, 2020

## LAST MODIFIED

Sep 17, 2020

## OWNERSHIP HISTORY

Aug 11, 2020  Julia Rossmanith protocols.io

Sep 02, 2020  Peter Sims

## PROTOCOL INTEGER ID

40288

## MATERIALS


NAME	CATALOG #	VENDOR
Gibco™ DPBS no calcium no magnesium	14190144	Thermo Fisher Scientific
Penicillin-Streptomycin-Glutamine (100X)	10378016	Thermo Fisher
UltraPure™ 0.5 M EDTA pH 8.0	15575020	Thermo Fisher Scientific
Thermo Scientific™ Nunc™ 50mL Conical Sterile Polypropylene Centrifuge Tubes	12-565-271	Fisher Scientific

NAME	CATALOG #	VENDOR
5mL Falcon™ Round-Bottom Polypropylene Test Tubes	14-959-11A	Fisher Scientific
Biotin anti-human CD235ab Antibody	306618	BioLegend
Biotin anti-human CD66b Antibody	305120	BioLegend
BioMag® Plus Streptavidin	BP628	Bangs Laboratories
Corning™ Externally Threaded Cryogenic Vials	09-761-71	Fisher Scientific
CryoStor CS10 100ML	NC9930384	Fisher Scientific
Gibco™ Fetal Bovine Serum qualified Australia	10-099-141	Fisher Scientific
Ficoll-Paque™ PLUS Media	45-001-749	Fisher Scientific
Human TruStain FcX™	422302	BioLegend
Gibco™ IMDM (Iscoves Modified Dulbeccos Medium)	12-440-053	Fisher Scientific
NC-Slide A8™ box with 25 Slides	942-0003	Chemometec
Solution 13 AO – DAPI	910-3013	Chemometec
Dead Cell Removal Kit	130-090-101	Miltenyi Biotec
MS Columns	130-042-201	Miltenyi Biotec

#### MATERIALS TEXT

##### Equipment


- Centrifuge
- Cell Counter - NC-3000
- EasyEights™ EasySep™ Magnet (Stemcell Technologies, Cat. No.: 18103)



EasyEights™ EasySep™ Magnet  
Magnet for column-free immunomagnetic separation

EasySep      18103      [↗](#)

- MACS Multistand (Miltenyi, Cat. No.: 130-108-934)



MACS MultiStand  
Separator for magnetic cell separation

MACS      130-108-934      [↗](#)

#### EQUIPMENT

NAME	CATALOG #	VENDOR
EasyEights™ EasySep™ Magnet	18103	Stemcell Technologies
MACS MultiStand	130-108-934	Miltenyi Biotec

#### SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

##### Biosafety Notes

- All materials required for sample processing are to be prepared in the biosafety cabinet before handling of Airway samples.
- All Airway sample manipulation takes place in a biosafety cabinet unless specifically stated.
- All centrifugation steps must take place in capped containers. Upon completion of a centrifugation step, return

the capped container to the biosafety cabinet, remove all tubes, and spray and wipe them with >70% ethanol before continuing to the next step.

### Preparing Buffer and Media

1 

Create the following **DPBS Solution-EDTA** in a bottle of DPBS by using the table below:

Component	Volume (mL)	Starting Conc.	Final Conc.
DPBS	474	-	
FBS	25	100%	5%
EDTA	1	0.5 M	1 mM

Table 1.

2 

Create **IMDM Layering Media** in a bottle of IMDM.



Component	Volume (mL)	Starting Conc.	Final Conc.
IMDM	470	-	
FBS	25	100%	5%
EDTA	5	100%	1%

Table 2.

### Preparation of Sample


- 3 Place sample box into biosafety cabinet.
- 4 Remove samples from box and from the containment bags, discard bags, spray sample containers with >70% ethanol and wipe down.



**NOTE:** If sputum traps come with tubing attached flush tubing with  **1 mL** -  **2 mL** (or more) before proceeding to optimize yield.

- 5 Record the total volume of Airway sample to be processed in mL.

6 

Transfer the Airway Sample to a 50 mL tube, add  **4 mL DPBS** to the airway sample.

7  

Spin the Airway for  **400 x g, 20°C, 00:10:00**, remove and save four aliquots of the supernatant in cryovials. **Add the same volume to all cryovials.**

Record number of vials: \_\_\_\_\_ and the volume per vial \_\_\_\_\_ mL.



**NOTE:** Ensure that cryovials are decontaminated prior to removal from the biosafety cabinet.

8



Add **IMDM Layering Media** to the Airway sample to bring the *total volume* to **25 mL**.

9



Add **40 µl Benzonase** to the sample, pipette up and down vigorously to mix and dissociate airway pellet, incubate at **Room temperature** for **00:30:00**.

10



Add **750 µl 0.5M EDTA**, and filter sample through a **100 µM** filter.

#### Ficoll-Paque

11



Aliquot **15 mL Ficoll-Paque Media PLUS** to a 50 mL tube.

12



Using the **slow** setting on the pipette gun, **gently** layer the airway/IMDM mixture on top of the

**15 mL Ficoll-Paque Media PLUS**.

**Take extra care not to disturb the interface while layering. Disturbing the interface excessively prevents the mononuclear cells from becoming a clean layer.**

13



Spin for **1200 x g, 20°C, 00:20:00** with **no brake**, 4 acceleration.



**NOTE:** Centrifuge should be pre-warmed to **20 °C**.

14



Remove the mononuclear cell layer from each tube and transfer to a new 50 mL tube. **Take extra care to avoid pulling cells from the ficoll layer (underneath the mononuclear cell layer) as this typically contains a lot of granulocytes.** Pulling from the plasma layer is not an issue.



15 

Top MNC with **cold** DPBS Solution-EDTA to **40 mL** (ensure at least 2-3 volumes are added) and centrifuge the cell suspension(s) for **400 x g, 4°C, 00:10:00**.

16 

(Platelet Spin) Discard the supernatant, top tube to **40 mL** with **cold** DPBS Solution-EDTA, and centrifuge the cell suspension for **120 x g, 4°C, 00:10:00**.

17 

Remove the supernatant (**caution: pellet may be loose**), and re-suspend the cell pellet in **4.5 mL Dulbecco's Phosphate Buffered Saline (DPBS)** (final volume should be about **5 mL**).

#### Cell Counting of COVID Samples

18 

Add **0.05 mL sample**, **0.05 mL DPBS**, and **0.005 mL Solution 13** to a 1.5 mL centrifuge tube, incubate for **00:02:00** at **Room temperature**.

19 

Add **0.1 mL BD Cytofix Fixation Buffer** to the samples and incubate **00:30:00**, **Room temperature**, and protect from light.

20 

Aliquot **0.01 mL sample** to the well of a NC-Slide A8 and count on the NC-3000.

Record number and viability below, calculate total cells:

cell number: \_\_\_\_\_ cells/mL, \_\_\_\_\_ % viable

#### Division of Sample for scRNAseq Analysis and Freeze-down

21 

Aliquot up to  $2 \times 10^7$  cells to a 5 mL Falcon Round-Bottom tube and place **On ice** for subsequent sample clean-up (next section).

22 Freeze down up to  $1 \times 10^8$  cells in approximately  $1 \times 10^7$  aliquots (**1 mL** each) using Cryostor CS10 Medium, a Mr. Frosty, and a **-80 °C** freezer.

Record the number of vials frozen: \_\_\_\_\_ and the cells per cryovial frozen: \_\_\_\_\_.

## Sample Clean Up for scRNAseq – CD66b and CD235ab removal

23



Centrifuge the single cell suspension for **400 x g, 4°C, 00:05:00**.

24



Discard the supernatant and resuspend the cell pellet in **50 µl DPBS Solution-EDTA**.

25



Add **10 µl Human TruStain FcX** to the single cell suspension and incubate for **00:10:00**, **4 °C**.

26



Add **10 µl biotinylated anti-CD66b and biotinylated anti-CD235ab** to the sample and incubate for **00:30:00**, **4 °C**.

27



While the single cell suspension is incubating add **0.2 mL BioMag Plus Streptavidin Beads** to a 5 mL Falcon Round-Bottom tube.

28



Add **2 mL DPBS Solution-EDTA** to the BioMag Plus Streptavidin Beads and place on a magnet for **00:05:00**.

29



Remove all the supernatant from the BioMag Plus Streptavidin Beads, remove from the magnet and resuspend the beads in **0.1 mL DPBS Solution-EDTA**.

30





Once step 26 ( **go to step #26** ) is complete, add **3 mL DPBS Solution-EDTA** to the single cell suspension and centrifuge for **400 x g, 4°C, 00:05:00**.

31



Resuspend the single cell suspension in the BioMag Plus Streptavidin Beads from step 29 ( **go to step #29** ), and incubate at **Room temperature** for **00:05:00**.

32 

Add  **3 mL DPBS** to the tube and place on a magnet for  **00:05:00**.

33 

Remove supernatant from tube and transfer to a separate 5 mL Falcon Round Bottom tube.

#### Sample Clean Up for scRNAseq – Dead Cell Removal


34  

Centrifuge the single cell suspension for  **400 x g, 4°C, 00:05:00**, and discard supernatant.

35  

Resuspend cell pellet in  **0.1 mL Dead Cell Removal Microbeads**, mix well, incubate,  **Room temperature**,  **00:15:00**.

36 

While the cell suspension is incubating, place an MS Column onto the MACS Multistand and rinse with  **0.5 mL 1x Binding Buffer Solution**.

37 Post incubation, apply cell suspension to the MS Column and capture the flow through in a 5 mL Falcon Round Bottom tube.

38 

Rinse with  **1.5 mL 1x Binding Buffer** and capture in the same tube.

39 






Centrifuge the single cell suspension for  **400 x g, 4°C, 00:05:00**, and discard supernatant.

40  




Resuspend cell pellet in  **0.5 mL DPBS**, and count cells.

#### Cell Counting of COVID Samples (10x)


41  

Add  **0.05 mL sample**,  **0.05 mL DPBS**, and  **0.005 mL Solution 13** to a 1.5 mL centrifuge tube, incubate for  **00:02:00** at  **Room temperature**.

42  

Add  **0.1 mL BD Cytofix Fixation Buffer** to the samples and incubate  **00:30:00** ,  **Room temperature** , and protect from light.

43 

Aliquot  **0.01 mL sample** to the well of a NC-Slide A8 and count on the NC-3000.

44 

Record number and viability below, calculate total cells:

cell number: \_\_\_\_\_ cells/mL, \_\_\_\_\_ % viable

#### 10X Encapsulation

45 Follow the appropriate 10X protocol (Chromium Next GEM Single Cell 3' Reagent Kits v3.1 User Guide – Rev D) for encapsulation of cells from the airway sample.