



MAY 13, 2023

NCEM Drop - Tissue Dounce Homogenisation (TM-014)

sandra.crameri¹

¹CSIRO ACDP



sandra.crameri

ABSTRACT

tissue dounce homogenisation

OPEN  ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.eq2ly7odw1x9/v1

Protocol Citation: sandra.crameri 2023. NCEM Drop - Tissue Dounce Homogenisation (TM-014).

protocols.io

<https://dx.doi.org/10.17504/protocols.io.eq2ly7odw1x9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 13, 2023

Last Modified: May 13, 2023

PROTOCOL integer ID:
81823




HEADER

1 SAN:

SPEC No:





OPERATOR:

Dounce Homogenisation

- 2 This method is often used for scabs, warts, or chunks of tissue such as liver. Place ~ 5mm x 5mm piece of scab or tissue into glass homogeniser tube.
- 3 Insert teflon plunger into cordless drill, set drill on max speed.
- 4 Add equal volume TC water to tube, homogenise(20 strokes) inside Class II BSC cabinet.
- 5 Clarify resultant solution at  13000 rpm in Eppendorf centrifuge for  00:01:00 1m
- 6 Stand for  00:05:00 to permit viruses to diffuse back into solution, from the debris, and provide an interface for the sampling of membrane associated viruses. Use supernatant for sample below. 5m

Conventional

11m

- 7 Adsorb  10 μ L sample to grid  00:10:00 , inspect to ensure sample does not dry out. 10m
- 8 Drain excess sample from grid using filter paper, leave wet.
- 9 Stain  nano-W Contributed by users Catalog ##2018-5ML  00:01:00 1m
- 10 Drain & dry using filter paper