



Oct 18, 2020

Colony PCR

Jiaxin Li¹¹South China University of Technology

1

Works for me

dx.doi.org/10.17504/protocols.io.bnihmcb6

Jiaxin Li

DOI

dx.doi.org/10.17504/protocols.io.bnihmcb6

PROTOCOL CITATION

Jiaxin Li 2020. Colony PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bnihmcb6>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 16, 2020

LAST MODIFIED

Oct 18, 2020

PROTOCOL INTEGER ID

43305

- 1 Draw grid on clean agar plate (use the appropriate antibiotic). You will use this to inoculate an overnight culture for a positive clone.
- 2 Select colonies to pick. Streak portion of colony to numbered sector and place the remainder in a correspondingly numbered PCR tube with 20 μ L of autoclaved water.
- 3 Heat at 95 °C for 10 minutes.



This can be done in PCR machine

- 4 Use 1 μ L of reaction for colony PCR (T5 Taq Polymerase)

5

Component	Volume/ μ L	Final Concentration

2XT5 Super PCR Mix	5ul	1x
10uM Forward Primer	0.4ul	0.4uM
10uM Reverse Primer	0.4ul	0.4uM
Template DNA	1ul	<1ug
ddH2O	Add to 10ul	∅
Total	10ul	∅

6

Temperature	time
98°C	03:00
98°C	00:10
Tm-5°C	00:10
72°C	20s/kb
72°C	02:00
4°C	∞