



Oct 24, 2020

# Measuring the amount of bacteria in a soil sample

Andreea S<sup>1</sup><sup>1</sup>University of Groningen

Other dx.doi.org/10.17504/protocols.io.bkrckv2w

iGEM Groningen 2020

r.rajadhyaksha

## ABSTRACT

Colony Forming Units (c.f.u) is a unit that is used in microbiology to estimate the number of viable bacteria or fungal cells in a sample. It also depends on their ability to multiply under controlled conditions. In the paper published by *El-Hassan and Gowen, 2006*, they analyzed various formulations of *Bacillus subtilis* by counting the CFU of *B. subtilis* present in every formulated product.



S. A. El-Hassan and S. R. Gowen (2006). Formulation and Delivery of the Bacterial Antagonist *Bacillus subtilis* for Management of Lentil Vascular Wilt Caused by *Fusarium oxysporum* f. sp. *lentis*. *Journal of Phytopathology*, Volume 154, Issue 3.

<https://doi.org/10.1111/j.1439-0434.2006.01075.x>

## DOI

[dx.doi.org/10.17504/protocols.io.bkrckv2w](https://dx.doi.org/10.17504/protocols.io.bkrckv2w)

## PROTOCOL CITATION

Andreea S 2020. Measuring the amount of bacteria in a soil sample. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bkrckv2w>

## LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 03, 2020

## LAST MODIFIED

Oct 24, 2020

## PROTOCOL INTEGER ID

41476

## ABSTRACT

Colony Forming Units (c.f.u) is a unit that is used in microbiology to estimate the number of viable bacteria or fungal cells in a sample. It also depends on their ability to multiply under controlled conditions. In the paper published by *El-Hassan and Gowen, 2006*, they analyzed various formulations of *Bacillus subtilis* by counting the CFU of *B. subtilis* present in every formulated product.



S. A. El-Hassan and S. R. Gowen (2006). Formulation and Delivery of the Bacterial Antagonist *Bacillus subtilis* for Management of Lentil Vascular Wilt Caused by *Fusarium oxysporum* f. sp. *lentis*. *Journal of Phytopathology*, Volume 154, Issue 3.  
<https://doi.org/10.1111/j.1439-0434.2006.01075.x>

#### CFU determination

- 1 Colony Forming Units (CFU) can be determined by estimating the OD of spore suspension using a tube-reading spectrophotometer adjusted at 1.978 [corresponding to  $8.5 \cdot 10^{10}$  CFU/ml] at 600nm absorbance wavelength
- 2 The formulation will be placed on sterile aluminum foil in pans and air-dried for **24:00:00** with occasional stirring in a laminar airflow cabinet.
- 3 Dried formulations (35% moisture content) of *B. mycoides* will be passed through a 250µm mesh sieve to attain the desired particle size.
- 4 Pack in sterilized polypropylene bags, seal and store at room temperature prior to use.
- 5 Count CFUs to estimate the number of viable propagules of *B. mycoides* using the standard dilution platin method described in step 6.

#### Standard dilution method

- 6 Take three **1 g** aliquots of the dried powder and place in **99 mL** sterile PBST solution (this will include PBS + **0.05 % (v/v)** Tween 20). Stir magnetically at high speed for **00:15:00** . Now dilute this suspension with approximately and take **0.2 mL** of this suspension and plate on Nutrient Agar (NA) media.