



Jan 18, 2021

PCR for SARS-CoV-2

Leonardo Caserta¹¹Cornell University

1

Works for me

dx.doi.org/10.17504/protocols.io.brk4m4yw

Diel Lab



Leonardo Caserta

DOI

dx.doi.org/10.17504/protocols.io.brk4m4yw

DOCUMENT CITATION

Leonardo Caserta 2021. PCR for SARS-CoV-2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.brk4m4yw>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 18, 2021

LAST MODIFIED

Jan 18, 2021

DOCUMENT INTEGER ID

46460

Multiplex PCR

Set up two multiplex PCR reaction mix as follows:

Q5 High-Fidelity 2X Master Mix	12.5 µl
H2O	6.75 µl
Primers pool (1 or 2) at 10 µM	1.75 µl
cDNA	4 µl
TOTAL	25 µl

Cycling conditions for **Q5 High-Fidelity 2X Master Mix (~1500bp)**

A	B	C
Temp.	Time	Number of cycles
98	30 sec	1
98	10 sec	4
65	2 min	
72	3 min	
98	10 sec	3
64	2 min	
72	3 min	
98	10 sec	2
63	2 min	
72	3 min	
98	10 sec	2
62	2 min	
72	3 min	
98	10 sec	2
61	2 min	
72	3 min	
98	10 sec	22
60	2 min	
72	3 min	
72	3 min	1
4	∞	