



Jul 06, 2021

T4PNK end-healing of RNA samples

Jonathan Howard¹¹University of California, Santa Cruz

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bwdspa6e

Jonathan Howard

ABSTRACT

Resolution of 5' and 3' ends of sample RNA prior to library preparation.

DOI

dx.doi.org/10.17504/protocols.io.bwdspa6e

PROTOCOL CITATION

Jonathan Howard 2021. T4PNK end-healing of RNA samples. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bwdspa6e>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 06, 2021

LAST MODIFIED

Jul 06, 2021

PROTOCOL INTEGER ID

51346

GUIDELINES

Do not let pellet from ethanol precipitation and wash over dry

Visually inspect top phase after phase lock tube spin to make sure there is no secondary phase at bottom of RNA phase. If there is, carefully remove top phase and leave bottom phase. This is the organic phase that has leaked into the RNA phase and should be avoided.

MATERIALS TEXT

5x PNK pH 6.5 buffer

350 mM Tris-HCl, pH 6.5(Calbiochem; CAS 77-86-1)

50 mM MgCl₂(Sigma; M8266-100G)

5 mM dithiothreitol(Sigma; D9779-5G)

Acid-Phenol:Chloroform pH 4.5(ThermoFisher; AM9720)

GlycoBlue Coprecipitant(ThermoFisher; AM9515)

Phase Lock gel heavy tubes, 2 mL(VWR; 10847-802)

T4 Polynucleotide Kinase (NEB; M0201S)

PNK treatment and cleanup

- 1 Set up the following reaction components in a sterile PCR tube:
 - 4 μL 5 \times PNK pH 6.5 buffer
 - 0.5 μL PNK
 - 0.5 μL RNasin
 - 5 μL Plus-AlkB/Minus-AlkB RNA sample
 - To 20 μL with Nuclease-free Water
- 2 Put reaction at 37°C for 30 mins. This can be done in a thermocycler
- 3 Put samples on ice, add to the reaction:
 - 5 μL of 10 \times PNK reaction buffer (NEB-supplied buffer)
 - 5 μL of 10 mM ATP
 - 1 μL of T4PNK (additional)
 - 1 μL of RNase inhibitor (optional)
 - Up to 50 μL with Nuclease-free Water
- 4 Put reaction at 37°C for 30 min.
- 5 Raise volume of T4PNK-treated RNA samples to 400 μL .
- 6 Collect the solution and add it together with 400 μL acid phenol/chloroform to a 2 ml Phase Lock Gel Heavy tube. Mix by inversion in hand for 30 sec.
- 7 Separate the phases by spinning for 5 min at full speed and room temperature.
- 8 Transfer the aqueous layer into a new tube (be careful not to touch the gel matrix with the pipette or remove organic phase that may have not migrated below the gel matrix).
- 9 Precipitate by adding 0.75 μL glycoblue and 40 μL 3 M sodium acetate pH 5.5. Mix, then add 1 ml 100% ethanol, mix again, and place at -20°C overnight.
- 10 Centrifuge at 15,000 rpm at 4°C for 20 min. Remove the supernatant and wash the pellet with 500 μL 80% ethanol and vortex for 30 sec. Spin again for 5 min. to pellet RNA again.

- 11 Remove ethanol wash completely (be careful to not remove pellet). Allow to air-dry at room temperature for 10 mins. Remove any excess ethanol in tube with additional pipetting.
- 12 Resuspend the pellet in 5 μ L of nuclease-free H₂O and transfer to new 0.2 mL RNase-free PCR tube.