




Apr 15, 2021

Metal-antibody MIBItag conjugation kit

Marc MB Bosse¹, Sean Bendall¹, Mike Angelo¹¹Department of Pathology, Stanford University**1** Works for me dx.doi.org/10.17504/protocols.io.bhyej7teHuman BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu Marc Bosse
Stanford University

ABSTRACT

This is a "protocol IO" extended version of MIBItag conjugation kit protocol from lonpath.

DOI

dx.doi.org/10.17504/protocols.io.bhyej7te

PROTOCOL CITATION

Marc MB Bosse, Sean Bendall, Mike Angelo 2021. Metal-antibody MIBItag conjugation kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bhyej7te>

KEYWORDS

null, lonpath, MIBItag, conjugation, metal, antibody

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 26, 2020

LAST MODIFIED


Apr 15, 2021

PROTOCOL INTEGER ID

38630

MATERIALS TEXT

MATERIALS

 Centrifugal Filter Unit: 50 kDa Amicon Ultra 500 µL V

[bottom Millipore Catalog #UFC505096](#) In 2 steps

Worksheet antibody information

1

A	B	C	D	E	F	G	H	I	J
Antibody (clone)	Isotope/Element	Pre-label Nanodrop (µg/mL)	Volume to label 200 µg	Post-label Nanodrop (µg/mL)	Labelled volume	Dilution for 200 µg/mL	Total volume	Volume of stabilizer	Conjug ID
ALK									

[☒ MIBItag conjugation](#)

[kit IonPath Catalog #600XXX](#)

Optional antibody concentration

10m

10m

2 

Critical: the 50-kDa MWCO micro-filter device is designed to hold up to 500 µL

If 100 µg of antibody requires addition of more than 200 µL in volume

Add up to  **500 µl** of antibody solution in the 50-kDa MWCO micro-filter device

Centrifuge  **12000 rcf, Room temperature**  **00:10:00** or  **21000 rcf, Room temperature** , **00:05:00**

[☒ Centrifugal Filter Unit: 50 kDa Amicon Ultra 500 µL V](#)


[bottom Millipore Catalog #UFC505096](#)

2.1 

Make sure to align the filters parallel to centrifuge force to ensure optimal tangential flow filtration of both filters

This can be done by centering the hinge of the collecting tube between the two filters

The 50-kDa MWCO micro-filter device is then placed in the centrifuge with the hinge outward




Antibody buffer exchange




- 3 Add 300 μ L of Buffer 2 (100 mM phosphate buffer, 2.5 mM EDTA, pH 7.2) to a 50-kDa MWCO micro-filter device

 [Centrifugal Filter Unit: 50 kDa Amicon Ultra 500 \$\mu\$ L V](#)

[bottom Millipore Catalog #UFC505096](#)

- 4 If step 2 not performed and antibody solution volume is less than  200 μ l

Add the 100 μ g antibody solution volume

- 5  **12000 rcf, Room temperature**  **00:10:00** or  **21000 rcf, Room temperature , 00:05:00**

- 6 Add 400 μ L of buffer 2 IonPath


- 7 

 **12000 x g, Room temperature**  **00:10:00** or  **21000 rcf, Room temperature , 00:05:00**

Make sure that the residual dead volume is below the bottom end of the filter (~20 μ L)

Do not measure with a pipette, just inspect visually if the level of the fluid is as shown in the pictures

(blue dye was used for demonstration only to help visualizing the fluid)



- 8 Dilute **8 µl** of TCEP (0.5M) in **992 µl** of buffer 2 (**4 Millimolar (mM)** final)

Tris(2-carboxyethyl)phosphine hydrochloride solution Sigma

Aldrich Catalog #646547-10X1ML

Note: **10 µl** TCEP aliquots are stored at **-20 °C**

Antibody reduction

- 9 Add **100 µl** of diluted TCEP to the concentrated antibody in the 50-kDa MWCO micro-filter device
- 10 Mix antibody and TCEP solution by flicking the 50-kDa MWCO micro-filter device gently

- 11

Incubate at **37 °C** for **00:30:00 ONLY**

Do not incubate longer than 30 min, this will result in full reduction of disulfide bonds necessary for the structural integrity of the antibody

Antibody wash

- 12 After incubation, add immediately **300 µl** of buffer 3 to the partially reduced antibody in the 50-kDa MWCO micro-filter device
- 13 Mix reduced antibody and wash buffer 3 (Tris buffered salt, 1 mM EDTA, pH 7.5) by flicking the 50-kDa MWCO micro-filter device gently
- 14 **12000 rcf, Room temperature , 00:10:00** or **21000 rcf, Room temperature , 00:05:00**

- 15 Discard the flow through
- 16 Add additional **400 µl** of buffer 3 [go to step #14 once](#)

Conjugation of metal and antibody


- 17 Reconstitute preloaded polymer in **100 µl** of buffer 3
- 18 Transfer the preloaded polymer solution in the reduced antibody/ 50-kDa MWCO micro-filter device
- 19 Mix reduced antibody and preloaded polymer by flicking the 50-kDa MWCO micro-filter device gently
- 20 Incubate **37 °C** for **01:30:00** to **02:00:00**

Final wash

- 21 Add 250 µL of buffer 4 to the antibody conjugation mixture in the the 50-kDa MWCO micro-filter device
- 22 Mix antibody conjugation mixture and buffer 4 (Tris buffed salt, pH 7.5) by flicking the 50-kDa MWCO micro-filter device gently
- 23 **12000 rcf, Room temperature , 00:10:00** or **21000 rcf, Room temperature , 00:05:00**
- 24 Discard the flow through
- 25 Add **400 µl** of buffer 4 in the the 50-kDa MWCO micro-filter device [go to step #22](#) to 24
- 26 Add **400 µl** of buffer 4 in the the 50-kDa MWCO micro-filter device [go to step #22](#) to 24
- 27 Add **400 µl** of buffer 4 in the the 50-kDa MWCO micro-filter device [go to step #22](#) to 24

28 Add  **400 µl** of buffer 4 in the the 50-kDa MWCO micro-filter device [↻ go to step #22](#) to 24


Labelled antibody retrieval

29 Add  **50 µl** of buffer 4 to the washed labelled antibody in the the 50-kDa MWCO micro-filter device

30 Mix washed labelled antibody and buffer 4 by flicking the 50-kDa MWCO micro-filter device gently

31 Invert the micro-filter into a new collection tube

32  **10000 rcf, Room temperature , 00:01:00**


33 Add  **80 µl** of buffer 4 to the washed labelled antibody in the the 50-kDa MWCO micro-filter device
[↻ go to step #30](#) to 33

34 Measure the volume retrieved with a pipette (200 µL) and register data in labelled volume on the worksheet
[↻ go to step #1](#)

35 Open Nanodrop using Protein 280 and chose type = IgG

35.1 Clean the cell with ultrapure (type 1) water

35.2 Blank with buffer 4

35.3 Take  **1.5 µl** of labelled antibody and measure concentration

35.4 Register data in labelled volume on the worksheet [↻ go to step #1](#)


36 Open **MIBitracker** and click on **CONJUGATES** tab

36.1 Click on **Add Conjugate** tab



36.2 Select the antibody reference used for the conjugation

Fill the following fields: **Isotope**; **Channel**; **Element**; **µg for labelling**; **Responsible**



36.3 Click on **SUBMIT**, this will generate a conjugation ID number

Click on back to **Conjugate** and filter by conjugation date

36.4 Record the ID number on your experiment worksheet [go to step #1](#)

36.5 Click on Edit Conjugate

Conjugate Details	Delete Conjugate	Edit Conjugate
view	delete	edit

36.6 Enter Nanodrop reading post-labelling and Volume of labeled antibody in the respective fields

Nanodrop post labelling (ug/ml)
<input type="text" value="495"/>
Volume of labeled antibody (ul)
<input type="text" value="140"/>
Labeled stock (ug/ml)
<input type="text" value="200"/>
Dilution factor
<input type="text" value="2.475"/>
Volume to add (ul)
<input type="text" value="206.50"/>
Total volume
<input type="text" value="347"/>

36.7 For a standard labeled stock concentration (200 µg/mL), the Dilution factor and the volume of Antibody Stabilizer is calculated automatically

If the Dilution factor is less than <1, change labeled stock concentration to 100 µg/mL

36.8 Add the volume of Antibody Stabilizer to the antibody labeled collection tube. Use Tris-based antibody stabilizer (Boca Scientific).

36.9 Transfer the diluted labeled antibody into a new screw cap microcentrifuge tube

36.10 Label the tube with printer Brady 21-Lab, (Label for 1.5 mL vial, Auto bold)

Line1: Target name, Channel

Line2: ID: number
Line3: date (mm/dd/yy)

36.11 Hand write the target name and the channel on a dot label and place on the screw cap

Short term storage at 4°C/ Long term storage Lyophilized

37 Labeled can be stored short term for <3 months at 4°C in solution

Preferably, the newly labeled antibody is immediately placed in the lyophilization box at 4°C

in queue to be lyophilized and then placed at 4°C in designated box for longterm storage management