



Aug 05, 2020

LC-MS/MS analysis

avinash.kale¹¹domnic colvin

1

Works for me

dx.doi.org/10.17504/protocols.io.bjdvki66

UM-DAE Centre for Excellence in Basic Sciences



avinash.kale

DOI

dx.doi.org/10.17504/protocols.io.bjdvki66

DOCUMENT CITATION

avinash.kale 2020. LC-MS/MS analysis. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bjdvki66>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 05, 2020

LAST MODIFIED

Aug 05, 2020

DOCUMENT INTEGER ID

40085

LC-MS/MS analysis

Nano-LCMS using 6540 UHD Accurate Mass QTOF (CAMS, Venture Center, NCL Pune, India) was used to carry out all our mass spectrometric analysis. The peptides were analyzed using reversed-phase nanoscale liquid chromatography coupled to tandem mass spectrometry. The system consisted of a column; initially, the sample would be subjected to Agilent Binary LC1290 using column Polaris, high-performance chip, 360 nl enrichment chip (150mm X 75µm), and a separation column G4240-62030 with a thermostat temperature of 4°C. The mobile phase was Water (0.1% Formic Acid) and Acetonitrile: Water (90:10 v/v). The chromatography was carried out in a gradient manner (3% to 90% Acetonitrile: Water). The mass spec analysis was carried out on Agilent 6540 UHD QTOF MS (Gas temp: 250, Gas flow (l/minute): 8 minutes,) MS analysis was carried out in a data-dependent manner with survey scan resolution. Mass range: 250-1700(m/z), MS/MS: 50-1700(m/z), Data type: Centroid and software: Masshunter workstation software. Obtained data were searched against the UniProt protein databases of all the respective four Bacillus sp. using default parameters.