



Version 2

Nov 30, 2020

Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 µg) (NEB #T1030) V.2

New England Biolabs¹¹New England Biolabs

1

Works for me

dx.doi.org/10.17504/protocols.io.bp9dmr26

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.comDanielle Freedman
NEB

ABSTRACT

Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 µg) ([NEB #T1030](#))

EXTERNAL LINK

<https://www.neb.com/protocols/2017/04/25/quick-protocol-for-oligonucleotide-cleanup-using-the-monarch-pcr-dna-cleanup-kit-5-g-neb-t1030>

ATTACHMENTS

[ProtocolCard_T1030.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bp9dmr26

EXTERNAL LINK

<https://www.neb.com/protocols/2017/04/25/quick-protocol-for-oligonucleotide-cleanup-using-the-monarch-pcr-dna-cleanup-kit-5-g-neb-t1030>

PROTOCOL CITATION

New England Biolabs 2020. Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 µg) (NEB #T1030). **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bp9dmr26>

Version created by Lenny Teytelman



EXTERNAL LINK

<https://www.neb.com/protocols/2017/04/25/quick-protocol-for-oligonucleotide-cleanup-using-the-monarch-pcr-dna-cleanup-kit-5-g-neb-t1030>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 30, 2020

LAST MODIFIED

Nov 30, 2020

OWNERSHIP HISTORY

Nov 30, 2020



Danielle Freedman

NEB

PROTOCOL INTEGER ID

45061

ATTACHMENTS

[ProtocolCard_T1030.pdf](#)

GUIDELINES

DNA Cleanup and Concentration: for the purification of up to 5 µg of DNA (ssDNA > 200 nt and dsDNA > 50 bp) from PCR and other enzymatic reactions. A [detailed protocol](#) and a [quick protocol](#) are available for your convenience.

Oligonucleotide Cleanup (steps): for the purification of up to 5 µg of DNA fragments ≥ 15 bp (dsDNA) or ≥ 18 nt (ssDNA). Expected recovery is > 70%. When purifying ssDNA of any size, recovery can be increased by using this protocol; however, it is important to note that this protocol shifts the cutoff for smaller fragments to 18 nt (rather than 50 nt for the [DNA Cleanup and Concentration Protocol](#)).

MATERIALS TEXT

MATERIALS

 **Monarch® PCR & DNA Cleanup Kit (5 µg) New England**

Biolabs Catalog #T1030

SAFETY WARNINGS

Please refer to SDS for safety warnings.

ABSTRACT

Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 µg) ([NEB #T1030](#))

BEFORE STARTING

- All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).
- Add 4 volumes of ethanol (≥ 95%) to one volume of DNA Wash Buffer.
- Please review the important information under the "Guidelines" & "Warnings" tabs before beginning.

- 1 Add 100 µl DNA Cleanup Binding Buffer to the 50 µl sample.

 **100 µl DNA Cleanup Binding Buffer**

We recommend a sample volume of 50 µl. For smaller samples, adjust the volume with nuclease-free water.

- 2 Add 300 µl ethanol (≥ 95%). Mix well by pipetting up and down or flicking the tube. Do not vortex.

 **300 µl Ethanol (≥ 95%)**

- 3 Insert column into collection tube and load sample onto column.

- 4 Spin for 1 minute, then discard flow-through.

 **00:01:00 Spinning**

- 5 Re-insert column into collection tube.

6 Add 500 µl DNA Wash Buffer and spin for 1 minute.

 **500 µl DNA Wash Buffer**

 **00:01:00 Spinning**

7 Discard flow-through.

8 **Repeat steps 5-7 (Optional).** Recommended for removal of enzymes that may interfere with downstream applications (e.g., Proteinase K).

 **Repeating steps 5-7**

9 **Transfer column to a clean 1.5 ml microfuge tube.** Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

 **00:01:00 Re-spinning**

10 Add ≥ 6 µl of DNA Elution Buffer to the center of the matrix.

 **6 µl DNA Elution Buffer**

11 Wait for 1 minute, then spin for 1 minute to elute DNA.

 **00:01:00 Waiting**

 **00:01:00 Spinning**

Typical elution volumes are 6–20 µl. Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated.