

Aug 21, 2024 Version 2

# 🌐 Imaging of cholinergic interneurons in post-mortem rodent tissue to identify striatal satellite astrocytes V.2

📁 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2](https://dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2)

Shinil Raina<sup>1,2,3</sup>, Stephanie J Cragg<sup>1,2,3</sup>

<sup>1</sup>Department of Physiology, Anatomy and Genetics, University of Oxford, Oxford OX1 3PT, UK;

<sup>2</sup>Oxford Parkinson's Disease Centre, University of Oxford, Oxford OX1 3PT, UK;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA


ASAP Collaborative Rese...

Team Cragg



Cláudia C. Mendes

University of Oxford



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2](https://dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2)

**Protocol Citation:** Shinil Raina, Stephanie J Cragg 2024. Imaging of cholinergic interneurons in post-mortem rodent tissue to identify striatal satellite astrocytes. [protocols.io https://dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2](https://dx.doi.org/10.17504/protocols.io.j8nlkkq31l5r/v2) Version created by Cláudia C. Mendes

**Manuscript citation:**

Stedehouder & Roberts et al. (2024) Rapid modulation of striatal cholinergic interneurons and dopamine release by satellite astrocytes, bioRxiv

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 19, 2022



**Last Modified:** August 21, 2024

**Protocol Integer ID:** 98956

**Keywords:** Histology, immunofluorescence, postmortem, immunohistochemistry, astrocytes, cholinergic interneurons, confocal imaging

**Funders Acknowledgement:**

**Medical Research Council  
(MRC)**

**Grant ID:** MR/V013599/1

**Aligning Science Across**

**Parkinson's (ASAP)**

**Grant ID:** ASAP-020370

## Abstract

This protocol details confocal image acquisition steps and the analysis pipeline used to calculate inter-soma distance between cholinergic interneurons/pMSNs and their nearest astrocyte. The partner **Protocol: Immunofluorescent Labelling of Post-Mortem Rodent Brain Tissue** describes how to label cholinergic interneurons/pMSNs along with striatal astrocytes in PFA perfusion-fixed, 50- $\mu$ m thick, post-mortem rodent brain tissue.

## Materials

**Equipment:**

- **ZEISS LSM 980 with Airyscan 2**

**Software:**

- **Zen Blue** (Zeiss, version 3.9)
- **FIJI** (version 1.54f)

**Other:**

- Slide with ~50  $\mu$ m thick brain section with fluorescent labels.

## Before start

The tissue used in the steps below were obtained after processing as described in **Protocol: Immunofluorescent Labelling of Post-Mortem Rodent Brain Tissue**.

## Confocal Imaging Acquisition

- 1 Turn on the Leica LSM980 confocal microscope and Zen Blue imaging acquisition software.
- 2 Secure a glass slide containing one or more fluorescently-labelled brain section into the stage of the microscope. Take care at this step to ensure the coverslip does not break or dislodge.
- 3 Using a 10x objective, identify anatomical landmarks using the DAPI or ChAT/NeuN signal in the region of interest. This would include the corpus callosum for the dorsal striatum, and the anterior commissure for the ventral striatum.

### Note

The DAPI signal would be visualised with a 358 nm laser. The ChAT/NeuN signal can both be visualised with a 568 nm or similar wavelength laser as appropriate for the secondary fluorescent antibody used.

- 4 Switch to a 63x/1.4 NA (oil immersion) objective and identify randomly selected ChAT/NeuN positive somata within the region of interest. Center the objective over this particular neuron.

### Note

Only include neurons for which the full soma resides within the z-axis of the tissue (*e.g.*, no cut neurons).

Take care to identify somata while blinded to the S100 $\beta$  (green) channel, to prevent bias to selecting ChAT/NeuN cells that have an astrocyte (as identified by a positive S100 $\beta$  signal) close to them.

- 5 Image the identified cell across the different channels using the following lasers:
  - ChAT/NeuN - 568 nm wavelength laser
  - S100 $\beta$  - 488 nm wavelength laser
  - DAPI - 358 nm wavelength laser

The laser wavelength can also be altered as appropriate for the secondary antibody that is used.

**Note**

The settings used were a 30  $\mu\text{m}$  thick z-stack (512 x 512 pixels) at 63x magnification with 1x digital zoom, 1  $\mu\text{m}$  optical thickness/optical plane, and a step size of 2  $\mu\text{m}$ . The z-stack centre was set to z-coordinate where the soma is largest (*e.g.* centre of the neuron).

- 6 Save the images with the appropriate details of the sample and imaging settings. A .czi format is recommended but other formats could also work.
- 7 Move to the next cell in the region of interest. Repeat image acquisition across the desired number of cells. Sample an equal number of cells from both hemispheres.

**Image Analysis**

- 8 Open the images in FIJI (ImageJ) and convert to a maximum projection across the z axis.
- 9 While remaining blind to the S100 $\beta$  (green) channel to prevent bias on astrocyte location, draw a region of interest (ROI) around the cell in the centre of the image with the polygon selection tool. This would either be the ChAT/NeuN positive cell, depending on the sample.
- 10 Open the S100 $\beta$  (green) channel. Draw a ROI around the nearest astrocyte by identifying a signal in the S100 $\beta$  (green) channel (>5  $\mu\text{m}$  in diameter) featuring a DAPI-positive nucleus residing closest to centre of the image.
- 11 From the ROIs, measure the surface area of the cells and the x and y coordinates of the centre point to calculate the intersoma distance.
- 12 Using the x and y coordinates obtained earlier, calculate the Euclidean distance between the centre point of the ChAT/NeuN cell and their nearest astrocyte.
- 13 Repeat analysis steps outlined above for all imaged cells.