



Oct 17, 2021

# Ligation and gel electrophoresis

Chia-Hsien Shih<sup>1</sup><sup>1</sup>Chung Shan Medical University[dx.doi.org/10.17504/protocols.io.by5dpy26](https://dx.doi.org/10.17504/protocols.io.by5dpy26)

CSMU\_Taiwan

CSMUTaiwan  
Chung Shan Medical University

This protocol is used for ligation of circular probe

DOI

[dx.doi.org/10.17504/protocols.io.by5dpy26](https://dx.doi.org/10.17504/protocols.io.by5dpy26)

Chia-Hsien Shih 2021. Ligation and gel electrophoresis. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.by5dpy26>



\_\_\_\_\_ protocol ,

Oct 17, 2021

Oct 17, 2021




54149

## Preparation

- 1 Dilute linear DNA into 10μM
- 2 Dilute T4 ligation primer into 10μM

## Protocol of Ligation





- 3 Add  11 μL of RNase-free water into a eppendorf

- 4 Add  **3  $\mu$ L** of 10X T4 Ligase buffer
- 5 Add  **9  $\mu$ L** of 10 $\mu$ M T4 ligation primer
- 6 Add  **3  $\mu$ L** of 10 $\mu$ M linear DNA
- 7 Add  **4  $\mu$ L** of T4 Ligase
- 8 Spin down after vortex
- 9 Incubate for  **04:30:00** at  **Room temperature**

4h 30m


#### Protocol of Gel Electrophoresis

30m

- 10 Put 2% agarose gel into the electrophoresis tank.
- 11 Prepare 0.5X TAE buffer ( add  **3.5 mL** 50X TAE buffer and  **350 mL** ddH<sub>2</sub>O and mix them), then pour into the electrophoresis tank. Make sure the gel is soaked into the buffer completely.
- 12 Take  **5  $\mu$ L** of marker and load it in the first well.
- 13 Take  **1  $\mu$ L** of 6X loading dye and drop on the parafilm separately, pipetting dye with

 **5 µL** of sample

14 Load samples into each well

15 Run gel at a constant voltage of 100V for  **00:30:00** at room temperature

30m

16 After gel electrophoresis, put the gel into gel reading machine and report the result.