

Aug 29, 2024

## Purification of Lambda Protein Phosphatase

DOI

**[dx.doi.org/10.17504/protocols.io.kqdg322bqv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg322bqv25/v1)**

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

---

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.kqdg322bqv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg322bqv25/v1)**

**Protocol Citation:** Elias Adriaenssens 2024. Purification of Lambda Protein Phosphatase. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.kqdg322bqv25/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 23, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 101120

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**Aligning Science Across  
Parkinson's (ASAP)**

**Grant ID:** ASAP-000350

**Marie Skłodowska-Curie  
MSCA Postdoctoral  
fellowship**

**Grant ID:** 101062916

## Abstract

This protocol details the purification of Lambda protein phosphatase.

## Materials

 Rosetta™(DE3)pLysS Competent Cells - Novagen **Merck Catalog #70956-4**

### Lysis buffer:

A	B
Tris-HCl, pH 7.4	50 mM
NaCl	300 mM
MgCl <sub>2</sub>	2 mM
Glycerol	5%
Imidazole	10 mM
β-mercaptoethanol	2 mM

### Wash buffer:

A	B
Tris-HCl, pH 7.4	50 mM
NaCl	300 mM
Imidazole	10 mM
β-mercaptoethanol	2 mM

### SEC buffer:

A	B
Tris-HCl, pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM



## Purification procedure

1d 1h 45m 30s

- 1 To purify Lambda protein phosphatase ( $\lambda$  PPase), fuse the protein phosphatase to a N-terminal 6xHis-tag through cloning into a pET-DUET1 vector (available from Addgene).
- 2 After the transformation of the pET-DUET1 vector encoding 6xHis-TEV- $\lambda$  PPase in E. coli Rosetta pLysS cells (Novagen Cat# 70956-4), grow cells in 2x Tryptone Yeast extract (TY)medium at 37 °C until an OD<sub>600</sub> of 0.4 and then continue at 18 °C .
- 3 Once the cells reached an OD<sub>600</sub> of 0.8, induce protein expression with 100 micromolar ( $\mu$ M) isopropyl  $\beta$ -D-1-thiogalactopyranoside (IPTG) for 16:00:00 at 18 °C . 16h
- 4 Collect the cells by centrifugation and resuspend in lysis buffer, complete EDTA-free protease inhibitors (Roche), CIP protease inhibitor (Sigma), and DNase (Sigma)).

### Lysis buffer:

A	B
Tris-HCl, pH 7.4	50 mM
NaCl	300 mM
MgCl <sub>2</sub>	2 mM
Glycerol	5%
Imidazole	10 mM
$\beta$ -mercaptoethanol	2 mM

- 5 Sonicate cell lysates twice for 00:00:30 . 30s
- 6 Clear lysates by centrifugation at 18000 rpm, 4°C, 00:45:00 in a SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific). 45m
- 7 Filter the supernatant through an 0.45  $\mu$ m filter and load onto a pre-equilibrated 5 ml His-Trap HP column (Cytiva).





- 8 After bind His-tagged proteins to the column, wash the column with three column volumes of wash buffer.

**Wash buffer:**



A	B
Tris-HCl, pH 7.4	50 mM
NaCl	300 mM
Imidazole	10 mM
β-mercaptoethanol	2 mM

- 9 Elute the proteins with a stepwise imidazole gradient (30, 75, 100, 150, 225, 300 mM).

- 10 Pool and incubate the fractions containing the 6xHis-TEV-λ PPase  Overnight with TEV protease at  4 °C .

8h



- 11 After the 6xHis tag was cleaved off, recapture 6xHis tag and His-tagged TEV protease with nickel beads for  01:00:00 at  4 °C .

1h

- 12 Pellet the beads by centrifugation and the supernatant, concentrate containing the λ PPase protein using a 30 kDa cut-off Amicon filter (Merck Millipore) and load onto a pre-equilibrated Superdex 200 Increase 10/300 GL column (Cytiva).

- 13 Elute the proteins with SEC buffer.

**SEC buffer:**

A	B
Tris-HCl, pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM


- 14 Analyse the fractions by SDS-PAGE and Coomassie staining.



- 15 Pool the fractions containing purified λ PPase.



16 After concentrating the purified protein, aliquot the protein and snap-frozen in liquid nitrogen.

17 Store the proteins at  -80 °C .