



MAY 31, 2023

LC3 Lipidation assay for ATG3 Mutants

Imstrong¹¹UC Berkeley

Imstrong

ABSTRACT

LC3 lipidation assay on liposomes

OPEN  ACCESS**DOI:**dx.doi.org/10.17504/protocols.io.e6nvwjxodlmk/v1

Protocol Citation: Imstrong 2023. LC3 Lipidation assay for ATG3 Mutants.

protocols.io<https://dx.doi.org/10.17504/protocols.io.e6nvwjxodlmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 31, 2023**Last Modified:** May 31, 2023**PROTOCOL integer ID:**
82662**Keywords:** ASAPCRN

- 1 Mix 1 μ M ATG3 WT or mutant, 1 μ M ATG7, 500 nM WIPI2d, 1 μ M ATG12–ATG5–ATG16L1, 5 μ M LC3B

and 0.5 mg/mL extruded liposomes in reaction buffer (25mM HEPES pH 7.4, 150 mM NaCl, 1 mM TCEP).

- 2** Incubate at room temperature for 20 minutes.
- 3** Quench with 4X LDS loading buffer. Boil samples for 10 minutes at 60°C. Load and run SDS PAGE gel
- 4** Visualize protein bands with coomassie blue