



Sep 16, 2020

# High-Performance Liquid Chromatography (HPLC)

In 1 collection

Yingchao Xue<sup>1,2</sup>, Xiping Zhan<sup>3</sup>, Shisheng Sun<sup>4</sup>, Senthilkumar S. Karuppagounder<sup>5,6,7</sup>, Shuli Xia<sup>2,5</sup>, Valina L Dawson<sup>5,6,7,8,9</sup>, Ted M Dawson<sup>5,6,7,8,10</sup>, John Laterra<sup>2,5,8,11</sup>, Jianmin Zhang<sup>1</sup>, Mingyao Ying<sup>2,5</sup>

<sup>1</sup>Department of Immunology, Research Center on Pediatric Development and Diseases, Institute of Basic Medical Sciences, Chinese Academy of Medical Sciences and School of Basic Medicine, Peking Union Medical College, State Key Laboratory of Medical Molecular Biology;

<sup>2</sup>Hugo W. Moser Research Institute at Kennedy Krieger; <sup>3</sup>Department of Physiology and Biophysics, Howard University;

<sup>4</sup>College of Life Sciences, Northwest University; <sup>5</sup>Department of Neurology, Johns Hopkins University School of Medicine;

<sup>6</sup>Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine;

<sup>7</sup>Adrienne Helis Malvin Medical Research Foundation; <sup>8</sup>Department of Neuroscience, Johns Hopkins University School of Medicine;

<sup>9</sup>Department of Physiology, Johns Hopkins University School of Medicine;

<sup>10</sup>Department of Pharmacology and Molecular Sciences, Johns Hopkins University School of Medicine;

<sup>11</sup>Department of Oncology, Johns Hopkins University School of Medicine

1 Works for me This protocol is published without a DOI.

## Neurodegeneration Method Development Community

Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



Anita Broellochs  
protocols.io

## ABSTRACT

This protocol describes High-Performance Liquid Chromatography (HPLC) for lines ND1014, N1, and ND27760 from *Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons*.

## EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

## PROTOCOL CITATION

Yingchao Xue, Xiping Zhan, Shisheng Sun, Senthilkumar S. Karuppagounder, Shuli Xia, Valina L Dawson, Ted M Dawson, John Laterra, Jianmin Zhang, Mingyao Ying 2020. High-Performance Liquid Chromatography (HPLC). **protocols.io**  
<https://protocols.io/view/high-performance-liquid-chromatography-hplc-9zdh726>

## MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

## EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

## COLLECTIONS ①



**Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons**

## KEYWORDS

ND1014, N1, ND27760, ipsc, SNCA, Atoh2, Ngn2, HPLC

#### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### CREATED

Dec 02, 2019

#### LAST MODIFIED

Sep 16, 2020

#### OWNERSHIP HISTORY

Dec 02, 2019  Liz Brydon [Protocols.io](#)

Sep 16, 2020  Anita Broellochs [protocols.io](#)

#### PROTOCOL INTEGER ID

30469

#### PARENT PROTOCOLS

Part of collection

[Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons](#)

#### MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">Potassium chloride</a>	P9333	<a href="#">Sigma Aldrich</a>
<a href="#">Perchloric acid</a>	244252	<a href="#">Sigma – Aldrich</a>
<a href="#">Savant SPD121P SpeedVac Kit</a>	SPD121P	<a href="#">Thermo Scientific</a>
<a href="#">Savant Refrigerated Vapor Traps</a>	RVT5105	<a href="#">Thermo Scientific</a>

#### SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.

#### BEFORE STARTING

Obtain approval to work with human stem cells from an appropriate Institutional Review Board.

#### HPLC Preparation

- 1 Replace neuron culture medium with Hanks' balanced saline solution buffer with the addition of **[M]56 Milimolar (mM) KCl**.



2

Incubate for **00:15:00** at **37 °C**.



3

Collect media and centrifuge to clear cell debris. Collect neuron pellet.



4


Freeze immediately and store at  $-80^{\circ}\text{C}$ .

#### HPLC Analysis

5 Thaw samples.

6 Concentrate using a vacuum (Savant SPD 121P) connected with a refrigerated vapor trap (Savant RVT 5105).

7 Resuspend freeze-dried samples in  $10\text{ mM}$  perchloric acid.

8 

Analyze monoamines using HPLC electrochemical detection by dual channel Coulochem III electrochemical detector (model 5300).

9 Separate monoamines using a reverse phase C18 column (3-mm  $\times$  150-mm C-18 RP-column) with a flow rate of 0.600 ml/minute.

10 

Quantify monoamine concentrations by comparison of the area under the curve to known standard dilutions.