



Sep 10, 2020

Freezing/Cell Banking

In 1 collection

Andrea Argouarch¹¹University of California, San Francisco

1 Works for me dx.doi.org/10.17504/protocols.io.8gmhtu6

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Andrea Argouarch

University of California, San Francisco

ABSTRACT

Protocol includes banking dural cells into cryovials for long term storage in liquid nitrogen.

DOI

dx.doi.org/10.17504/protocols.io.8gmhtu6

PROTOCOL CITATION

Andrea Argouarch 2020. Freezing/Cell Banking. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8gmhtu6>

COLLECTIONS ⓘ

Dural Cell Isolation and Culturing - Collection

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 19, 2019

LAST MODIFIED

Sep 10, 2020

PROTOCOL INTEGER ID

28909

PARENT PROTOCOLS

Part of collection

[Dural Cell Isolation and Culturing - Collection](#)

STEPS MATERIALS

NAME	CATALOG #	VENDOR
Mr. Frosty Freezing Container, 2mL tubes, Nalgene		
Mr. Frosty Freezing Container for 1-2mL cryogenic tubes, PC, clear w/ blue lid, 1/Cs.	5100-0001	Thermo Fisher
Cryo-Tags	9187-1100	USA Scientific
Trypan Blue Solution, 0.4%	15250061	Thermo Fisher
Penicillin-Streptomycin	15140122	Gibco - Thermo Fisher
Fetal Bovine Serum	97068-091	Vwr

NAME	CATALOG #	VENDOR
DMSO Bio-Max, Cell Culture Grade	40470005-2	bioworld
C-Chip [®] ; Disposable Hemacytometers, Improved Neubauer, 50/Pk.	22600100	Thermo Fisher
DPBS, no calcium, no magnesium	14190250	Thermo Fisher
DMEM, high glucose, pyruvate	11995073	Thermo Fisher
Trypsin-EDTA (0.05%), phenol red	25300062	Thermo Fisher

Observations

1 At 100% confluency, take pictures of cell morphology, freeze down for banking (~0.5 Million per vial), collect remaining cells in flask for gDNA, and seed cells into 1xT25 for karyotyping (~30,000 cells, 60 ul)

2 Turn off UV lights and clean hood with 70% ethanol

3 Clean items with 70% ethanol and bring into hood
a. DPBS -/-



DPBS, no calcium, no magnesium

by Thermo Fisher

Catalog #: 14190250

b. Sterile Filtered Media



DMEM, high glucose, pyruvate

by Thermo Fisher

Catalog #: 11995073



Penicillin-Streptomycin

by Gibco - Thermo Fisher

Catalog #: 15140122



Fetal Bovine Serum

by Vwr

Catalog #: 97068-091

c. Label 2xT150 flasks with ID, date, and p2, add 28 mls 28 mL of media per flask

d. 0.05% Trypsin



Trypsin-EDTA (0.05%), phenol red

by Thermo Fisher

Catalog #: 25300062

Preparation

- 4 Turn off UV lights and clean hood with 70% ethanol

- 5 Labeling cryovials
 - a. Create labels using Microsoft Word, following the formatting instructions



Cryo-Tags

by USA Scientific

Catalog #: 9187-1100

- b. Print out estimated number of labels for the both flasks (~15) with ID, date, p2, and ~0.5 mill cells.
 - c. Arrange cryovials in 96, 1.5ml tube holder with spaces in between each vial
 - d. Spray label sheets thoroughly with 70% ethanol, dry, and bring into hood to adhere onto each vial
 - e. Label top of each vial with ID
 - f. Unscrew cap in order to easily lift up when pipetting in cell suspension
- 6 Clean items with 70% ethanol and bring into hood
 - a. DPBS ^{-/-}



DPBS, no calcium, no magnesium

by Thermo Fisher

Catalog #: 14190250

- b. Sterile Filtered Media



DMEM, high glucose, pyruvate

by Thermo Fisher

Catalog #: 11995073



Penicillin-Streptomycin

by Gibco - Thermo Fisher

Catalog #: 15140122



Fetal Bovine Serum

by Vwr

Catalog #: 97068-091

c. 0.05% Trypsin

- i. Aliquot to and warm before use



Trypsin-EDTA (0.05%), phenol red

by Thermo Fisher

Catalog #: 25300062

d. Freezing media (made fresh)

- i. 100% DMSO to 10%



DMSO Bio-Max, Cell Culture Grade

by bioworld

Catalog #: 40470005-2

- ii. 90% Media

e. Incyto disposable hemocytometers



C-Chip™ Disposable Hemacytometers, Improved Neubauer; 50/Pk.

by Thermo Fisher

Catalog #: 22600100


f. Trypan blue, add 40 ul 40 µl into a 1.5 ml tube











Trypan Blue Solution, 0.4%











by Thermo Fisher





Catalog #: 15250061

- g. T25 labeled to send to Cell Line Genetics for karyotyping (KT)
- h. 15 ml conical with 10 mls  10 mL of media to dilute KT cells
- i. 15 ml conical to collect cells for gDNA

Culturing

- 7 Take 3-4 BF images of morphology. Save TIFF under ID and passage number
- 8 Aspirate old media
- 9 Rinse by adding 5 mls  5 mL of DPBS per flask and gently swirling
- 10 Aspirate DPBS
- 11 Add 3 ml  3 mL of trypsin per flask and place in incubator  37 °C for 3 mins  00:03:00 until cells are detached
 - a. Can also gently tap the side of the flask to detach cells, check under microscope to confirm
- 12 Add 9 mls  9 mL of fibroblast media per flask to inactivate trypsin
- 13 Collect cell suspension in a 50 ml conical
- 14 Rinse flask with 12 mls  12 mL of PBS to collect remaining cells, place in 15ml conical for gDNA
- 15 Measure final cell suspension volume with serological and mix
- 16 Dilute 10 ul  10 µl of cell suspension into trypan blue, tap to mix
- 17 Load 10 ul  10 µl of trypan blue and cell dilution into the hemocytometer
- 18 Count non blue cells in 4 grind corners during spin and then calculate total cells and resuspension volume for 0.5E6 cells per ml

- 19 Total cells for 2 flasks = (Cell Count/4 grinds) x 5 dilution factor x 10,000 x ~25mls suspension volume
a. Resuspension volume for 0.5E6 cells per ml = Number of vials
- 20 Spin both conicals at 1000 rpm  **1000 rpm** for 5 mins  **00:05:00** during counting
- 21 Aspirate supernatant for both pellets (freezing and gDNA), being careful to not aspirate the pellet
- 22 Tap the pellet to resuspend for freezing
- 23 Store gDNA pellet at -80C  **-80 °C** until mycoplasma assay
- 24 Add freezing media dropwise to cell pellet and add 1 ml  **1 mL** per cryovial. Spin KT cell pellet during this time.
- 25 Prepare for KT in T25 flask
 - a. Take 60-100 ul  **60 µl** of DMSO cell suspension and dilute into 10 mls  **10 mL** of media
 - b. Spin for 1000 rpm  **1000 rpm** for 5 mins  **00:05:00**
 - c. Aspirate, resuspend in 5 ml  **5 mL** of media, and add to 1xT25 for KT
 - d. Place flask in incubator (5% CO2, 37°C)  **37 °C**
 - e. Send KT to Cell Line Genetics at 40-50% confluency in about 5-6 days
- 26 Secure each cryovial cap tightly and place in Mr. Frosty

 Mr. Frosty Freezing Container, 2mL tubes, Nalgene Mr. Frosty Freezing Container for 1-2mL cryogenic tubes, PC, clear w/ blue lid, 1/Cs.
by Thermo Fisher
Catalog #: 5100-0001
- 27 Place Mr. Frosty in -80°C  **-80 °C** for 24-48 hours  **48:00:00** and then transfer to liquid nitrogen  **-190 °C**

- 28 Freeze back up cells from 6 well plate, place in Mr. Frosty for 24-48 hours ⌚ **48:00:00** , then transfer to LN2
🧊 **-190 °C**

Clean Up

- 29 Throw away biohazard materials properly
- 30 Clean and sterilize hood with 70% ethanol and turn on UV
- 31 Update cell culture notes in lab notebook