

FEB 08, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io. 81wgbxbonlpk/v1

Protocol Citation: james.evans 2024. ASO transfection of iPSC-derived cells. **protocols.io** https://dx.doi.org/10.17504/protocols.io.81wgbxbonlpk/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Feb 08, 2024

ASO transfection of iPSC-derived cells

james.evans¹

¹University College London, University of London



james.evans

University College London, University of London, The Francis...

ABSTRACT

A protocol for transfecting iPSC-derived midbrain dopaminergic neurons with antisense oligonucleotides (ASOs).

Oct 8 2024



Last Modified: Feb 08, 2024

PROTOCOL integer ID: 94881

Keywords: transfection, ASO

ASO transfection of iPSC-derived mDA neurons

- 1 Make up transfection mix in N2B27. The transfection mix should be 1/5 of the final volume in the well. Transfection mix should have:
 - 1. 0.48 % <u>DharmaFECT</u> transfection reagent
 - 2. ASO (adjust concentration depending on cell type, ASO chemical modification, and knockdown required). Calculate final concentration required in the well not in the transfection mix.

Example - 1000 ul final volume in the well with 300 nm of ASO = 200 ul of N2B27 + 0.96 ul of DharmaFect + 3 ul of 100 uM ASO Stock

- 2 Vortex transfection mix and leave at room temperature for 30 minutes.
- Aspirate media from cells and replace with the transfection mix. After adding the transfection mix add the rest of the media to the cells to make up the final volume.