



Upload image

Jul 27, 2020

# Phenol-based RNA extraction from polycarbonate filters

Robbie M. Martin<sup>1</sup>, Steven W. Wilhelm<sup>1</sup><sup>1</sup>University of Tennessee, Knoxville

1

Works for me

[dx.doi.org/10.17504/protocols.io.bivuke6w](https://dx.doi.org/10.17504/protocols.io.bivuke6w)[The Aquatic Microbial Ecology Research Group - AMERG \(The Buchan, Zinser and Wilhelm labs\)](#)

Robbie Martin

## ABSTRACT

This protocol is an acid-phenol-based method for extracting RNA from samples collected onto polycarbonate filters. Bead-beating is included to increase yields from difficult to lyse cells, for example, in some species of cyanobacteria. We have used it successfully in RNA-sequencing projects involving lab cultures of cyanobacteria and from freshwater and marine environmental samples.

## DOI

[dx.doi.org/10.17504/protocols.io.bivuke6w](https://dx.doi.org/10.17504/protocols.io.bivuke6w)

## PROTOCOL CITATION

Robbie M. Martin, Steven W. Wilhelm 2020. Phenol-based RNA extraction from polycarbonate filters.

**protocols.io**[dx.doi.org/10.17504/protocols.io.bivuke6w](https://dx.doi.org/10.17504/protocols.io.bivuke6w)

## KEYWORDS

RNA extraction, polycarbonate filters, bead beating, cyanobacteria

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 22, 2020

## LAST MODIFIED

Jul 27, 2020

## PROTOCOL INTEGER ID

39572

## GUIDELINES

This protocol assumes that sample biomass has been appropriately collected onto polycarbonate filters and properly stored at -80 °C.

## MATERIALS TEXT

frozen polycarbonate filters containing sample biomass

ice

Solution A (lysing buffer) (0.5% SDS, 20 mM sodium acetate, 10 mM EDTA)

Lysing Matrix E tubes (bead beating tubes)

basic phenol with 8-hydroxyquinoline (pH 7.8-8.2)

acid phenol:chloroform:isoamyl alcohol (125:24:1; pH 4.8)

chloroform

100% ethanol

70% ethanol  
RNase-free water

#### SAFETY WARNINGS

Phenol is a very hazardous chemical. Use extreme caution and follow recommended manufacturer and institutional safety guidelines.

#### Phenol-based RNA Extraction

- 1 Add 700  $\mu$ L of Solution A (lysing buffer) to an MP Biomedical Lysing Matrix E tube (or equivalent bead beating tube). Solution A consists of 0.5% SDS, 20 mM sodium acetate, and 10 mM EDTA in RNase-free molecular biology grade water.
- 2 Pre-chill tube containing Solution A on ice to prevent tube from overheating during bead beating.
- 3 Unfold polycarbonate filter (either frozen, or thawed but maintained on ice) and place into pre-chilled bead-beating tube.
- 4 Agitate tube on maximum speed for 40 s in Biospec Products Mini-beadbeater (or equivalent), then place on ice.
- 5 Add 500  $\mu$ L basic phenol with 8-hydroxyquinoline (pH 7.8-8.2) plus 100  $\mu$ L chloroform to tube.
- 6 Agitate tube on maximum speed for 40 s, then place briefly on ice to chill.
- 7 Centrifuge tube at maximum speed in benchtop centrifuge for 5 min.
- 8 Transfer as much of aqueous supernatant as possible from bead beating tube to new 2 mL centrifuge tube.
- 9 Re-extract aqueous supernatant twice with 500  $\mu$ L of acid phenol:chloroform:isoamyl alcohol (125:24:1; pH 4.8).
- 10 Re-extract aqueous supernatant once with 500  $\mu$ L basic phenol with 8-hydroxyquinoline plus 100  $\mu$ L chloroform.
- 11 Re-extract twice with 500  $\mu$ L pure chloroform to remove last traces of residual phenol.

## Ethanol Precipitation

- 12 Add sodium acetate to aqueous supernatant to a concentration of 0.3 M and mix.
- 13 Add 2-2.5x volume of ice-cold 100% ethanol and mix well.
- 14 Place in -80 °C freezer for at least 1 h, or preferably overnight.
- 15 Centrifuge at maximum speed at 0 °C for at least 30 min in refrigerated benchtop centrifuge. If you expect low yields, or wish to increase recovery, spin for 1 h at maximum.
- 16 Carefully decant ethanol supernatant without disturbing RNA pellet.
- 17 Wash pellet by adding 1 mL of ice-cold 70% ethanol and gently pipetting to dislodge RNA pellet.
- 18 Centrifuge at maximum speed at 0 °C for at least 15 min in a refrigerated benchtop centrifuge.
- 19 Gently and carefully aspirate as much of the 70% ethanol from the tube as possible taking caution not to accidentally aspirate the RNA pellet.
- 20 Place the open tube in a 37 °C heating block for 3-5 min to evaporate all ethanol. Do not dry pellet to completion, as this can make RNA difficult to redissolve in water.
- 21 Dissolve the RNA pellet in an amount of RNase-free water appropriate for your expected yield and needed final concentration, typically about 50-100 µL.