



VERSION 3

APR 19, 2023

OPEN  ACCESS

Protocol Citation: Haim Barr, Noa Lahav 2023. MERS-CoV Main Protease (Mpro) Fluorescence Dose Response. **protocols.io** <https://protocols.io/view/mers-cov-main-protease-mpro-fluorescence-dose-resp-cs48wgzw> Version created by [ASAP Discovery](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Apr 19, 2023

Last Modified: Apr 19, 2023

PROTOCOL integer ID:
80768

MERS-CoV Main Protease (Mpro) Fluorescence Dose Response V.3

Haim Barr^{1,2}, Noa Lahav^{1,2}

¹The Weizmann Institute of Science; ²ASAP Discovery



ASAP Discovery

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This is a functional, biochemical assay used to identify treatments for viral infectious disease in MERS-CoV 3C-like protease.

Utilizing a direct enzyme activity measurement method, the experiment was performed in a 384-well plate reading the fluorescence intensity. This assay tested the mode of action of inhibition.

Final Experiment Concentrations

A	B	C
Reagent	Concentration	Units
MERS Mpro	50	nM
MERS Substrate peptide	550	nM
HEPES pH=7.3	20	mM
NaCl	50	mM
BSA	0.1	mg/ml
Triton X-100	0.01	% (v/v)
TCEP	1	mM

For more information, please see the Materials section.

GUIDELINES

Plate Information:

Total Assay Volume: 20 μ L

Compounds Top Assay Concentration: 100 μ M

Dilution Factor: 2

Dose Response Points: 12

Number of Replicates: 2

Backfill with DMSO: Yes

MATERIALS

Assay Buffer Reagents (Concentration listed are the stock concentrations)

- [M] 40 millimolar (mM)
⊗ HEPES Buffer (pH 7.3) Fisher Scientific Catalog #BP299-1 (or similar)
- [M] 100 millimolar (mM)
⊗ Sodium Chloride Fisher Scientific Catalog #S271 (or similar)
- [M] 10 mg/mL
⊗ BSA-Molecular Biology Grade - 12 mg New England Biolabs Catalog #B9000S
(or similar)
- [M] 10 % volume
⊗ Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML
(or similar)
- [M] 1000 millimolar (mM) | ⊗ TCEP HCl P212121 Catalog #SV-TCEP (or similar)

(all components are added fresh to the assay buffer before each experiment)

Additional Reagents

- [M] 507000 nanomolar (nM) **MERS Mpro Enzyme***

***Note:** The original MERS Mpro stock enzyme had a concentration of [M] 507000 nanomolar (nM) but when a new stock solution was made/delivered on 2023-02-21 the new stock had a concentration of [M] 478000 nanomolar (nM). Both stock solutions were diluted with fresh assay buffer to create a [M] 100 nanomolar (nM) solution before each experiment.

- [M] 750000 nanomolar (nM) **MERS Substrate***

***Note:** MERS Substrate (5-FAM)-GVLQSGLV-K(Dabcyl)-K-NH₂ Stock was purchased from Peptide 2.0 and dissolved in DMSO with an original concentration of [M] 750000 nanomolar (nM) however it was diluted with the same assay buffer used in the experiment before conducting each experiment to yield a concentration of [M] 1100 nanomolar (nM).

SAFETY WARNINGS



Please be sure to wear proper Personal Protective Equipment (PPE) while performing this experiment.

BEFORE START INSTRUCTIONS

Note: Inhibitor compounds stock concentration is **20 mM**. Compounds are pre-dispensed into 384 plates and stored at -20 C until use.

Prepare Reagents

1 PREPARE all of the reagents/buffers required for this experiment.

A	B	C	D	E
Reagent	Stock Concentration	Concentration Loaded into Combi	Final Concentration	Units
MERS Mpro Enzyme (original stock)	507000	100	50	nM
20230221 MERS Mpro Enzyme	478000	100	50	nM
MERS Substrate	750000	1100	550	nM
Assay Buffer				
HEPES (pH 7.3)	40	20	20	mM
Sodium Chloride	100	50	50	mM
BSA	10	0.1	0.1	mg/mL
Triton X-100	10	0.01	0.01	% by volume
TCEP	1000	1	1	mM

For more information, please see the Materials sec

Prepare 384 Well Plate

15m

2 PRIME with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely.

2.1 DISPENSE  10 µL Assay Buffer to Columns **1** and **23** of assay plate

- **Note:** These will represent the ***inhibitor control columns*** (Contain: Substrate, Assay Buffer, DMSO, **no experimental compounds**)

2.2 EMPTY Multi-Drop Combi Tube Dispensing Cassette (by selecting the **EMPTY** button on the Combi Dispenser until the tubes of the cassette are emptied).

- Discard Assay Buffer discharged from the cassette.

3 PRIME with **100 nanomolar (nM) MERS MPro Enzyme** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely.

- **Note:** Be sure to cycle dispensing several times on a clean plate lid (This confirms there are no bubbles in the Dispensing Cassette).

3.1 DISPENSE **10 µL** **100 nanomolar (nM) MERS MPro Enzyme** to Columns **2 through 22** and Column **24**

Note:

- **100 nanomolar (nM) MERS MPro Enzyme** is two times the final concentration for the assay. It is diluted to be a final concentration of **50 nanomolar (nM) MERS MPro Enzyme**.
- Column 2 and Column 24 are **neutral control columns** (Contain: Enzyme, Substrate, DMSO, no compounds)

3.2 EMPTY Multi-Drop Combi Tube Dispensing Cassette (by selecting the **EMPTY** button on the Combi Dispenser until the tubes of the cassette are emptied.)

- Discard the **100 nanomolar (nM) MERS MPro enzyme** discharged from the cassette.

4 CENTRIFUGE **15000 rpm, Room temperature, 00:01:00** plate to remove bubbles

1m

5 INCUBATE plate for **00:15:00** at **Room temperature**


15m

⚠ Make sure the plate is protected from light!

6 PRIME with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely. Then, **EMPTY** the Multi-Drop Combi Tube Dispensing Cassette.


7 PRIME with **1100 nanomolar (nM) MERS Substrate** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely.

- **Note:** Be sure to cycle dispensing several times on a clean plate lid (This confirms there are no bubbles in the Dispensing Cassette).


8 DISPENSE  10 µL 1100nM MERS Substrate into Columns **1 through 24** (the full plate)

Note:

- **[M] 1100 nanomolar (nM) MERS Substrate** is two times the final concentration for the assay. It is diluted to be a final concentration of **[M] 50 nanomolar (nM) MERS Substrate**

9 CENTRIFUGE  15000 rpm, Room temperature, 00:01:00 plate in plate centrifuge to remove bubbles

1m

10 INCUBATE plate for  01:00:00 at  Room temperature

1h

⚠ Make sure the plate is protected from light!

Recommended: Clean the Multi-Drop Combi Reagent Dispenser during this incubation step

Read Plate Fluorescence

11 READ and RECORD the plate Relative fluorescence units (RFU) via the **"MERS Protocol"** on the **PHERASTAR FS Control Software**.

- Software is a standard Fluorescence Assay set for Optimal excitation wavelength 485 nm, emission wavelength 528 nm, and a Gain of 300.

Equipment

PHERASTAR FS	NAME
Microplate reader	TYPE
BMG LABTECH	BRAND
0471B0001A	SKU
https://www.bmglabtech.com/en/pherastar-fsx/?utm_term=pherastar%20plate%20reader&utm_campaign=usa.roi.products&utm_source=adwords&utm_medium&gclid=Cj0KCQjw8qmhBhCIARIsANAtbodGRjigZtEYwcoMXUtxsLn25xp4gjKra3ZNt9jLh9-FwOoFR_5EUHUaAlkREALw_wcB	LINK

Expected result

Gain 300 should yield ~10,000 RFU in full reaction and ~6,000 RFU in Buffer Control