



Mar 25,
2020

Protocol for mouse perfusion with dye, DAPI staining, and slide preparation

Forked from [Protocol for mouse perfusion with dye, DAPI staining, and slide preparation](#)

In 1 collection

Boaz Mohar¹, Monique Copeland¹

¹HHMI Janelia Research Campus

1 Works for me dx.doi.org/10.17504/protocols.io.59jg94n



Boaz Mohar
HHMI Janelia Research Campus






MATERIALS

| NAME | CATALOG # | VENDOR |
|---|-----------|------------------------------|
| Phosphate Buffered Saline | 28374 | Thermo Fisher Scientific |
| Paraformaldehyde fixative: 4% paraformaldehyde in phosphate buffered saline (PBS) | | |
| Falcon Tube (50 mL) | | Fischer Scientific |
| DMSO | | |
| Disposable gloves, nitrile | | |
| Surflo Winged Infusion Set 21Gx3/4 | SV*21BLK | Terumo |
| Surgical Scissors | 14090-09 | Fine Science Tools |
| Standard Pattern Forceps | 11001-12 | Fine Science Tools |
| Forceps | 11223-20 | Fine Science Tools |
| Isoflurane | CED1360 | Penn Veterinary Supply, Inc. |

STEPS MATERIALS

| NAME | CATALOG # | VENDOR |
|--|--------------|------------------------------|
| Falcon Tube (50 mL) | | Fischer Scientific |
| Dimethyl sulfoxide (DMSO) | D2650 | Sigma Aldrich |
| PBS | BP24384 | Fisher Scientific |
| PBS | BP24384 | Fisher Scientific |
| Isoflurane 1-3% | 07-893-1389 | Patterson Veterinary |
| Agarose LE | 87046 | Electron Microscopy Sciences |
| PBS | BP24384 | Fisher Scientific |
| Silver blue double edge razor blade (Gillette) | 884718268996 | Amazon |
| cell culture plate 24 well | View | Sigma-aldrich |
| DAPI | D1306 | Thermo Scientific |

| NAME  | CATALOG #  | VENDOR  |
|--|---|--|
| Superfrost Plus Microscope Slies | 4951PLUS4 | Thermo Fisher Scientific |
| VECTASHIELD® Hardset™ Antifade Mounting Medium | H-1400 | Vector Laboratories |

Prepare dye and fix


- 1 Add  **100 μ l** of DMSO to lyophilized dye (100nmol) and Pipette up and down ~20 times to get 1mM concentration stock.



Dimethyl sulfoxide (DMSO)

by Sigma Aldrich

Catalog #: D2650

- 2 Mix further using a vortex machine to ensure the dye is evenly distributed in the DMSO
- 3 Prepare 4% Paraformaldehyde in 0.1M PB [see here](#)
- 4 Add 50-100 ul of 1mM dye stock to  **50 ml per mouse** of 4% paraformaldehyde in .1M PB

Set-up hood for perfusion

5 Set up fume hood with:



Peristaltic pump-102R
Watson-Marlow 013.7101.000



Dissection tray
VWR 470006-956



Forceps
Fine Science tools 11223-20



Surgical Scissors
Fine Science tools 14090-09



Falcon Tube (50 mL)
by Fischer Scientific

Perfusion and dissection


6

Perfusion pump preparation

Attach needle:



Surflo winged infusion set 21G 3/4
Terumo TER 3SV-21BLK

to pump tubing and flush with ~  20 ml of 1x PBS




PBS
by Fisher Scientific
Catalog #: BP24384

to clear tubing of previous solution

7 Retrieve animal per vivarium protocol

8 Anesthetize mouse in 2000ml glass beaker

8.1 pour  5 ml isoflurane down the side of the container.



Isoflurane 1-3%

by Patterson Veterinary

Catalog #: 07-893-1389

- 8.2 Place the animal inside the container on a raised platform and cover.
(Platform can be square slides box lid or styrafoam from 50ml tube

packaging)

(DO NOT allow the liquid isofluorane to come into contact with the animals' skin)

- 8.3

As soon as the animal stops breathing, remove from contain, check for pain response by gently pinching the foot with forceps and make sure no reflexive response and place on dissection tray.

Timing is extremely important. If the animals is over-exposed to the isofluorane than the heart will stop beating and cell death can occur.

- 9 Begin perfusion

- 9.1 Lay mouse on its back and place a pin in each front paw, others as needed.

- 9.2 Use scissors to cut parallel to the sternum, opening the chest cavity and exposing the heart.

- 9.3 Insert winged needle into the right ventricle.

- 9.4 Cut left atrium.

- 9.5 Start perfusion pump and flush 50 ml 4% PFA with 100 µl of Dye at a rate of 7.5ml/min.

(Animal will "twitch" when the fixative circulates through the bloodstream. If no reaction occurs check postion of the needle in the heart.)

- 9.6 "Stop" perfuison pump when all fixative is flushed through the animal.

- 9.7 Unpin mouse from dissection tray and decapitate behind ears.
(Place carcass in small biohazard bag for disposal.)

- 9.8 Remove the fur from the head and cut small slits into the skull by placing scissor tip on eithier side of the opening where the spinal cord connect to the brain.

- 9.9 Use forceps to peel off skull cap and scope out brain.

- 9.10 Post-fix brain in 20 ml 4% PFA 16:00:00 at 4 °C

9.11 Rinse 3 times 🕒 00:15:00 🕒 00:15:00 🕒 00:15:00 in PBS



PBS

by Fisher Scientific

Catalog #: BP24384

at 🌡 Room temperature

Prep brain for sectioning

10 Prepare 🧴 50 ml 4% agarose

10.1 Add 🧴 2.5 g



Agarose LE

by Electron Microscopy Sciences

Catalog #: 87046

into 🧴 50 ml 1x PBS



PBS

by Fisher Scientific

Catalog #: BP24384

in 500ml glass bottle with cap.
(do not fully tighten cap onto bottle so steam can vent while heating)

10.2 Heat in microwave for 🕒 00:01:20 .

Check to see if all agarose is dissolved. If not re-heat at 🕒 00:00:10 intervals till completely dissolved.

10.3 Pour agarose into



Peel-A-Way® Embedding Mold (Square - S22)
polysciences 18646A-1

and orient brain in mold.

10.4 Allow agarose to set about 🕒00:10:00

(will change from transparent to slightly cloudy and will become firm to the touch.)

11 Prep Vibratome



VT 1200s Vibratome
Leica 1491200S001 [🔗](#)

with



Silver blue double edge razor blade
(Gillette)
by Amazon
Catalog #: 884718268996

11.1 Adjust vibratome settings:

Amplitude 1.00mm
Speed .70mm/s
Feed (thickness) 100um

11.2 Fill buffer tray with 🧴300 ml 1x PBS

Section brain

12 Section desired region at 50um or 100um.

12.1 Float sections in



cell culture plate 24 well

by Sigma-aldrich

[View](#)

(Single section per well) with each well filled with **500 µl 1x PBS**

(Label plate and wrap with aluminum foil)

DAPI stain

13 DAPI stain

13.1 Add **2.1 µl DAPI stock** to **50 ml 1xPBS**



DAPI

by Thermo Scientific

Catalog #: D1306

13.2 Select desired sections and remove PBS from wells.

13.3 Add **500 µl per well** of DAPI stain to each well for **02:00:00**

13.4 Rinse 2 times **00:05:00** **00:05:00** with **500 µl 1x PBS**

Mount and coverslip with Vecatshield

14 Mount sections on



Superfrost Plus Microscope Slies

by Thermo Fisher Scientific

Catalog #: 4951PLUS4

14.1 Allow sections to dry for at least **02:00:00**

14.2 Rehydrate slides in  **250 ml 1x PBS** for  **00:30:00**

14.3 Cover-slip slides with  **130 µl vectashield** and allow to fully dry for  **12:00:00**



**VECTASHIELD® Hardset™ Antifade
Mounting Medium**

by Vector Laboratories

Catalog #: **H-1400**



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited