

Oct 22, 2020

Illumina TruSeq Library quantification with qPCR probe method

Kentaro Itokawa¹¹National Institute of Infectious Diseases, Japan

1 Works for me This protocol is published without a DOI.

Kentaro Itokawa
National Institute of Infectious Diseases, Japan

ABSTRACT

A homemade solution for quantification of Illumina library.

Because the method uses a quenched-probe (not intercalator), it is expected to be relatively stable for the difference of size distribution of libraries.

Currently, we confirmed this protocol with a probe for TruSeq type adapter.

PROTOCOL CITATION

Kentaro Itokawa 2020. Illumina TruSeq Library quantification with qPCR probe method. **protocols.io**
<https://protocols.io/view/illumina-truseq-library-quantification-with-qpcr-p-bnpamdje>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 21, 2020

LAST MODIFIED

Oct 22, 2020

PROTOCOL INTEGER ID

43458

MATERIALS

NAME	CATALOG #	VENDOR
UltraPure [®] Salmon Sperm DNA Solution	15632011	Thermo Fisher
PhiX Control v3	FC-110-3001	Illumina, Inc.
PrimeTime [®] Gene Expression Master Mix	1055770	Integrated DNA Technologies

ABSTRACT

A homemade solution for quantification of Illumina library.

Because the method uses a quenched-probe (not intercalator), it is expected to be relatively stable for the difference of size distribution of libraries.

Currently, we confirmed this protocol with a probe for TruSeq type adapter.

Preparing calibrators and primer/probe

- 1 Prepare low-TE buffer (10 mM Tris-HCl, 0.1 mM EDTA, pH 8.0) containing 5 ng/μl salmon sperm DNA.

1 M Tris-HCl, pH8.0	100 μl
0.5 M EDTA	10 μl

UltraPure™ DNA Solution 10 mg/mL	5 µl
milli-Q water	9889 µl
Total	10 mL

Store in -20 °C.

- 2 Prepare dilution series of PhiX control with low-TE buffer containing 5 ng/µl salmon sperm DNA in 8-strip PCR tubes. These are used as calibrators.

50 pM	100 µl
5 pM	100 µl
0.5 pM	100 µl
0.05 pM	100 µl

Store in -20 °C.

- 3 Prepare 20x primers & probe mix solution as below.

Probe	5 µM
P5 primer	5 µM
P7 primer	5 µM

Store in -20 °C.

Double-quenched probe ([IDT](#)): /56-FAM/ACACTCTTT/ZEN/CCCTACACGACGCTCTTC/3IABkFQ/

P5 primer: 5'-AATGATACGGCGACCACCGA-3'

P7 primer: 5'-CAAGCAGAAGACGGCATACGA-3'

- 4 Prepare the following qPCR master mix per reactions for the number of your samples (+ calibrators).

Prime Time PCR master mix	10 µl
20x primers & probe mix	1 µl
H2O	8 µl

Distribute above to each well.

Add 1 µl of template (diluted library or calibrator) to each well.

- 5 Conduct PCR with the following condition.

95 °C for 2 min

30 cycles of



97 °C for 5s

62 °C for 10s

68 °C for 15s (photo with a filter appropriate for SYBR-green or FAM)

Result

6



An example of expected result

Calculate the original concentration of your library.