



May 05, 2021

Static Glucose-stimulated Insulin Secretion (GSIS) Protocol: Mouse Islets

Aliya F Spigelman¹, Jocelyn E Manning Fox¹, Patrick E Macdonald¹¹University of Alberta**1** Works for me dx.doi.org/10.17504/protocols.io.sp7edrn

CIRTNR2FIC

Aliya Spigelman
University of Alberta

ABSTRACT

Static glucose stimulated insulin secretion (GSIS) protocol for mouse islets.

DOI

dx.doi.org/10.17504/protocols.io.sp7edrn

PROTOCOL CITATION

Aliya F Spigelman, Jocelyn E Manning Fox, Patrick E Macdonald 2021. Static Glucose-stimulated Insulin Secretion (GSIS) Protocol: Mouse Islets. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.sp7edrn>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 17, 2018

LAST MODIFIED

May 05, 2021

PROTOCOL INTEGER ID

14815

MATERIALS TEXT

MATERIALS

[BSA Sigma](#)

Aldrich Catalog #A7906

[Sodium bicarbonate Sigma](#)

Aldrich Catalog #S5761

[HEPES Fisher](#)

Scientific Catalog #BP310-500

[Penicillin-Streptomycin Gibco - Thermo](#)

Fisher Catalog #15140122

[RPMI 1640 Gibco - Thermo](#)

Fischer Catalog #11875

[FBS \(Canadian Origin\) Gibco - Thermo](#)

Fischer Catalog #12483-020

[Sodium Chloride Fisher](#)

Scientific Catalog #BP358-212

[Potassium Chloride Sigma](#)

Aldrich Catalog #P9541

[Calcium Chloride Sigma](#)

Aldrich Catalog #C4904

[Magnesium Chloride Hexahydrate Emd](#)

Millipore Catalog #MX0045

[D-\(-\)-Glucose Sigma](#)

Aldrich Catalog #G8270

[STELLUX® Chemi Rodent Insulin ELISA](#)

Jumbo Alpco Catalog #80-INSMR-CH10

Day Before Secretion

- 1 Mouse islets isolated as described in [Mouse Islet Isolation](#) protocol.
- 2 Pick mouse islets into Mouse Islet Culture Media until as close as possible to 100% purity.

Mouse Islet Culture Media

500ml RPMI 1640 (11.1mM glucose)	Gibco 11875-119
50ml FBS Canadian Origin	Gibco 12483-020
5ml Pen/strep	Gibco 15140-122

- 3 Culture islets overnight in incubator at 37°C, 5% CO₂.

Solution Preparation

4 Prepare solutions as follows:

A	B	C
	mM Final	per 100mL total
NaCl	115	5.75mL (2M)
KCl	5	0.5mL (1M)
NaHCO₃	24	0.2g
CaCl₂	2.5	0.25mL (1M)
MgCl₂	1	0.1mL (1M)
HEPES	10	1mL (1M)
BSA	0.1% w/v	0.1g

KRBH is made fresh on the day of the experiment

Mix chemicals from above table in milliQ water (approximately **80 mL**). Warm KRBH solution to **37 °C** in incubator 37°C, 5% CO₂ (Approximately 30 min). Once solution is warmed, pH KRBH solution to **pH7.4** with NaOH and bring to volume (**100 mL**). KRBH Solution should be kept in incubator throughout the experiment.

Add glucose and/or additional treatments as required.

A	B	C
per 50mL total	from 1M stock	from powder
2.8mM	140μL	0.025g
16.7mM	835μL	0.150g

1M glucose stock should be made fresh on day of experiment

Acid Ethanol:

A	B
95% Ethanol	150mL
Acetic Acid	47mL
Concentrated HCL	3mL

This solution can be made in advance

Experimental Protocol

- Pick islets into 35mm non-tissue cultured coated (NTCC) dish and 'wash' islets with **2 mL** of KRBH with low glucose (2.8mM).
- Pick islets into new 35mm NTCC dish in **2 mL** of low glucose KRBH and pre-incubate in incubator at **37 °C** ^{1h} 5% CO₂ for **01:00:00**.
- Transfer islets into a new 35mm NTCC dish and add **2 mL** of low glucose (2.8mM) KRBH and pre-incubate for ^{1h} **01:00:00** at **37 °C**, 5% CO₂.
- Pick 15 islets into eppendorf tubes. Each treatment group should be done in triplicate.
 - Control group is typically 2.8mM glucose for low glucose and 16.7mM glucose for high glucose
 - Number of islets can be increased/decreased depending on ELISA kit sensitivity

- 9 Gently add **500 µl** of low glucose KRBH to the islets, and incubate for **01:00:00** at **37 °C**, 5 % CO₂. Leave tube lids open. 1h
- 10 Close lids, gently invert tubes, and centrifuge at **1000 rpm, 00:01:00** to pellet islets. 1m
- 11 Collect as much of the **500 µl** supernatant as possible without disturbing the pellet. Store supernatant at **-20 °C** until insulin assay.
- 12 Gently add **500 µl** of high glucose (16.7mM) KRBH to islets, and incubate for **01:00:00** with conditions at **37 °C**, 5 % CO₂. Leave lids open. 1h
- 13 Close lids, gently invert tubes, and centrifuge at **1000 rpm, 00:01:00** to pellet islets. 1m
- 14 Collect as much of the **500 µl** supernatant as possible without disturbing the pellet. Store supernatant at **-20 °C** until insulin assay. **-20 °C**
- 15 Add **500 µl** of acid ethanol to the islets. Store tube at **-20 °C** until insulin assay.

ELISA

- 16 Samples are assayed using ALPCO Stellux Rodent Insulin ELISA kit (Cat # 80-INSMR-CH10). Content samples are diluted with zero buffer 1:400.