



DEC 27, 2022

OPEN ACCESS

DOI:

[dx.doi.org/10.17504/protocols.io.ewov1oodylr2/v1](https://doi.org/10.17504/protocols.io.ewov1oodylr2/v1)

Protocol Citation: Sarah M Prostack, Edgar M Medina, Erik Kalinka, Lillian Fritz-Laylin 2022. Protocol 4: Creating depressions in induction media plates. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1oodylr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 02, 2022

Last Modified: Dec 27, 2022

PROTOCOL integer ID: 72193

Keywords: Spizellomyces, induction media

Protocol 4: Creating depressions in induction media plates

In 1 collection

Sarah M Prostack¹, Edgar M Medina¹, Erik Kalinka¹, Lillian Fritz-Laylin¹

¹University of Massachusetts at Amherst



fritzlaylinlab.umass

ABSTRACT

Transformation is carried out by co-culturing Agro and Spizellomyces on induction media ("IM"). Once its virulence genes are induced, Agro will infect Spizellomyces and the plasmid the transfer DAN is randomly incorporated into the Spizellomyces genome. The IM media and agar plates (see recipes) should be no more than 1 month old, as the hormone Acetosyringone is likely to degrade, decreasing its effectiveness. For each plasmid, four ratios of Agro and Sp will be made and plated onto one plate (See Protocol 5.3). Prepare IM plates with small depressions using a glass culture tube to keep the co-culture liquids from touching each other.

ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

MATERIALS TEXT

Materials

- Induction Media agar plates (1.5% agar w/v; see recipe)
- Large diameter glass culture tube (>10mm, such as



PYREX™ Reusable Borosilicate Glass Tubes with Plain End **Fisher Scientific Catalog #14-957N**

)

- 50 mL conical (such as



Centrifuge Tubes-Bag 50mL Centrifuge Tube - Bag Sterile **Cell Treat Scientific Products Catalog #229421**

)

- 70% ethanol to maintain sterility
- Open flame source to maintain sterility

SAFETY WARNINGS



Protective eye goggles must be worn during this procedure, as the potential for broken or cracked glass is a hazard.

BEFORE START INSTRUCTIONS

The protocol below outlines how to prepare the depressions and should be practiced until the agar is no longer scarred. Prepare these plates within a month of your intended Spizellomyces transformation time.

Steps

1 On the bottom, divide an IM plate into quadrants.

2 Draw circles 2.5 mm in diameter in each quadrant.

Note

Ensure the circles do not overlap.

3 Place a 20-30 mm diameter glass culture tube into a 50 mL conical tube.

4 Add enough 70% ethanol to cover the first inch of the culture tube.



5 Burn off the ethanol on the tube using a flame, let the glass cool in the sterile area around the open flame.


6 Dip the culture tube back into the ethanol.

7 Repeat steps 3-5 a total of 3 times.

8 Working within the sterile area around the flame, gently press the warm (but not enough to melt the agar), sterile culture tube into the agar at the center of one circle drawn earlier.

- 9 Gently press on the agar and move the culture tube in circles until a depression is formed roughly the size of the circle drawn.

Note

- Ensure the four circles of one plate do not overlap and can hold about  200 µL of liquid.
- Depressions can be seen by holding the plate to light at an angle; the light should bend around the circle.

- 10 Repeat for all quadrants of all IM plates needed.