



AUG 01, 2023

Agarose gel for RNA and DNA

kaylee Beine¹

¹University College of Dublin



kaylee Beine

ABSTRACT

General method for running agarose gels for RNA and DNA samples

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.5qpvo3bm7v4o/v1

Protocol Citation: kaylee Beine 2023. Agarose gel for RNA and DNA. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo3bm7v4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited









Protocol status: Working
We use this protocol and it's working

Created: Aug 01, 2023



Last Modified: Aug 01, 2023

Preparing the agarose gel and plate

4m 30s


- 1  1 g agarose in  100 mL of  TAE (Tris-Acetate-EDTA) buffer, 1x Contributed by users
- 2 Microwave solution in a flask for  00:03:00 until the agarose is dissolved. Microwave for  00:01:00 thereafter  00:00:30 4m 30s
- 3 Allow to cool before adding  1 μ L  Ethidium bromide, 10mg/mL, 10mL Amresco Catalog #X328-10ML in the gel medium before completely cool. Then pour into the gel mold
- 3.1 Make sure that your plate is balanced and the teeth are in. Decide if you need a long ladder or a short weather.

Preparing the chamber

- 4 Add  5 μ L  Ethidium bromide, 10mg/mL, 10mL Amresco Catalog #X328-10ML in the gel medium within the chromatography chamber

Preparing RNA samples

5m

- 5  250-300 ng /RNA per sample $C1.V1 = C2.V2$

V1 = (10)(50)/RNA conc. (from nanodrop)

Add the V1 sample of RNA to 1:1 ratio of



RNA Loading Dye (2x) - 4.0 ml New England Biolabs Catalog
#B0363S

to each sample

6

Heat up the ladder



RiboRuler High Range RNA Ladder Thermo Fisher Catalog
#SM1821

5m

and RNA samples to 70 °C for 00:05:00 for the gel

Preparing DNA samples

7



SyberSafe DNA Gel Stain Invitrogen - Thermo Fisher Catalog
#S33101

is already in the

samples for the PCR; then use the



GeneRuler 1 kb Plus DNA Ladder Thermo Fisher Catalog
#SM1331

Setting up the sample gels

1h 50m

8

Add 5 µL of ladder and samples to gel wells

9

For RNA: 80V for 01:00:00

1h 40m

For DNA: 80V for 00:40:00

10

*Note: can check gel after 00:10:00 to see if there is product showing on the gel

10m