



Feb 25, 2022



CIDC_S16_NMR_Celegans_Extraction_Protocol

Brianna M Garcia¹, Amanda Shaver¹, Goncalo Gouveia¹¹University of Georgia

2

dx.doi.org/10.17504/protocols.io.b2rbqd2n

Brianna Garcia
University of Georgia

A sample preparation protocol for lyophilized *C. elegans* samples to be analyzed by NMR spectroscopy

DOI

dx.doi.org/10.17504/protocols.io.b2rbqd2n

Brianna M Garcia, Amanda Shaver, Goncalo Gouveia 2022.
CIDC_S16_NMR_Celegans_Extraction_Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b2rbqd2n>






protocol ,

Dec 09, 2021

Feb 25, 2022

55811

 zirconia beads **BioSpec****Products Catalog #11079110zx** Isopropanol **Fisher****Scientific Catalog # A4614** Methanol **Fisher****Scientific Catalog #A4564** Water **Fisher****Scientific Catalog #W64**

- 1 Remove lyophilized *C. elegans* samples from  **-80 °C** freezer
- 2 Add approximately  **200 µL** of  **1.0 mm** Zirconia beads (BioSpec Cat. No. 11079110zx) to each sample.

3

1m 30s



FastPrep-96





High-throughput bead beating grinder and lysis system

MP Biomedicals**116010500****QuickFlex™ Sample Holder**

Adapter for 96 x 2 mL tube holder on FastPrep-96

MP Biomedicals**116010570**

Using a FastPrep-96™ instrument (MP Biomedicals) equipped with QuickFlex™ Sample Holder (or similar equipment) homogenize samples at  **420 rcf** for  **00:01:30**

- 4 Place samples on dry ice for  **00:01:30** to avoid overheating due and sample degradation. 1m 30s
- 5 Using a FastPrep-96™ instrument (MP Biomedicals) equipped with QuickFlex™ Sample Holder (or similar equipment) homogenize samples at  **420 rcf** for  **00:01:30** 1m 30s
- 6 Place samples on dry ice for  **00:01:30** to avoid overheating due and sample degradation. 1m 30s

- 7 Using a FastPrep-96™ instrument (MP Biomedicals) equipped with QuickFlex™ Sample Holder (or similar equipment) homogenize samples at **420 rcf** for **00:01:30** 1m 30s
- 8 Store on dry ice and proceed with Extraction (Section 2) or store at **-80 °C** until ready to **00:00:00** process.

Sequential non-polar (i) and polar(ii) extraction 7m 30s

- 9 Add **500 µL** of 100% IPA chilled to **-20 °C** to each homogenized sample containing: (1) lyophilized *C. elegans* material and (2) zirconia beads
- 10 Vortex each sample for **00:00:30** to **00:01:00** 1m 30s
- 11 Let sample sit at **Room temperature** for **00:15:00** to **00:20:00** 35m
- 12 Add an additional **500 µL** of 100% IPA chilled to **-20 °C** to each homogenized sample containing: (1) lyophilized *C. elegans* material, (2) zirconia beads, and (3) **500 µL** IPA.
- 13 Vortex each sample for **00:00:30** to **00:01:00** 1m 30s
- 14 Let sample sit at **Room temperature** for **00:15:00** to **00:20:00** 35m
- 15 Store samples **Overnight** (~12 hours) at **-20 °C**
- 16 Centrifuge samples for **00:30:00** at **20800 rcf, 4°C** 30m

17 Transfer via Pasteur pipette the supernatant of each centrifuged sample to a new 2 mL tube for the analysis of non-polar molecules.

18 Place all non-polar extracts in a CentriVap at **Room temperature** and monitor until completely dry.

CentriVap
Benchtop Centrifugal Vacuum Concentrator
Labconco 7810010

19 Add **1 mL** of 80:20 Methanol:Water (MeOH:H₂O) chilled to **4 °C** to each tube containing (1) the remaining worm pellet and (2) zirconia beads.

20 Shake samples at **4 °C** for **00:30:00** 30m

21 Centrifuge samples for **00:30:00** at **20800 rcf, 4 °C** 30m

22 Transfer via Pasteur pipette the supernatant of each centrifuged sample to a new 2 mL tube for the analysis of polar molecules.

23 Place all polar extracts in a CentriVap at **Room temperature** and monitor until completely dry.

24 Remaining worm pellet and zirconia beads can be stored at **-80 °C** for future protein and/or carbohydrate analysis or discarded appropriately.