



FEB 08, 2024

Descending Platform

daniel.dautan^{1,2}, Per Svenningsson^{1,2}¹Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

Eileen Ruth Torres

Weill Cornell Medicine

ABSTRACT

Behavioral test to assess motor function.

MATERIALS

Horizontal grid: 0.5cm space between grid, 45cm long x 5cm wide (Custom made)

Clean cage

OPEN  ACCESS**Protocol Citation:** daniel.dautan,
Per Svenningsson 2024.

Descending Platform.

protocols.io

<https://protocols.io/view/descending-platform-c8wxzxfn>**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited**Protocol status:** Working

We use this protocol and it's working

Created: Feb 08, 2024**Last Modified:** Feb 08, 2024**PROTOCOL integer ID:** 94903

Keywords: ASAPCRN, behavior, motor

Funders Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: 020608

- 1 Place the custom-made horizontal grid (0.5cm space between grid, 45cm long, 5cm wide) to have an angle of 45 degrees with the floor.
- 2 A clean cage was placed at the bottom of the grid.
- 3 Place a camera on top of the grid to allow recording of the time to descend.
- 4 Naïve mice were placed on the horizontal end of the grid and gently pushed to descend the grid towards the cage.
- 5 Quantification:
 - Time to descend the first and second half of the grid
 - Using the video and a proper scale, the distance between the hindlimbs of the mice was determined at 5 randomly selected frames and used as average.