



Jul 09, 2021

# Chromosomal DNA extraction from Gram-positive bacteria

Anders Kiledal<sup>1</sup>, Julia A Maresca<sup>2</sup><sup>1</sup>University of Delaware, Department of Biological Sciences;<sup>2</sup>University of Delaware, Department of Civil and Environmental Engineering

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.bs7knhkw](https://dx.doi.org/10.17504/protocols.io.bs7knhkw)

CivilMicroLab

Julia Maresca

## ABSTRACT

Extraction of high-molecular-weight DNA from Gram-positive bacterial species, with optional steps for removing surfactants. This DNA is suitable for sequencing and the protocol can be scaled up at least 5-fold. Modified from a protocol by Tina Wecke, [LMU-Munich](#).

## DOI

[dx.doi.org/10.17504/protocols.io.bs7knhkw](https://dx.doi.org/10.17504/protocols.io.bs7knhkw)

## PROTOCOL CITATION

Anders Kiledal, Julia A Maresca 2021. Chromosomal DNA extraction from Gram-positive bacteria.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bs7knhkw>

## KEYWORDS

Gram-positive, chromosomal DNA, Firmicutes, Actinobacteria, sequencing

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## IMAGE ATTRIBUTION

Julia Maresca, University of Delaware

## CREATED

Mar 10, 2021

## LAST MODIFIED

Jul 09, 2021

## PROTOCOL INTEGER ID

48076

## GUIDELINES

Recommend wearing gloves throughout and working in a biosafety cabinet if possible to prevent contamination. This protocol can be scaled up at least 5-fold. If the isopropanol precipitation step is used, subsequent steps do not have to be scaled up unless the culture volume is substantially larger.

## MATERIALS TEXT

### SOLUTIONS

#### TEN

- 10 mM Tris-HCl, pH 8.0
- 10 mM EDTA
- 150 mM NaCl

#### TEN\*

- 10 mM Tris-HCl, pH 8.0
- 1 mM EDTA
- 50 mM NaCl

#### RNAse A

- 20 mg/mL in water

#### Lysozyme

- 20 mg/mL in water

#### SDS

- 10% (w/v) in water

#### Other reagents:

Isopropanol, ethanol (100% and 70%), phenol, chloroform:isoamyl alcohol (24:1), sterile water.

Enzyme solutions should be stored at -20 between uses or prepared freshly. Other solutions can be stored at room temperature.

### CONSUMABLES

- microcentrifuge tubes (or larger centrifuge tubes, depending on volume)
- pipetment (P1000, P200, P20)
- pipet tips (P1000, P200, P20)
- Glass Pasteur pipet with tip bent

### SAFETY WARNINGS

Phenol and chloroform:isoamyl alcohol should be handled in a fume hood and the liquid waste and contaminated tubes should be disposed of in accordance with the institution's rules for handling organic solvent waste.

### BEFORE STARTING

Grow culture to high cell density and prepare all solutions.

#### Grow culture

- 1 Inoculate 10 mL rich medium from a fresh overnight culture, and incubate at appropriate temperature on shaker. At OD600 of ~0.8-1.0, harvest cells by centrifugation (10 min., 5000 rpm).

#### Cell lysis

- 2 Resuspend cell pellet in 2 mL TEN (10 mM Tris-HCl, pH 8.0, 10 mM EDTA, 150 mM NaCl).
- 3 Add 100 µL lysozyme (20 mg/mL) and incubate for 20 min at 37°C.
- 4 Add 20 µL RNAse (10 mg/mL) and incubate for 3 min at 65°C.

#### Remove surfactants (optional)

- 5 IF THE STRAIN PRODUCES A SURFACTANT THAT INTERFERES WITH THE PHASE SEPARATION, Add 0.1 volume 3 M sodium acetate and 1 volume cold isopropanol, mix, and incubate on ice for 20 min. Centrifuge for 10 min at 5000 rpm and decant the supernatant. Then resuspend in 400  $\mu$ L TEN and 550  $\mu$ L TEN\* and transfer to a microcentrifuge tube.

#### Phenol & chloroform:isoamyl alcohol extractions

- 6 Add 900  $\mu$ L phenol, mix by inversion. Centrifuge for 5 min. at 13000 rpm and transfer the upper phase to a clean microcentrifuge tube.
- 7 Re-extract once with phenol (1 volume) and twice with chloroform: isoamyl alcohol (24:1 v/v, 1 volume)

#### DNA precipitation

- 8 Collect DNA by coiling on the end of a glass Pasteur pipet.
- 9 Air dry, then resuspend DNA in 100  $\mu$ L sterile water overnight at 4°C.