


May 29, 2024

## Cryo-EM sample preparation for RCKW:DARPin complex

 Forked from [Preparation of LRRK2 RCKW cryo-EM grids](#)

DOI

**[dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1)**

Marta Sanz Murillo<sup>1</sup>

<sup>1</sup>University of California, San Diego

ASAP Collaborative Rese...



Marta Sanz Murillo

University of California, San Diego

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1)**

**Protocol Citation:** Marta Sanz Murillo 2024. Cryo-EM sample preparation for RCKW:DARPin complex. **protocols.io**  
**<https://dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 31, 2021

**Last Modified:** May 29, 2024

**Protocol Integer ID:** 100824

**Keywords:** cryo-EM, LRRK2, structural biology, ASAPCRN

**Funders Acknowledgement:**

**ASAP**

**Grant ID:** ASAP-000519



## Abstract

This is Leschziner's Lab protocol for making cryo-EM grids for RCKW:DARPin complex.

## Materials

LRRK2 Buffer:

- [M] 20 millimolar (mM) HEPES pH 7.4
- [M] 150 millimolar (mM) NaCl
- [M] 0.5 millimolar (mM) TCEP
- [M] 5 % volume Glycerol
- [M] 2.5 millimolar (mM) MgCl<sub>2</sub>
- [M] 20 micromolar (μM) GDP

Note: please change salt as needed to maintain final salt of 150 mM NaCl

## Safety warnings

- ! For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).  
Take proper precautions while freezing grids.


## Before start

Decide which protein concentration to use, and create the proper LRRK2 buffers in order to obtain the right salt concentration (150 mM NaCl).



## Freezing Grids

20s

- 1 Plasma clean grids.  
We used UltrAuFoil Holey Gold 1.2/1.3 300 mesh grids and plasma cleaned them in the Solarus II (Gatan) using the QuantiFoil Au preset.
- 2 Dilute samples to desired concentration in the **LRRK2 buffer**. Make sure final salt is at 150 mM NaCl.  
For best results, make  8  $\mu$ L samples, good for freezing 2 grids. This is to minimize time spent outside of storage buffer, reducing aggregation.
- 3 Apply protein to grids and plunge freeze (3-4  $\mu$ L)  
We used a Vitrobot (FEI) to blot away excess sample and plunge freeze
- 4 Store grids in liquid nitrogen until ready for imaging.