



Dec 09, 2020

Fluorescent Western Protocol

Steven J Burgess¹¹University of Illinois at Urbana-Champaign

In Development

This protocol is published without a DOI.

Ag SynBio Lab UIUC

Steven Burgess
University of Illinois at Urbana-Champaign

ABSTRACT

Analysis of proteins using fluorescent immunoblot.

Note:

- The choice of secondary antibody depends on the choice of primary antibody, whether it is derived from a mouse (monoclonal) or a rabbit (polyclonal).
- It is advisable to stick to the 800CW wavelength to avoid problems with chlorophyll autofluorescence encountered with the 680CW antibodies.

Literature: <https://www.licor.com/documents/fxc6evxvxbub4srkqy6i9yg46l7i0xz5>

PROTOCOL CITATION

Steven J Burgess 2020. Fluorescent Western Protocol. **protocols.io**
<https://protocols.io/view/fluorescent-western-protocol-bqhwmt7e>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2020

LAST MODIFIED

Dec 09, 2020

PROTOCOL INTEGER ID

45334

MATERIALS TEXT

- Intercept^R PBS Blocking Buffer (LI-COR Biosciences; [927-70001](#)) - milk powder with PBS could equally be used
- IRDye[®] 800CW Donkey anti-Rabbit IgG Secondary Antibody (LI-COR Biosciences; [926-32213](#))
- IRDye[®] 800CW Donkey anti-Mouse IgG Secondary Antibody (LI-COR Biosciences; [926-32212](#))
- Black Western Blot Incubation Box (LI-COR Biosciences; [929-97110](#))
- 10x PBS buffer
- TweenTM 20 (Fisher Biosciences; [BP337-100](#))
- Methanol
- Odyssey CLx Imager (LI-COR Biosciences)
- Primary antibody ([various](#))

ABSTRACT

Analysis of proteins using fluorescent immunoblot.

Note:

- The choice of secondary antibody depends on the choice of primary antibody, whether it is derived from a mouse (monoclonal) or a rabbit (polyclonal).

- It is advisable to stick to the 800CW wavelength to avoid problems with chlorophyll autofluorescence encountered with the 680CW antibodies.

Literature: <https://www.licor.com/documents/fxc6evxvxbub4srkqy6i9yg46l7i0xz5>

1 

Keep membranes in the black Western Blot incubation box for all steps, this is important after adding the secondary antibody because the signal is light-sensitive and will become bleached if not exposed to light for a long enough period.

2 Wet with 1x PBS for  **00:02:00** min

2m

3 Rinse membrane with dH₂O

4  **Intercept (PBS) Blocking Buffer LI-**




1h

Discard PBS and incubate with **COR Catalog #927-70001**

 **01:00:00**

at  **Room temperature**


As an alternative to Intercept Blocking Buffer you can use PBS-T 5% w/v milk powder. This is to prevent unspecific binding of antibody and lowers background signal

5 Incubate with primary antibody (appropriate dilution in Intercept Blocking Buffer) at  **4 °C**  **Overnight** (in cold room) with gentle agitation on a platform shaker  **50 rpm, 4 °C** .

6 Prepare  **1 L** PBS-T solution

7 Pour off the primary antibody and rinse the membrane with PBS-T.

Some primary antibodies can be re-used multiple times depending on the concentration used, in this instance collect the primary antibody in a tube and store the solution at -20 °C before re-use.

8 Cover the membrane with PBS-T, shake vigorously on a platform shaker at  **50 rpm, Room temperature , 00:10:00** . Repeat 3 times


9 Create a working dilution of secondary antibody using



[☒ Intercept \(PBS\) Blocking Buffer LI-](#)

COR Catalog #927-70001

with 0.2 % (v/v) Tween, 0.01 % SDS (w/v). IRE

Dye secondary antibody (1:20,000).

 **10 mL** is sufficient to cover one blot and 0.01% SDS is only to be used in conjunction with PVDF membranes.

10 Incubate for  **01:00:00**  **Room temperature** with gentle agitation on a platform shaker. 1h

11 Pour off the secondary antibody

12 Cover the membrane with PBS-T, agitate  **80 rpm, Room temperature , 00:10:00**

13 Discard PBS-T. Repeat step 12 three times  **go to step #12**

Note: more washes (x5) and for longer can be done to reduce background

14 Rinse then cover the membrane with 1x PBS,  **80 rpm, Room temperature , 00:10:00** . Repeat (2X)

15 Proceed to imaging blot on LI-COR Odyssey CLx imaging system

Odyssey CLx
Imaging System

LI-COR

Odyssey CLx

