



SEP 23, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocol  
s.io.e6nvwd6x2lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwd6x2lmk/v1)

**Protocol Citation:** Elias  
Adriaenssens 2023. GFP pull  
down assay. **protocols.io**  
[https://dx.doi.org/10.17504/p  
rotocols.io.e6nvwd6x2lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwd6x2lmk/v1)

**License:** This is an open  
access protocol distributed  
under the terms of  
the [Creative Commons  
Attribution License](#), which  
permits unrestricted use,  
distribution, and reproduction  
in any medium, provided the  
original author and source  
are credited

**Protocol status:** Working  
We use this protocol and it's  
working

**Created:** Jun 27, 2023

**Last Modified:** Sep 23,  
2023

**PROTOCOL integer ID:**  
84086

**Keywords:** GFP pull down  
assay

## GFP pull down assay

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

### ABSTRACT

This protocol describes GFP pull down assay.

### ATTACHMENTS

[754-1922.pdf](#)

### MATERIALS

#### Materials

- GFP-Trap agarose beads (Chromotek)
- dH<sub>2</sub>O
- Protein Loading dye


#### Bead assay buffer

A	B
Tris-HCl pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM



## GFP pull down assay

2h 5m

- 1 Mix GFP-tagged TBK1 with 20  $\mu$ L of equilibrated GFP-Trap agarose beads (Chromotek) at a final concentration of 1 micromolar ( $\mu$ M). Make sure to wash beads 2x in dH<sub>2</sub>O before washing with bead assay buffer to equilibrate the beads adding the protein to the beads.

2 To this end, wash  20  $\mu\text{L}$  of beads twice with  $\text{dH}_2\text{O}$  and equilibrate with bead assay buffer.



3 Resuspend the beads in  40  $\mu\text{L}$  bead assay buffer, to this add GFP-TBK1 at a final concentration of  5 micromolar ( $\mu\text{M}$ ) .



4 Incubate the beads with GFP-TBK1 for  01:00:00 at  4  $^{\circ}\text{C}$  at a horizontal tube roller.





1h

5 Wash the beads three times to remove unbound GFP-tagged bait protein.






6 Prepare protein master mixes with prey protein in bead assay buffer at the following concentrations:

- mCherry-OPTN (1  $\mu\text{M}$ ),
- mCherry-NDP52 (1  $\mu\text{M}$ ),
- GST-NAP1 (1-10  $\mu\text{M}$ ).

7 Add the protein master mixes to the beads and incubate for  01:00:00 at  4  $^{\circ}\text{C}$  at a horizontal tube roller.



1h

8 Wash the beads three times to remove unbound proteins, remove any supernatant from the beads and resuspend the beads in  60  $\mu\text{L}$  of 1x Protein Loading dye, and heat-inactivate at  95  $^{\circ}\text{C}$  for  00:05:00 .



5m

9 Analyze the samples by SDS-PAGE and Coomassie staining as described above.

