



Version 2

Jan 22, 2021

cDNA synthesis using SuperScript™ IV V.2

Roey Angel¹, Eva Petrova¹¹Soil and Water Research Infrastructure

1

Works for me

dx.doi.org/10.17504/protocols.io.brqum5ww

SoWa RI Anaerobic and Molecular Microbiology (public)

Tech. support email: eva.petrova@bc.cas.cz

Roey Angel

Soil and Water Research Infrastructure

ABSTRACT

The following protocol is intended as a downstream application for our [Purification of RNA from a DNA/RNA Extract](#) protocol. This protocol describes how to synthesise a first-strand non-specific complementary DNA (cDNA) from a purified RNA extract using [SuperScript IV Reverse Transcriptase](#). The second strand synthesis is usually not required for most downstream applications. This protocol is a simplified and condensed version of the full protocol provided by the manufacturer.

ATTACHMENTS

[SSIV_First_Strand_Synthesis_System_UG.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.brqum5ww

PROTOCOL CITATION

Roey Angel, Eva Petrova 2021. cDNA synthesis using SuperScript™ IV. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.brqum5ww>

Version created by Roey Angel

KEYWORDS

cDNA, RNA, mRNA, SuperScript, reverse transcriptase

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 22, 2021

LAST MODIFIED

Jan 22, 2021

PROTOCOL INTEGER ID

46580

MATERIALS TEXT

MATERIALS

[Bovine Serum Albumin \(BSA\)](#) **Thermo Fisher**

Scientific Catalog #B14

Step 4

[SuperScript™ IV First-Strand Synthesis System](#) **Thermo Fisher**

Scientific Catalog #18091050

Step 4

[RNaseOUT™ Recombinant Ribonuclease Inhibitor](#) **Thermo Fisher**

Scientific Catalog #10777019

Step 4

[Ribonuclease H \(RNase H\)](#) **Thermo Fisher**

Scientific Catalog #18021071

In 2 steps

[DNA Polymerase I \(10 U/μL\)](#) **Thermo Fisher**

Scientific Catalog #EP0041

Step 7

[Random hexamers](#) **Thermo**

Scientific Catalog #N8080127

Step 1

STEP MATERIALS

[Random hexamers](#) **Thermo**

Scientific Catalog #N8080127

Step 1

[Nuclease-free autoclaved DEPC-treated water](#) **Carl**

Roth Catalog #T143.1

In 2 steps

[USB Dithiothreitol \(DTT\) 0.1M Solution](#) **Thermo Fisher**

Scientific Catalog #707265ML

Step 4

[RNaseOUT™ Recombinant Ribonuclease Inhibitor](#) **Thermo Fisher**

Scientific Catalog #10777019

Step 4

[SuperScript™ IV First-Strand Synthesis System](#) **Thermo Fisher**

Scientific Catalog #18091050

Step 4

[Bovine Serum Albumin \(BSA\)](#) **Thermo Fisher**

Scientific Catalog #B14

Step 4

[Ribonuclease H \(RNase H\)](#) **Thermo Fisher**

Scientific Catalog #18021071

In 2 steps

[DNA Polymerase I \(10 U/μL\)](#) **Thermo Fisher**

Scientific Catalog #EP0041

Step 7

[Ribonuclease H \(RNase H\)](#) **Thermo Fisher**

Scientific Catalog #18021071

In 2 steps

[Nuclease-free autoclaved DEPC-treated water](#) **Carl**

Roth Catalog #T143.1

In 2 steps

BEFORE STARTING

Make sure your RNA is pure and contains no traces of DNA. A simple and very sensitive way to ensure that is to use the purified RNA as a template for a PCR reaction targeting a gene that should be present in the sample. A negative result indicates a lack of DNA template in the sample.

Primer annealing



1 

Prepare the following mixture in a PCR tube:

1. **1 µl to 4 µl purified RNA** (**10 pg - 5 µg** ; usually **200 ng** for soil extract)
2. **1 µl random hexamers (50 µM)** or a gene-specific primer (**2 µM**)
3. **9.8 µl RNase free water**

 [Random hexamers Thermo](#)

Scientific Catalog #N8080127

 [Nuclease-free autoclaved DEPC-treated water Carl](#)

Roth Catalog #T143.1

2 

Mix gently and spin down the solution.

3 

Incubate the mixture at **65 °C** for **00:05:00** in a thermocycler and chill **On ice** (or in the cyclor at **> 4 °C**) for at least **00:01:00** .

Reverse transcription

4 


Prepare the following mixture and add to each tube:

1. **4 µl 5x Reaction buffer**
2. **1 µl dNTP mix, 10 mM**
3. **1 µl 0.1M DTT**
4. **1 µl RNaseOUT™ (40 U/µl) ***
5. **0.2 µl BSA (20 µg/µl)**
6. **1 µl SuperScript™ IV RT (200 units/µl)**

* Optional

 [USB Dithiothreitol \(DTT\) 0.1M Solution Thermo Fisher](#)

Scientific Catalog #707265ML

 [SuperScript™ IV First-Strand Synthesis System Thermo Fisher](#)


Scientific Catalog #18091050

 [RNaseOUT™ Recombinant Ribonuclease Inhibitor Thermo Fisher](#)

Scientific Catalog #10777019

 [Bovine Serum Albumin \(BSA\) Thermo Fisher](#)

Scientific Catalog #B14

 dNTP Mix (10 mM each) Thermo Fisher

Scientific Catalog #R0191

5 

Incubate the mixture in a thermocycler at **23 °C** for **00:10:00** (only if using random hexamers, skip if using a specific primer) followed by **50 °C** for **01:00:00** to **03:00:00** and **80 °C** for **00:10:00**. Chill **On ice**.

6 

For PCR templates > 1kb remove the RNA by adding **1 µl** (2 units) of E. coli RNase H and incubate at **37 °C** for **00:20:00**.

 Ribonuclease H (RNase H) Thermo Fisher

Scientific Catalog #18021071

Optional: Second strand synthesis

7 

Prepare the following mixture and add to each tube:

1. **1 µl DNA Polymerase I reaction buffer**
2. **0.75 µl DNA Polymerase I**
3. **0.2 µl RNase H**
4. **3.05 µl RNase-free water**
5. **5 µl Template cDNA**

 DNA Polymerase I (10 U/µL) Thermo Fisher

Scientific Catalog #EP0041

 Ribonuclease H (RNase H) Thermo Fisher

Scientific Catalog #18021071

 Nuclease-free autoclaved DEPC-treated water Carl

Roth Catalog #T143.1

8 

Incubate for at **15 °C** for **02:00:00** followed by **00:10:00** at **75 °C** for deactivation.

9 3. Purify the reaction through phenol/chloroform purification followed by ethanol precipitation or using a PCR purification kit.